



**100Plus remote patient monitoring devices assist physicians and our staff to monitor patients by providing informative data virtually. The data from these devices has enabled our practice to care for our patients when they are at risk and that's what matters most to our practice.**

— Dr. Liberoni



## Regulatory And Safety Notices

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

Modifications not expressly approved by 100Plus could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### Radiofrequency radiation exposure information

This router has been tested and meets the FCC RF exposure guidelines when used with the 100Plus accessories supplied or designated for this product. Use of other accessories may not ensure compliance with FCC RF exposure guidelines.

**100plus**

## Health Hub

### Quick start guide



### We're here to help

+1 (844) 445-8267

support@100plus.com

## Intended Use

The intended use of the device is to relay data from certain supported Bluetooth enabled medical sensors to our secure portal for physician review. The user is expected to plug this device into a power outlet somewhere in their home where it is likely to get good cell reception (e.g. near a window) while still close to the sensors they use (1-2 non structural walls in between is usually OK).

## Included Accessories

Power adapter

## Visual Alerts

### Slow flashing yellow

The hub is trying to connect to the network

### Steady green

The hub is connected to the network and idle

### Slow flashing blue

The hub is communicating with one or more Bluetooth sensors

### Quick flashing green

The hub is attempting to send data to the network

### Steady red

An irrecoverable hardware failure has occurred (SIM not working/ cell module not working).

## Audio Alerts

### One short beep

The hub has obtained a measurement from a sensor

### Three short beeps in succession

The hub has successfully transmitted a measurement to the network

## Contact 100Plus Customer Support

If you need help with your 100Plus Health Hub, please give us a call or send us an email. We may also send you emails periodically with reports from your doctor to help you manage your chronic condition with 100Plus.

### We're here to help

+1 (844) 445-8267

[support@100plus.com](mailto:support@100plus.com)

**100plus**

