



OR scan the QR code to download the app



3

1

2

The screenshot shows the XLife app interface. At the top, the 'XLife' logo is visible. Below it, there are three icons: a person, a syringe, and a sun. The sun icon is highlighted with a red box. Below these icons, there are two status indicators: 'Remaining Battery' at 99% and 'Remaining Time' at 35 min. At the bottom, there is a home button icon (a house) which is circled in red.



1

2

3

Trouble	Solution
starting up and no response	Please use the attached charging cable for quick charging. The charging time is 35 minutes
Camera does not turn on	Close and open the cover to try a few times
APP quit quickly	Delete the existing APP and re-scan the QR code for installation
Installation of tools in accessory barrel is not smooth	Unscrew the ear spoon head and pull out the silicone ear spoon for replacement. Screw and plug the ear spoon head well
Rod's power indicator is red	Please use the attached charging cable for charging. The charging time is 35 minutes
The display on the app is not reading properly.	Disconnect the phone with product and open cover to restart it, then set phone to reconnect product
WiFi connection is unstable, sometimes it is disconnected in seconds	Close and open the cover to restart rod and restart app in phone to connect rod
If you still have problems with the above procedures, please contact our Customer Service Representatives directly for further assistance	

Part name	Hazardous material					
	(Pb)	(Hg)	(Cd)	(Cr(VI))	(PBB)	(PBDE)
Printed Circuit Components	×					
Case	○					
Plastic rubber accessory						

Seller: Telephone number



Enable everyone to enjoy visible ear cleaning

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

RF Exposure Information

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.