


# BACtrack Skyn Operator Manual

Item# SK-15  
FCC ID: S8B-SK-15  
IC: 10935A-SK15

## Powering ON/OFF the BACtrack Skyn



On the top of the BACtrack Skyn, press and hold the power button for two seconds. The button will turn blue and remain on for about a minute, indicating the device is ON.

To turn the device off, press and hold the power button for two seconds. The button will fade from blue to gray, indicating the device is OFF.

The BACtrack Skyn can still be ON when the button is not blue. To check if the device is still ON, firmly tap the top of BACtrack SKYN. If the button pulses blue, the device is still ON and recording data. If the button remains gray, the device is OFF.

## Wearing the BACtrack Skyn



1. Slide the BACtrack Skyn over the fingers and palm of your non-dominant hand.
2. Place the device approximately 1 inch below your wrist joint, making sure that the white sensor circle is facing the underside of your wrist (ie. the opposite position of a normal watch).
3. Tighten the strap. Make sure it is tight, but not uncomfortable.

## Connecting to the BACtrack Skyn App



1. Enable Bluetooth on your Apple iOS device.
2. Download and install the BACtrack Skyn App from the Apple iTunes App Store.
3. Open the BACtrack Skyn App and power your device ON.
4. The App will automatically connect to your Skyn and you will see a confirmation message that you've successfully connected your device.

Additional product information, instructions and troubleshooting can be found in the App.

## Charging the BACtrack Skyn



1. Make sure the BACtrack Skyn is powered OFF by pressing and holding the power button for two seconds.
2. Attach the charging cable to the bottom of the device by lining up the two magnetic squares on the cable and device.
3. Insert the USB connector into a USB port to begin charging. When the BACtrack Skyn is connected to the charger the light will begin pulsing green.

It takes approximately 3 hours to fully charge, at which point the light will stop pulsing and remain green, indicating the device is fully charged. The BACtrack Skyn should not be worn while it is being charged.



## Regulatory Compliance

This device (will) comply with Part 15 of the FCC Rules. Operation is subject to the following conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation

NOTE: This equipment (will be) tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected

Consult the dealer or an experienced radio/TV technician for help. The FCC requires the user to be notified that any changes or modifications made to this device that are not expressly approved by KHN Solutions Inc may void the user's authority to operate the equipment.

### RF EXPOSURE WARNING STATEMENT

This equipment (will) comply with FCC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance. The information in this guide may change without notice. The manufacturer assumes no responsibility for any errors that may appear in this guide.

### RSS-GEN, ISSUE 4, NOVEMBER 2014

This device complies with ISED license-exempt RSS standard(s). Operation is subject to the following two conditions:

1. This device may not cause harmful interference;
2. This device must accept any interference received, including interference that may cause undesired operation of the device.

Radiation Exposure Statement: This equipment complies with the IC RSS-102 radiation exposure limits set forth for an uncontrolled environment.

Cet appareil est conforme à Industrie Canada une licence standard RSS exonérée. Son fonctionnement est soumis aux deux conditions suivantes:

1. Cet appareil ne doit pas provoquer d'interférences.
2. Cet appareil doit accepter toute interférence reçue, y compris les interférences pouvant provoquer un fonctionnement indésirable de l'appareil.

Énoncé d'exposition aux rayonnements : Cet équipement est conforme aux limites d'expositions aux rayonnements ioniques RSS-102 pour un environnement incontrôlé.

To access product regulatory information, please visit <https://skyn.bactrack.com/regulatory>