CMS50F User Manual Pulse Oximeter

CONTEC Contec Medical Systems Co., Ltd.

Address:No.112 Qinhuang West Street, Economic & Technical Development Zone, Qinhuangdao, Hebei Province, PEOPLE' S REPUBLIC OF CHINA Tel: +86-335-8015430 Fax:+86-335-8015588 Technical support:+86-335-8015431 E-mail:cms@contecmed.com.cn

Website:http://www.contecmed.com CMS2.782.078ESS/1.1 1.4.01.01.822 2021.10

Instructions to Lise

Dear users, thank you very much for purchasing the Pulse Oximeter (hereinafter refe It is a medical device, which can be used repeatedly,

The Manual describes, in accordance with the Pulse Oximeter's features and requirements, main structure functions, specifications, correct methods for transportation, installation, usage, operation, repair, maintenance and storage, etc. as well as the safety procedures to protect both the user and equipment. Refer to the respective chapters for details.

Please read the User Manual carefully before using this product. The User Manual which describes the preating procedures should be followed strictly. Failure to follow the User Manual may cause measuring abnormality, equipment damage and human injury. The manufacturer is NOT responsible for the safety,

reliability and performance issues and any monitoring abnormality, human injury and equipment damage due to users' negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

Owing to the forthcoming renovation, the specific products you received may not be totally in accordance with the description of this User Manual. We would sincerely regret for that.

Our company has the final interpretation to this manual. The content of this manual is subject to change without prior notice. Warnings:

Remind that it may cause serious consequences to tester, user or environmen

- Explosive hazard-DO NOT use the device in environment with inflammable gas such as anesthetic. DO NOT use the device while examining by MRI and CT, as the induced current may cause burn.
- Do not take the information displayed on the device as the sole basis for clinical diagnosis. The device is only used as an auxiliary means in diagnosis. And it must be used in conjunction with doctor's advice,
- clinical manifestations and symptoms. be performed by qualified service personnel specified by manufacturer, dangers (such as over-temperature, fire or explosion) may occur when replacing the battery by the personnel not fully trained. Users are not permitted to maintain or refit the device by themselves or replacement of the
- Uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the nicrocirculation disturbance users. It is not recommended that the sensor is used on the same finger for more than 2 hours.
- For some special users who need a more careful inspection on the test site, please don't place the device n the edema or tender t
- Please do not stare at the red and infrared light emitter (the infrared light is invisible) after turning on First or the device, including the maintenance staff, as it may be harmout on the device, including the maintenance staff, as it may be harmful to the eyes.
 The device contains silicone, PVC, TPU, TPE and ABS materials, whose biocompatibility has been
- tested in accordance with the requirements in ISO 10993-1, and it has passed the recommended biocompatibility test. The person who is allergic to silicone, PVC, TPU, TPE or ABS can not use this
- The disposal of scrap device, its accessories and packaging should follow the local laws and regulations. to avoid polluting to the local environment. And the packaging materials must be placed in the region where the children are out of reaching.
- The device can not be used with the equipment not specified in the Manual. Only the accessories appointed or recommended by the manufacturer can be used, otherwise it may cause injury to the tester nd operator or damage to the device.
- The SpO2 probe accompanied is only suitable for using with the device. The device can only use the SpO2 probe described in the Manual, so the operator has the responsibility to check the compatibility between the device and the SpO2 probe before using, incompatible accessories may cause device performance degradation, device damage or user injury.
- Do not reprocess the accompanying SnO: probe-
- Check the device before use to make sure that there is no visible damage that may affect user's safety
- and device performance. When there is obvious damage, please replace the damaged parts before use. When the message "Sensor Off" or "Sensor Fault" appears on the screen, it indicates that the SpO₂ probe is disconnected or line fault occurs. Check the connection of the SpO₂ probe and whether there is damage for the probe, if necessary, please replace the probe to avoid risks. The probe fault will not result in a safety hazard.
- Functional testers can not be used to assess the accuracy of the SpO2 probe and Pulse Oximeter.
- Some functional testers or national simulators can be used to verify whether the device works normally or example, INDEX-2LFE Simulator (software version: 3.00), please refer to the Manual for the detailed operation steps.
- Some functional testers or patient simulators can measure the accuracy of the device copied calibration curve, but they can not be used to evaluate the device accuracy.
- When using the device, please keep it away from the equipment which can generate strong electric field or strong magnetic field. Using the device in an inappropriate environment may cause interference to he surrounding radio equipment or affect its working
- When storing the device, keep it away from children, pets and insects to avoid affecting its performance
- Do not place the device in places exposed to direct sunlight, high temperature, humidity, dust, cotton wool or easy to splash water, to avoid affecting its performance.
- The measured accuracy will be affected by the interference of electrosurgical equipment When several products are used on the same people simultaneously, danger may occur which is arisen
- from the overlap of leakage current.
- Solution of the operation of the set of t
- The intended operator of the device may be a user.
- Avoid maintaining the device during using
 Users should read the product manual carefully before use and operate according to the requirements

1.Overview Insert the finger when measuring the device will directly display the SnO+ value measured, it has a high

aracy and repeatability 1.1 Features

- Easy to use
- Small in volume, light in weight, convenient to carry.
- Low power consumption.

1.2 Intended purpose

The Pulse Oximeter can be used in measuring the pulse oxygen saturation and pulse rate through finger. The product is suitable for being used in family, hospital, oxygen bar, community healthcare, physical care in sports (It can be used before or after doing sports, and it is not recommended to use the device during the process of having sport) and etc. 1.3 Environment Requirements

- Storage Environment
- Temperature : -40°C~+60°C
- b) Relative humidity : ≤95%
 c) Atmospheric pressure : 500hPa~1060hPa
- Operating Environment

- a) Temperature: +10 °C~ + 40 °C b) Relative Humidity : <75 %
- Atmospheric pressure: 700hPa-1060hPa 1.4 Precautions
 - 1.4.1 Attention Point out conditions or practices that may cause damage to the device or other properties.
 - Before using the device, make sure that it locates in normal working state and operating environment.
- In order to get a more accurate measurement, it should be used in a quiet and comfortable environment. When the device is carried from cold or hot environment to warm or humid environment, please do not use it
- immediately, wait four hours at least is recommended.
- If the device is splashed or coagulated by water, please stop operatin DO NOT operate the device with sharp things. High temperature, high pressure, gas sterilizing or immersion disinfection for the device is not permitted. Refer
- to User Manual in the relative chapter (6.1) for cleaning and disinfection. Please turn off the device before cleaning and disinfection. The device is suitable for children and adult.
- The device may not be suitable for all users, if you can't get a satisfactory result, please stop using it
- The device may not be suitable (of an above, n you can get a substactory result, prease stop using it. Data averaging and signal processing have a delay in the upgrade of SpO₂ data values. When the data update period is less than 30 seconds, the time for obtaining dynamic average values will increase, which is arisen from signal degradation, low perfusion or other interference, it depends on the PR value,
- The device has 3-year service life, date of manufacture sees the label. The expected service life of the attached parts or accessories of the equipment is two year
- If the shelf life is less than the expected service life, the shelf life of the attached parts or accessories of the
- equipment is two year. This device has the function of prompting, users can check on this function according to chapter 5.3.1 as a
- The device has the function of limits prompting, when the measured data is beyond the highest or lowest limit
- The device has the function of minis prompting, which the inclusive data is beyond the ingress on lowest minit, the device would start prompting automatically on the premise of the prompting function is on. The device has the function of prompting, this function can either be paused, or closed (default setting) for good This function could be turned on through menu operation if you need. Please check the chapter 5.3.1 as a
- The maximum temperature at the SpO₂ probe -tissue interface should be less than 41°C which is meas the temperature tester
- E During measuring, when abnormal conditions appear on the screen, please pull out your finger and reinsert it to measure again.
- If some unknown error appears during measuring, press "RESET" button to reset it.
- Do not contort or drag the wire of the device. The plethysmographic waveform is not normalized, as a signal inadequacy indicator, when it is not smooth and stable, the accuracy of the measured value may degrade. When it tends to be smooth and stable, the measured value read is the optimal and the waveform at this time is also the most standard.
- The device can not be used during charging. If necessary, please visit our official website to get the information about SpO_2 probe that can be used with this device. If the device or component is intended for single-use, then the repeated use of these parts will pose risks on the
- parameters and technical parameters of the equipment known to the manufacturer. f necessary, our company can provide some information (such as circuit diagrams, component lists,
- illustrations, etc.), so that the qualified technical personnel of the user can repair the device components designated by our company. The measured results will be influenced by the external colouring agent (such as nail polish, colouring agent or
- color skin care products, etc.), so don't use them on the test site.
- As to the fingers which are too cold or too thin or whose fingernail is too long, it may affect the measured results, so please insert the thicker finger such as thumb or middle finger deeply enough into the probe when measuring
- The finger should be placed correctly (see Attached figure 4), as improper installation or improper contact The high should be placed correctly (see Allached figure 4), as improper instantion or improper contact position for sensor will influence the measurement. The light between the photoelectric receiving tube and the light-emitting tube of the device must pass through
- the subject's arteriole. Make sure the optical path is free from any optical obstacles like rubberized fabric, to void inaccurate results
- Excessive ambient light may affect the measured results, such as surgical light (especially xenon light sources), bilirubin lamp, fluorescent lamp, infrared heater and direct sunlight, etc. In order to prevent interference from ambient light, make sure to place the sensor properly and cover the sensor with opsque material. Frequent movement (active or passive) of the subject or severe activity can affect the measured accuracy
- The SpO₂ probe should not be placed on a limb with the blood pressure cuff, arterial ductus or intraluminal
- The measured value may be inaccurate during defibrillation and in a short period after defibrillation, as it has not defibrillation function.
- The device has been calibrated before leaving factory.
- The device is calibrated to display functional oxygen saturation. The device is calibrated to display functional oxygen saturation. The equipment connected with the Oximeter interface should comply with the requirements of IEC 60601-1 Please select medical nower adapter to charge it when connecting the special adapter with the socket, make sure there is no shelter near the socket and it is easy to plug and unplug, otherwise the power will not be cut off in time when necessary, causes damage.

1.4.2 Clinical restriction A. As the measure is taken on the basis of arteriole pulse, substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to

- nterference. B. The measurement will be influenced by intravascular staining agents (such as indocvanine green or methylene blue), skin pigmentation.
- C. The measured value may be normal seemingly for the tester who has anemia or dysfunctional hemoglobin(such as carboxyhaemoglobin (COHb), methaemoglobin (MetHb) and sulfhaemoglobin (SuHb)), but the tester may appear

hypoxia, it is recommended to perform further assessment according the clinical situations and symptoms. D. Pulse oxygen only has a reference meaning for anemia and toxic hypoxia, as some severe anemia users still show better pulse oxygen measured valued. E. Contraindication:

- a The person who is allernic to silicone. PVC TPU TPE or ABS can not use this device
- b. The damaged skin tissue can't be measured.
- c. During cardiopulmonary resuscitation
- d. When the user is hypovolemic.
- e. For assessing the adequacy of ventilatory support.
 f. For detecting worsening lung function in users on a high concentration of oxygen.
- 1.5 Clinical indications
- The Pulse Oximeter can be used in measuring the pulse oxygen saturation and pulse rate through finger.

2 Principle

Principle of the Oximeter is as follows: An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive Hemoglobin (Hb) and Oxyhemoglobin (HhO₂) in glow & near-infrared zones. Oneration principle of the device is: Photoelectric Oxyhemoglobin Inspection (inc) in how the mathematical operation principle on the device is a how receive of mathematical operation of the provided of the second principle of measured signal can be obtained by a photosensitive element, information acquired through which will be shown on screen through treatment in electronic circuits and microprocessor.



Figure 1. Operating principle



- B. PR value and bar graph display
- C. Pulse waveform display D. Low-battery indication: low-battery indication appears when the battery voltage is too low to work

. Voice prompt for over-limit, sensor off /finger-out and low battery.

K. The data can be unloaded to the terminal equipment by wired mode

E. Automatic standby function

The data can be uploaded to the terminal equipment (Bluetooth wired equipment) by wireless mode

Figure 2 Appearance

Phaneniphoist file moniphi

01:28.8 6.0

PRbom

Figure 3 Measurement interface

Open the USB plug of the device, insert the micro end of USB cable into the SpO2 probe interface , the other end

B. Accessories: one SpO2 probe, one USB cable, one power adapter (optional), one User Manual, Bluetooth adapter

ليهج لاده

Serve Serve

Dot1

Figure 5 Main menu

"Direction": direction, "up": increase the value, "down": decrease the value

The values displayed in Figure 6 are the initial values of over-limit prompt.

1.....

Ward and

Dated with

"SpO2 HI": upper limit prompt for SpO2 over-limit

"PR HI": upper limit prompt for PR over-limit

"Pulse Sound": PR sound, "off": close, "on": open,

and return to Main Menu interface

measurement interface flickers.

turn on / off memory by "Record".

last time.

main menu

space is full.

switching mode

5.3.3 Clock setting

shown in Figure 8.

display "Recording".

"Seo": data segment

button again to enter the measurement interface.

power, pulse sound indication will turn off automatically.

Incl. Time

Total - Manual

device adopts 24-hour clock.

"set year": set the year "set month": set the month

"set minute": set the minute

"set day": set the day

"set hour": set the hour

as shown in Figure 9.

"ID": user name

"Hand Ver." hardware version

"Soft.Ver.": software version

that Minute.

Table

Figure 8 Time setting interface

"set time": set the time. "yes": to set the time. "no": not set the time

menu: if the time is not set, it directly return to the main menu.

"Brightness": set the screen brightness, adjustable range: 1~4

5.3.4 System setting and other options introduction

"Delete All": delete all records (auto record mode is shown as Figure 7).

5.3.2 Data storage

SpO2 LO": lower limit prompt for SpO2 over-limit

"PR LO": lower limit prompt for PR over-limit "Prompt Sound": prompt for over-limit, "off": close, "on": open.

Acres 1 Nov.

(99 group of data at most), the total duration does not exceed 72 hours Manual record: store up to 24-hour data.

5.3.1 Sound prompt setting

interface shown in Figure 6.

-11

Phi lithuai

Walter Tenned

Figure 6 Setting interface for sound prompt

Under main menu, short press the button to select "Sound", then long press the button to enter its setting

Table Sound 111 control of 1000cr and 1000cr and 1000cr limit can not be lower than the lower limit when adjusting the values. SpO₂ range: $0 \% \sim 100 \%$, PR range: $0 \sim 254$ bpm

After setting, short press the button to select "Exit", then long press the button to exit sound setting interface,

Under the main menu, short press the button to select "Record", then long press the button to enter the

Record Menu interface as shown in Figure 7. It indicates that the device is storing when the red dot "Ro" in

Figure 7 Record menu interface

"Mode": record mode selection, including: "Auto" and "Manual" mode. Under "Manual" mode, select to

Auto record: start recording after stable data annear, null out the finger to finish recording a group of data

When the memory is full, it will display "Memory is full!", then it will enter the standby mode after several

seconds. When exiting the standby mode, it will display "Memory is full!" to give user a prompt, press the

Note: under manual mode, when "Record" is "ON", the device will prompt to clear the data stored

It will display "Recording" when there is no operation under record state for 30s, then it will enter energy

saving mode after several seconds, long press the button to exit this mode; short press the button, it will

Note: under data recording state, after the display screen turns off automatically, in order to save

After setting, short press the button to select "Exit", long press the button to exit record menu and return to

Note: please upload data in time after recording, otherwise the data may be covered when the storage

Note: the historical data will be deleted once switching the record mode. Under recording state, the

record mode can not be switched; under manual mode, the "Record" should be turned off before

Under main menu, short press the button to select "Clock", long press the button to enter its sub-menu as

Short press the button to select the option to be adjusted, then long press the button to change the value. The

After setting, short press the button to select "Exit", then long press the button to exit clock setting interface

If user changes the setting, current set time will be displayed, then press the button to return to the main

Under main menu, short press the button to select "System", then long press the button to enter the interface

Short press the button to select the option to be adjusted, then long press the button to change the value.

"Domo": set the Demo mode, "on": turn on the Demo mode, "off": turn off the Demo mode. "Sound Volume": set the sound volume, adjustable range: 1 ~ 3

1000

10

100

Section Water

Anthe IEE

Figure 9 System menu interface

10000110001

front fighter

the Later same

Short press the button to select the option to be adjusted, then long press the button to change the value

Short press the button to select the option to be adjusted, then long press the button to change the value

Annual brand will

Tritt

Acres 1 Name

24.14

-

- 14

F. display mode can be changed 3. Adjustable screen brightness H. PR sound indication J. Memory function

M. With clock function

and in such

4.2 Interface introduction

USB interface: connect with USB cable or SpO2 probe

RESET hole: reset button inside of it, to reset the device.

facts only

Number of Contraction of Contract, Name

by their store.

the second

4.3 SpO₂ probe installation

4.4 Connection of USB cable

into computer or power adapter.

C.Software description Software name: CMS50F embedded software

Purpose: be used to measure SpOy, pulse rate, etc

A. Insert the finger into the probe as shown in Figure 4.

4.5 Structure and accessories

Software specification: r

Release version: 2.0

5 Operating

fingernail.

it will resume automatically after about 60s.

5.1 Measurement

(optional).

Button: power on, pause sound prompt, display clock, enter menu, menu operation.

4.40

O

5.5eDt

Open the USB plug of the device, then insert the SpO2 probe interface into the USB port of the device

4.5 Structure: main unit, SpO₂ probe, USB cable, power adapter (optional) and Bluetooth adapter (optional).

Please check the device and accessories according to the list to avoid that the device can not work normally

linical function: calculate SpO₂ and pulse rate values by collecting and processing the testee's pulse signa

Figure 4 Sketch map for finger placement

(The appearance of actual probe may be different with the one shown as Figure 4, please refer to the actual probe.)

Note: when inserting the finger, the light emitting from the sensor must be directly irradiated to the side of the

5.2Measurement interface A.In sound prompt "ON" state, when the sound prompt occurs, short press the button to pause the sound prompt, and

Under the measurement interface, long press the button to enter the main menu interface as shown in Figure 5, sound,

B-Press button to turn on the device, it displays the measurement interface. C.Wait a few seconds, the device directly shows measurement result on the screet

Note: during measuring, do not shake the finger and keen quiet, not move

B.If you want to turn off the sound prompt permanently, please set it in menu,

record, clock,system and Bluetooth, etc. can be set, methods are as followings

ound prompt includes over-limit prompt, low-battery prompt, probe off or finger out prompt.

If there is no sound prompt, short pressing the button will directly enter clock interface, press again to exit 5.3 Menu operation

Naming rule for version: V <Maior enhancive software upgrade>.<Minor enhancive software

upgrade>.<Improvement software upgrade> Involved algorithm: name: plethysmography; type: mature arithmetic

....

N. Charging function

4 Installation

4.1 Appearance

After setting, short press the button to select "Exit" ,then long press the button to exit the system menu and return to main me 5.3.5 Bluetooth setting

Under main menu, short press the button to select "Bluetooth" then long press the button to enter its election interface as shown in Figure 10 and Figure 11. When the Bluetooth is "ON" jf there is no data to be transmitted, the Bluelooih will turn off automatically after a few seconds. Note: under data transmitting state by Bluetooth, the Bluetooth can not be turned off.



Figure 10 Bluetooth "ON" interface Figure 11 Bluetooth "OFF" interface

5.3.6 Shutdown Under the main menu interface, short press the button to select "Power off" then long press the button to n off the device-

5.3.7 Exil main menu Under main menu, short press the button to select "Exit", then long press the button to exit the main menu

and return to the measurement interface. 5.4 Data unload

5.4.1 Wired transmission

Connect the device to computer by the USB cable, upload the data after connecting with the PC software

properly, refer to "Software operating instruction" for details. Note: the PC software can be downloaded from our official website

5.4.2 Bluetooth wired transmis

Turn ou the device Bluetooth and the PC software to upload data, refer to "Software operating instruction" for details

5.5 Charging Power adapter can be selected to charge for lhe device.

It indicates that the device is charging when the indicator is orange, the charging is finished when the indicator turns to green.

5.6 Reset

Use a pointed and hard object (for example, a paper clip) to press lhe reset button inside of the RESET hole, to reset the device.

stream, 10 per

The device must be turned off before cleaning, and it should not be immersed into liquid. Please take out the internal battery before cleaning ,do not immerse it into liquid. Use 75% alcohol to wipe the device enclosure, and use liquid soap or isopropanol to wipe the watchband for disinfection, nature dry or clean it with clean and soft cloth. Do noJ spray any liquid on the device directly-

and avoid liquid penetrating into the device. A.Check the main unit and all accessories periodically to make sure that there is no visible damage that may

affect user's safety and monitoring performance. It is recommended lhai the device should be inspected weekly at least. When there is obvious damage, stop using it.

B.Please clean and disinfect the device before/after using it according to the User Manual (6.1).

D Recharge the battery in lime when low battery appears. D Recharge the battery soon after over-discharge. The device should be recharged every three months when

it is not used for some time. It can extend the battery life following this guidance

E.The device need not to be calibrated during maintenance. 6.3 Transportation and Storage

A. The packed device can be transported by ordinary conveyance or according io transport contract. During traiLsportation, avoid strong shock, vibration and splashing with rain or snow, and it can not be transported mixed with toxic harmful material

B. The packed device should be stored in room with no corrosive gases and good venlilation. Temperature: -40°C~60°C; Humidity: 595

Contraction in the local division of the loc				_		
Trouble		Possible Reason		Solution		
The values can not be displayed normally or stably		1. The finger is not properly positioned. 2. The finger is shaking or the user is moving. 3. The device is not used in environment required by the manual. 4. The device works abnormally.		 Please insert the finger properly and measure again. Lel the user keep calm Please use the device in normal Please contact the after-sales. 		
The device can not be turned on		 Low battery or the battery is drained away. The device works abnormally. 		1-Please charge the battery. 2-Please contact the after-sales.		
The display disappears suddenly.		1.The device enters into the energy saving mode. 2.Low battery. 3.The device works abnormally.		1.P 2.P 3.P	1. Please contact the after-sales. 2. Please charge the battery 3. Please contact the after-sales.	
The device can not be used for full time after charge.		1. The battery is not charged fully. 2. The device works abnormally.		1.Please charge the battery. 2.Please contact the after-sales.		
The battery can not be full charged even after 10 hours charging lime.		The battery works abnormally.		Please contact the after-sales.		
The data can not be stored.		1. The device is not operated according to the manual. 2. The device works abnormally.		 Please operate the device according to the manual. Please contact the after-sales. 		
Road Contractor	-	_	_			
Symbols		Meaning		ls	Meaning	
©	Cauti accor	on, consult npanying documents	PRbpn	1	Pulse rate (bpm)	
(b)	29	The Women of		E	Pulse oxygen saturation (%)	
.	Manufacturer		i		Fully charged	
SN	Serial mimheir		12		Use-by date	

R	Recycling garbage WEEE (2012/19/EU)	•4	USB	
+	Battery anode	-	Battery cathode	
IP22	It means this pulse is protected against harmful effects of dripping water when tilled at 15°.	RST	RESET hole	
1	Temperature limitation	15	Humidity limitation	
0	Atmospheric pressure limitation	II	This way up	
r	Fagin lack of an	Æ	Internet Second	
Û	Locates.	8	Close the sound prompt	
8	Pause lhe sound prompt	0	Open the sound prompt	
O- 山	Menu/ Power button	dx	Close ihe PR sound	
۵	Recyclable	3	Open the PR sound	
0'	Bluetooth icon (Bluetooth wired device)	E.	Manufaclure Dale	
Sensor Off	The probe is disconnected.	Finger Out	The finger is not inserted.	
R*	Recording	Sensor Fault	Probe failure	
P/N	Material code	LOT	Batch No.	
	1. The finger clip falls off 2. Probe error 3. Signal inadequacy		Alarm inhibit	

idicator

sPo2 see note 1

Display range

Resolution

Measured range

Resolution

note 41

Light interference

Pulse Intensity

Optical sensor see note 51

SnO -

Red ligh

Memory

Safety class

International Protection Working voltage

Working current

Power supply

Infrared light

PR

Accuracy see note 31

Accuracy under low perfusion[see

Upper and lower limit of measured values

PR Display range

Measured range

Accuracy see note 2]

Note : Your device may not contain all the following symbols,

0%_00%

0%~100%

1%

70%-100%; ±2%

30bpm 250bpra

30bpm~250bpm

Low perfusion 0.4%

Same and

%

stronger pulse.

0% 10.0%

0 bpm - 254 bpm

not exceed 72 hours

DC 36 V = 42 V

1P22

< 100 mA

Up to 24-liour data under manual mode

Internally powered equipment, type BF applied part

±2 bom during the pulse rate range of 30 bom 99 bom and

PR: ±2 bpm during the pulse rate range of 30 bpm = 99 bpm and ±2% during the pulse rate range of 100 bpm 250 bpm.

der normal and ambient light conditions, the SpO+ deviation

Continuous bar-graph display, the higher display indicates the

Wavelength: about 660 nm, optical output power 1< 6.65 mW

Wavelength: about 905 nm, optical output power: < 6.75 mW

Un to 99 group of data under auto mode, total duration does

A rechargeable lithium battery (3.7V) (The red wire on the

battery denotes anode, the black wire on the bancry denotes

±2% during the pulse rate range of 100 bpm - 250 bpm.

0%-69%: unspecified

cathode) Charge and discharge:no less than 500 times Battery life Output voltage: DC 5V Adapter specification Output current: 1000 mA Dimensions and Weight Dimension 61 mm(L) > 356 mni(W) x 24 mm(H) Weight About 60 g (including a lithium battery)

Note 1: the claims of SpO2 accuracy shall be supported by clinical study measurements taken over the full range. By artificial inducing, get the stable oxygen level to the range of 70 % to 100 % SpO 2compare the SpO palues collected by the secondary standard pulse oximeter equipment and the tested equipment at the time, to form paired data, which used for the accuracy analysis (ii is applicable for the probes equipped.

There are 12 healthy volunteers (male: 6, female: 6; age: 18-50; skin color: black: 2, light: 8, while: 2) data in the clinical report.

Note 2: because pulse oximeter equipment measurements are statistically distributed, only about hvo-thirds of pulse oximeter equipment measurements can be expected to fall within ±Amis of the value measured by a CO-OXIMETER.

Note 3: Patient simulator has been used to verily the pulse rate accuracy, it is stated as the root-mean-square difference between the PR measurement value and the value set by simulator-

Note 4: percentage modulation of infrared signal as the indication of pulsating signal strength patient simulator has been used to verify ils accuracy under conditions of low perfusion. SpO₂ and PR values are different due to low signal conditions. compare Ihem with the known SpO: and PR values of input signal. Note 5: optical sensors as the light-emitting components, will affect other medical devices applied the

wavelength range. The information may be useful for the clinicians who carry out the optical treatment. For example, photodynamic therapy operated by clinician.

State	Prompt condition delay	Prompt signal generation delay
Low voltag cprompt	Is	20ms
SpO: prompt	330ms	20ms
Pulse rate prompt	330ms	20ms
Probe error prompt	16ms	20ms

EMC

abic 1:		
Guidance a	nd manufacturer's declaration -ele	ctromagnetic emission
The Pulse Oximeter is intended for user of the device should assure th	or use in the electromagnetic envir at it is used in such environment.	onment specified below. The purchaser or Ihe
Emission test	Compliance	
RF emissions CISPR 11	Group :	
able 2: Guidance a	nd manufacturer's declaratiou-elec	tromagnetic immunity
The Pulse Oximeter is intended fo user of ihe Pulse Oximcler should	or use in the electromagnetic envir assure lhat it is used in such enviro	onment specified below. The purchaser or lhe nment.
Immunity test	IEC60601 test level	Compliance level
Electrostatic discharge (ESD)	±8kV contact	±8kV contact
IEC 61000-4-2	±15 kV air	±15kV air
Power frequency (50 / 60Hz) magnetic field IEC 61000-4-8	30 A/m	30A/m

Table 3:

The Pulse Oximeter is	intended for use in the electromagnetic enviro	onment specified below. The customer
user of the Pulse Oxim	eter should assure that it is used in such environ	nent.
Immunity test	IEC 60601 test level	Compliance level
Radiated RF IEC61000-4-3	10 V/m 80 MHz= 2.7 GHz	10 V/m80 MHz- 2.7 GHz
NOTE : AI 80 MHz	and 800 MHz, the higher frequency range applic	s.

IOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by

absorption and reflection from structures, objects and people. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios amateur radio. AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic emironment due to fixed RF transmitters, ait electromagnetic site survey should be considered. If the measured field strength in the location in which the Pulse Oximeter is used exceeds the applicable RE compliance level above the Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Pulse Oximeter.

h Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3V/m.

he [Code ser of the	SI] is inten Pulse Oxin	ded for use neter should	in the electromag assure that i(is u	netic environme sed in such an e	ent specified nvironment	i below. The c	ustomer or the
Radiated	Test Freque ncy (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modul alien b) (W)	Distance (m)	IMM UNI TY TEST LEVEL (V/m)
	385	380 390	TETRA 400	Pulse modulation b) 18 Hz	K8	0,3	27
	450	430 -470	GMRS 460, FRS 460	FMc) ±5 kHz deviation 1 kHz sine	2	0,3	28
F	710	704 -787	LTE Band 13,17	Pulse			9
C610	745			modulation b) 217 Hz	02	0.3	
0-4-3 Test	780				02	0,3	
ped flea	810		GSM 800/900, TETRA 800, iDEN 820. CDMA 850,				
ons for	870	800		Pulsa			
ENCLO SURE PORT IMMUN	930			modulation b) 18 Hz	2	0,3	28
F	1720		GSM 1800.	Pulse modulation b) 217 Hz		0.3	28
ireless	1945	4U I	CDMA 1900; GSM 1900; DECT; LTE Band 1, 3,4,25; UMTS				
mmun	1845	4					
ications equipme nt)	1970	1700 -1990			2		
	2450 5240 5500	2400 -2570 5100 -5800	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7 WLAN 802.11 a/n	Pulse modulation b) 217 Hz Pulse modulation b)	2	0,3	28
_	5785			217 Hz			-
OTE If atenna an cnnitled l	necessary t d the ME by IEC 6100	o achieve t EQUIPMEN 10-4-3.	he IMMUNITY	TEST LEVEL, EM may be rea	the distant	The I m t	e transmitting est distance is
For som The carr As an al	e services, o ier shall be ternative to	modulated FM modula	nk frequencies at using a 50 % duly ttion. 50 % pulse	e included. cycle square w modulation at 1	ave signal. 8 Hz may b	e used becaus	e while it does
ot rq/rese	nt actual m	odulation, it	would be worst o	ase.			
ne MAN	NAGEME!	SK should c	onsider reducing ing higher IMM	the minimum se UNITY TEST	paration dis	tance, based o hat are appro	opriate for the
educed EVELS s	se hall be ailc	paration dis ulated using	stance. Minimun the following eq	uation: $E = \prod_{i=1}^{n}$	stances for	higher IMM	UNITY TEST
/here P	is the max	imum pow	er in W, d is Il	a ne minimum se	paration dis	stance in ra,	and E is the

arning

Table 4

- Don't near active HFSURGICAL EQUIPMENT and the KF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.
- ÷ Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or -
- decreased electromagnetic immunity of this equipment and result in improper operation." Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device including cables specified by Che manufacturer. Otherwise- degradation of Ihe performance of this equipment could result.
- Active medical devices are subject to special EMC precautions and likey must be installed and . used in accordance vilh these guidelines

A Note :

When the device is disturbed, (he data measured may fluctuate ,please measure repeatedly or in another environment to ensure its accuracy.

Bluetooth Specification Working frequency: 2402 MHz 2480 MHz Modulation mode: GFSK Transmitting power:-6 dBm, +4 dBm Receiving sensitivity: -93 dBm

FCC Caution

§ 15.19 Labeling requirements.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

§ 15.21 Information to user.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

§ 15.105 Information to the user.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

-Reorient or relocate the receiving antenna.

-Increase the separation between the equipment and receiver.

-Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

-Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable/mobile exposure condition without restriction

.....