

尺寸68x181mm

正面



Promote the era of intelligent ear cleaning

Open bird C3


1

Scan QR code to download app



2

Open Bebird app and select Bebird C3



3

Press and hold the power key to turn on the ear cleaning rod




The blue light flashes after starting

Press and hold the power key for 1 second for power on

Connect product WiFi


1

Connect Product Wi-Fi "bebird-XXXX"



2

Back to Bebird app Enter and start to use it



Accessory replacement

Bayonet type ear spoon's end It can be replaced by side folding

Directly-inserted end of accessory Insert the buckle



Clipped end of accessory It can be removed by side-folded

Accessory

1

Q-elastic ear spoon × 1
The world's first version which is combined PC with silica gel



2

Transparent ear spoon × 2
Dig the ear canal deeply and take out a large amount of earwax



Product size

PC + silicone
two color injection
Outer diameter of lens tube
4.5mm



Hazardous substance

Part name	Hazardous material					
	(Pb)	(Hg)	(Cd)	(Cr(VI))	(PBB)	(PBDE)
Printed Circuit Components	×	○	○	○	○	○
Case	○	○	○	○	○	○
Plastic rubber accessory	○	○	○	○	○	○

This form is formulated in accordance with SJ/T 11364
○ : represent the hazardous material ratio is below GB/T26572-2011
× : represent the hazardous material ratio is over GB/T26572-2011

背面

Troubleshooting

Trouble	Solution
starting up and no response	Please try several more times according to the correct boot-up method
	Please use your own charging equipment for charging. The charging time is about 1h
Camera does not turn on	Shut down and restart the ear-cleaning rod
APP quit quickly	Delete the existing APP and re-scan the QR code for installation
Rod's power indicator is red	Please use your own charging equipment for charging. The charging time is about 1h
The product failed to display the real-time image after successfully connecting the Bebird APP	Disconnect the phone from WiFi, restart the product and reconnect the WiFi
WiFi connection is unstable, sometimes it is disconnected in seconds	Restart ear clean rod, reconnect Bebird APP
If the above operation still fails to solve the problem, please contact Bebird Customer Service. We will provide you practical solutions.	

Attention

- When cleaning lens, please wipe it carefully with professional alcohol cotton swab.
- Before use it, please check the surrounding environment. Do not use the product in the environment where someone runs, in order to avoid strike.
- This product is not suitable for children under one year old.
- It is forbidden for children to use this product alone in order to avoid accidental harm.
- Do not use the endoscope in liquid to avoid damage caused by immersion.
- Product which is built-in rechargeable lithium batteries, if not often used, must also be charged once a month.
- Avoid product exposure to sunshine directly, especially the accessories in the barrel, in order to avoid softening.
- During using it, the temperature of the product will rise slightly (up to 35 C). Please rest assured that it can be used.
- In order to ensure the best using experience, please keep the child in a sitting position and operate the product horizontally.
- For your ear cleaning experience, please replace the spoon's end after more than 6 months of use.



Product parameters

Brand	bebird
Product Name	Smart visual ear-clean Rod
Model	C3
Rod weight	16g
Network standard	IEEE 802.11b/g/n
Antenna	Internal FPC antenna
Working frequency	2.4Ghz
Image transmission rate	30fps
image sensor	CMOS
Working temperature	-10 ~ 50℃
Battery	250mAh lithium battery
Battery life	about 60 minutes
Charging time	1h
Input current	DC 5V 300mA
Lens diameter	4.5mm
Optimum focal length	1.5 ~ 2cm
Pixel	3M
Gravity Sensor	3 axis



Warranty

User name:	Telephone number:
Address:	
Seller:	Telephone number:



bebird

Promote the era of intelligent ear cleaning

FCC Statement
This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.