

GAUSSRING

User Manual



Please read the Quick Use Guide carefully before use.

1. Products Introduce



Sleep



Body Temperature



Heart Rate



Blood Oxygen



Menstrual cycle




Exercise Record

Heart Rate Monitoring	Resting and Ambulatory Heart Rate Monitoring with Periodic Analysis
Menstrual Cycle Management	Menstrual Cycle Tracking and Personalized Health Insights
Blood Oxygen monitoring	Continuous Monitoring and Analysis of Blood Oxygen
Sleep monitoring	Scientific recording of sleep parameters to help improve sleep quality
Temperature monitoring	Track and Analyze Your Body Temperature to Gain Insights into Your Overall Health
Exercise Record	Step Counting and Exercise Tracking

2. Jimi Care App: Getting Started

- Scan the QR code to download the APP (or search for it in Apple App Store and Android App Market).
- Connect the ring to the power supply through the charging cable, the charging status is detected and the ring will automatically enter the working state.
- Open JIMI CARE APP, click 'My - To add', enter the binding interface, and search the Bluetooth information automatically.



Please scan the QR code to download the Jimi Care APP



2. Jimi Care App: Getting Started

Pairing Your Ring:


1. **Open the app** and look for the Bluetooth devices list.
2. **When you see your ring's name appear**, tap on it to start pairing.
3. **If the ring isn't found**, make sure it's in the charging cradle. The app will automatically try to connect and pair with the ring after a few moments.
4. **During pairing**, the ring will blink green for 3 seconds.

Unpairing Your Ring (iPhone only):

1. Open the iPhone - **Settings**.
2. Tap on **Bluetooth**.
3. Find your ring's name in the list of devices.
4. Tap on the **information icon (i)** next to your ring's name.
5. Tap on **"Forget This Device"** to unpair the ring.

Unpairing Your Ring (Android devices):

1. Open the app and go to **"My Settings"** (or similar section).
2. Look for the **"Unbind Ring"** option and tap on it to complete the unpairing process.



3. Product Usage Guidance


- How to wear:

Wearing:

1. Find your perfect fit: The ring comes in various sizes. Choose the one that fits comfortably on your chosen finger. It shouldn't be too loose or too tight.

SOS Feature:


1. Emergency? Quickly tap the small light sensor on the inner side of the ring, facing your finger, 10 times in a row. This will trigger the SOS function.



3. Product Usage Guidance

- **How to Charge:**

1. Connect the charging cradle to the USB port and place the ring on the charging cradle as shown in the illustration. The ring blinks red to indicate that charging is in progress, and the ring blinks green to indicate that the battery is fully charged.
2. If the ring is charging as shown and the ring shows no light indication, please check the charging contacts of the ring and charger for stains and wipe the contacts clean before charging.



Ring Charging Contacts



Cable Charging Contacts



Align the charging probe and magnet to start charging

4. Safety & Precautions

Before using and operating the equipment, please read and observe the following precautions to ensure optimal performance and to avoid dangerous or illegal situations.

- This device has a built-in non-removable battery, and the battery is not replaceable, please do not disassemble or modify the battery.
- Do not expose the device to high temperatures or around heat-generating equipment, such as heaters, microwave ovens, etc. Do
- not throw the device into a fire, as this may cause the battery to catch fire and explode. If the battery leaks, do not allow your skin or
- eyes to come into contact with the leaked liquid. If it comes into contact with your skin or eyes, flush it with water immediately and go to the hospital for medical treatment. Please dispose of batteries in accordance with local regulations and do not dispose of batteries
- as household garbage, as improper disposal of batteries may cause them to explode. Do not drop, hit, crush or puncture the battery to avoid subjecting the battery to large external pressures
- that could cause internal short circuits and overheating. Children or pets are strictly prohibited from swallowing this product or accessories to avoid injury.
-

4. Safety & Precautions

- This device has a built-in non-removable battery, and the battery is not replaceable, please do not disassemble or modify the battery.
- Some people are allergic to materials such as plastic and stainless steel, and long-term exposure can cause symptoms such as redness, swelling, and inflammation. If similar situations occur, please stop using and consult a doctor. Radio waves generated by
- the device may interfere with the proper functioning of implanted medical devices or personal medical equipment, such as pacemakers, hearing aids, etc. If you use any of these devices, please consult their manufacturers for restrictions on the use of this equipment. Do not dispose of this equipment and its accessories as
- ordinary household garbage. Please observe local ordinances for the disposal of this equipment and its accessories and support recycling initiatives.



4. Safety & Precautions

- NOTE:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.



5. Names and contents of hazardous substances in products

Component Name	Harmful Substance					
	Pb	Hg	Cd	Cr(VI)	PBB	PBDE
Charging Cable	○	○	○	○	○	○
Housings	X	○	○	○	○	○
PCB Assembly	○	○	○	○	○	○
Batteries	○	○	○	○	○	○

- This table is based on the provisions of SJ/T113641.

○: Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit requirements stipulated in GB/T26572-2011.

X: Indicates that the content of the hazardous substance in all homogeneous materials of the part exceeds the limit requirements stipulated in GB/T26572-2011.

6. Product Specification

Brand	JimiloT
Product Name	GAUSSRING
Package includes	GAUSSRING *1 & Charging Cable *1
Charging method	Magnetic wired charging
Batteries	Lithium Polymer Battery
Charging time	≤1.5 hours
Waterproof Level	IP68
Connection Method	Bluetooth LE5.0

JimiloT

7. Warranty Statement

- This product is warranted for six months from the date of purchase for product quality problems caused by manufacturing, materials, design, etc. during normal use. For after-sales service, please
- contact the seller or service provider of the corresponding purchase channel. Free warranty is not provided for faults caused by the user's
- personal reasons, as follows: (1) Faults caused by unauthorized disassembly, modification, etc.; (2) Faults caused by accidental falls during use. (3) All man-made damages or damages caused by third party's negligence, external force cracking, external parts scratching, etc. are not covered by the warranty.

Jimi IoT