

# Medtronic

## Envision™ Recorder

User Guide





EN

English

1

FR

DE

ES

NL

IT

NO

FI

SV

DA

PT

SL

CS

HU

PL

TR

SK

EL

RU

RO

BG

HR

SR

ZH

BP

AR

HE

## Introduction

The Envision™ recorder is a component of the Envision™ Pro Continuous Glucose Monitoring (CGM) system. It is a retrospective CGM system, therefore data is not available to patients in real time. The recorder is compatible with the Envision™ sensor (MMT-7080). The recorder connects to the sensor, and receives and sends data to the Envision™ Pro application (app) through a Bluetooth™ wireless connection. For detailed information on system components consult the Envision™ Pro Continuous Glucose Monitoring System User Guide for Healthcare Professionals.

## Indications

The recorder is intended for single-patient, single-use in patients with diabetes mellitus. The recorder is a component of the Envision™ Pro CGM system.

## Contraindications

None known.

## Warnings

- Always refer to the Envision™ sensor (MMT-7080) user guide for all contraindications, warnings, precautions, and instructions relating to the sensor. Not referring to the sensor user guide can result in serious injury to the patient or damage to the sensor.
- No modification of this equipment is allowed.
- This product contains small parts and may pose a choking hazard for young children.
- Attempting to send data from the recorder when the recorder is near other medical devices that emit radio frequency should be avoided due to possible interference. If you have communication issues, try moving away from such devices.
- Do not expose your recorder to x-ray, ultrasound, or diathermy devices as the performance of the recorder has not been evaluated under those conditions and may be unsafe. If your recorder is exposed to any of these, discontinue use and contact your local country representative for further assistance.
- Do not expose your recorder to MRI equipment or other devices that generate strong magnetic fields as the performance of the recorder has not been evaluated under those conditions and may be unsafe.

If your recorder is inadvertently exposed to a strong magnetic field or ionizing radiation, discontinue use and contact your local country representative for further assistance.

- Do not expose your recorder to temperatures exceeding those listed in the specifications table for Storage Conditions as this may deplete the battery and result in a non-functional recorder.

## **Precautions**

Do not reuse recorders. The recorder is designed to be used for one patient, and one evaluation only. Once the recorder is activated for a patient, it cannot be used for another evaluation or patient. The recorder will not function and no data will be gathered.

## **Assistance**

Please contact your local country representative using the Medtronic Diabetes International Contacts list in this user guide.

## **IEC60601-1-2; Special EMC Precautions for Medical Electrical Equipment**

1. Special precautions regarding Electromagnetic Compatibility (EMC):  
This body worn device is intended to be operated within a reasonable residential, domestic, public or work environment, where common levels of radiated “E” (V/m) or “H” fields (A/m) exist; such as cellular phones, WiFi, Bluetooth™\*, electric can openers, microwave and induction ovens. This device generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the provided instructions, may cause harmful interference to radio communications.
2. Portable and mobile RF communications equipment can affect Medical Electrical Equipment as well. If you encounter RF interference from a mobile or stationary RF transmitter, move away from the RF transmitter that is causing the interference.

## **Installing the Envision™ Pro app**

For information on installing the Envision™ Pro app, consult your Envision™ Pro System User Guide for Healthcare Professionals.

## **Inserting the Envision™ sensor (MMT-7080) and connecting the Envision™ recorder**

For information about inserting the sensor and connecting the recorder,

consult your Sensor User Guide.

For information on pairing the recorder with the app, use the Envision™ Pro app and follow the instructions on your screen.


**Note:** *The graphical image, a type of barcode, on the back of the recorder is for manufacturing purposes only.*

## Ending the evaluation and uploading data

End the evaluation and upload the data according to the Envision™ Pro System User Guide for Healthcare Professionals. The patient or HCP can upload data either before or after removing the recorder and sensor from the patient. Be sure to upload data before recorder disposal.

## Removing the sensor and the recorder

1. Put on gloves.
2. Peel the sensor and recorder off of the body as one unit.
3. Separate the sensor and the recorder.
4. Dispose of the sensor in a sharps container. Dispose of the recorder according to local regulations for battery disposal (non-incineration).



**Note:** *Do not discard the recorder in a medical waste container or receptacle in which it would be exposed to extreme heat, above 55 °C (131 °F).*

## Bathing and swimming

After the recorder and sensor are connected, they form a waterproof seal to a depth of 2.4 meters (8 feet) for up to 30 minutes. Your patient can shower and swim without removing them. No additional tape is required.

## Help

The Envision™ Pro app provides the best source of information for help with the recorder. To access the Help screen, tap **Help**. The app will walk you through the various Help topics.

## Storing and transporting the devices

Store the recorder in a clean and dry location at room temperature

between 15 °C (59 °F) and 30 °C (86 °F). Do not transport the recorder at temperatures above 55 °C (131 °F) or below -30 °C (-22 °F). Temperatures outside this range can damage components

## Recorder use life

The recorder has a maximum life of 170 hours of glucose recording, plus an additional five days of battery life immediately following the glucose recording to allow for data upload. The life span of the recorder begins when it is connected to the sensor. After 170 hours the recorder will stop recording and no further data will be gathered. When the battery dies, any data not uploaded from the recorder will be lost.

## Specifications

Biocompatibility	Recorder: Complies with EN ISO 10993-1
Applied parts	Envision™ Sensor (MMT-7080)
Operating conditions	Temperature: 5 °C to 45 °C (41 °F to 113 °F) Relative humidity: 10% to 95% with no condensation Pressure: 57.6 kPa to 106.0 kPa (8.4 psi to 15.4 psi)
Shipping conditions	Temperature: -30 °C to 55 °C (-22 °F to 131 °F) Relative humidity: 10% to 95% with no condensation Pressure: 57.6 kPa to 106.0 kPa (8.4 psi to 15.4 psi)
Storage conditions	Temperature: 15 °C to 30 °C (59 °F to 86 °F)
Recorder communication frequency	Bluetooth™* version 4.0 (2.4 GHz band)
Modulation	G1D
Maximum output power	-11.5 dBm effective radiated power (ERP)
Operating range	Up to 2.4 meters (8 feet)



## Recorder wireless communication

### Quality of service

The recorder and mobile app connect via Bluetooth™\* Low Energy (BLE). The recorder sends data and related alerts to the app. The recorder and the app verify the integrity of received data after wireless transmission. Quality of the connection is in accordance with the Bluetooth™\* Specification v4.0.

### Data security

The recorder is designed to only accept BLE communications from recognized and linked devices. You must program the app to accept information from a specific recorder. Transmitted sensitive data is encrypted to prevent unauthorized receipt or communication.

## Guidance and Manufacturer's declaration

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The recorder is intended for use in the electromagnetic environment specified below. The customer or the user of the recorder should make sure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The recorder must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class B	The recorder is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

### Guidance and Manufacturer's Declaration - Electromagnetic Immunity


**The recorder is intended for use in the electromagnetic environment specified below. The customer or the user of the recorder should make sure that it is used in such an environment.**

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±2 kV, ±4 kV, ±8 kV, ±15 kV Air ±2 kV, ±4 kV, ±6 kV, ±8 kV Contact	±2 kV, ±4 kV, ±8 kV, ±15 kV Air ±2 kV, ±4 kV, ±6 kV, ±8 kV Contact	For use in a typical domestic, commercial, or hospital environment.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable	Requirement does not apply to this battery powered device.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Not applicable	Requirement does not apply to this battery powered device.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle	Not applicable	Requirement does not apply to this battery powered device.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	400 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical domestic, commercial, or hospital environment.

**Note:**  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Guidance and Manufacturer's Declaration - Electromagnetic Immunity

**The recorder is intended for use in the electromagnetic environment specified below. The customer or the user of the recorder should make sure that it is used in such an environment.**

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 V/m 150 kHz to 80 MHz	Not applicable	Not applicable
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m 80 MHz to 6 GHz	<p>Portable and mobile RF communications equipment should be used no closer to any part of the recorder, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the recorder.</p> <p>Refer to the recommended separation distance table for more information.</p> <p><b><math>d = 0.35 \sqrt{P}</math></b> 80 MHz to 800 MHz</p> <p><b><math>d = 0.70 \sqrt{P}</math></b> 800 MHz to 6 GHz</p> <p>Where <math>P</math> is the maximum output power rating of the recorder in watts (W) according to the recorder manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF recorders, as determined by an electromagnetic site survey<sup>3</sup>, should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The recorder is intended for use in the electromagnetic environment specified below. The customer or the user of the recorder should make sure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
<p><b>Note:</b> At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p><b>Note:</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption, and reflection from structures, objects and people.</p> <p><sup>a</sup>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcasts and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF recorders, an electromagnetic site survey should be considered. If the measured field strength in the location in which the recorder is used exceeds the applicable RF compliance level above, the recorder should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the recorder.</p> <p><sup>b</sup>Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the recorder			
The recorder is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. If you have communication issues when attempting to send data from the recorder, try maintaining a separation distance between the recorder and portable or mobile communications equipment as per the following table:			
Rated maximum output power of recorder (W)	Separation distance according to the frequency of recorder (m)		
	150 kHz to 80 MHz Not applicable	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 6.0 GHz $d = 0.70 \sqrt{P}$
0.01	Not applicable	0.035	0.07
0.1	Not applicable	0.11	0.22
1	Not applicable	0.35	0.7
10	Not applicable	1.1	2.2
100	Not applicable	3.5	7
For recorders rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be estimated using the equation applicable to the frequency of the recorder, where $p$ is the maximum output power rating of the recorder in watts (W) according to the recorder manufacturer.			

**Recommended separation distances between portable and mobile  
RF communications equipment and the recorder**









The recorder is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. If you have communication issues when attempting to send data from the recorder, try maintaining a separation distance between the recorder and portable or mobile communications equipment as per the following table:












Rated maximum output power of recorder (W)	Separation distance according to the frequency of recorder (m)		
	150 kHz to 80 MHz Not applicable	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 6.0 GHz $d = 0.70 \sqrt{P}$

**Note:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

**Icon Table**

	Serial number
	Catalogue or model number
(1x)	One recorder per container/package
(5x)	Five recorders per container/package
	Manufacturer
	Refer to instruction manual before every use (appears blue on label)
	Manufactured in
	Non-ionizing electromagnetic radiation (RF communication)
	Configuration or unique version identifier
	Degree of protection against electric shock: Type BF applied part
<b>IP48</b>	Recorder: 4 is the level of protection against solid objects with a diameter above 1 mm. 8 is the level of protection against the effects of continuous immersion in water 2.4 meters (8 feet) immersion for 30 minutes

Icon Table	
	Humidity limitation
	This product conforms to Australia Radio Requirements
<b>C</b> <b>CE</b> 0459	Signifies European technical conformity
<b>EC</b> <b>REP</b>	Authorized representative in the European community
	Do not reuse
	Fragile, handle with care
	Keep dry
	Recycle cardboard, paper, plastic packaging supplies and unwanted written material
	WEEE Initiative: DO NOT THROW IN TRASH. Recycle device according to local disposal requirements
	Magnetic Resonance (MR) unsafe: keep away from magnetic resonance imaging (MRI) equipment
<b>IC</b>	Complies with Industry Canada Radio Communication requirements
	Use by Date
	Storage temperature
	Transit Temperature

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.™\* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

**Africa:**

Medtronic Africa (Pty) Ltd.  
Tel: +27 (0) 11 677 4800

**Albania:**

Net Electronics Albania  
Tel: +355 697070121

**Argentina:**

Corpomedica S.A.  
Tel: +(11) 4 814 1333  
Medtronic Directo 24/7: +0800 333 0752

**Armenia:**

Exiol LLC  
Tel: +374 98 92 00 11 or +374 94 38 38 52

**Australia:**

Medtronic Australasia Pty. Ltd.  
Tel: 1800 668 670

**Azerbaijan:**

Isomed  
Tel: +994 (12) 464 11 30

**Bangladesh:**

Sonargaon Healthcare Pvt Ltd.  
Mobile: (+91)-9903995417  
or (+880)-1714217131

**Belarus:**

Zarga Medica  
Tel: +375 29 625 07 77 or:  
+375 44 733 30 99  
Helpline: +74995830400

**België/Belgique:**

N.V. Medtronic Belgium S.A.  
Tel: 0800-90805

**Bosnia and Herzegovina:**

Novopharm d.o.o. Sarajevo  
Tel: +387 33 476 444  
Helpline: 0800 222 33  
Epsilon Research Intern. d.o.o.  
Tel: +387 51 251 037  
Helpline: 0800 222 33

**Brasil:**

Medtronic Comercial Ltda.  
Tel: +(11) 2182-9200  
Medtronic Directo 24/7:  
+0800 773 9200

**Bulgaria:**

RSR EOOD  
Tel: +359 888993083  
Helpline: +359 884504344

**Canada:**

Medtronic of Canada Ltd.  
Tel: 1-800-284-4416 (toll free/sans frais)

**Chile:**

Medtronic Chile  
Tel: +(9) 66 29 7126  
Medtronic Directo 24/7: +1 230 020 9750  
Medtronic Directo 24/7 (From Santiago):  
+(2) 595 2942

**China:**

Medtronic (Shanghai) Ltd.  
24 Hour Help (Cell): +86 400-820-1981  
24 Hour Help (Land): +86 800-820-1981

**Colombia:**

Medtronic Latin America Inc. Sucursal  
Colombia  
Tel: +(1) 742 7300  
Medtronic Directo 24/7 (Landline):  
+01 800 710 2170  
Medtronic Directo 24/7  
(Cellular): +1 381 4902

**Croatia:**

Mediligo d.o.o.  
Tel: +385 1 6454 295  
Helpline: +385 1 4881144  
Medtronic Adriatic d.o.o.  
Helpline: +385 1 4881120

**Česká republika:**

Medtronic Czechia s.r.o.  
Tel: +420 233 059 111  
Non-Stop Helpline (24/7): +420 233 059 059  
Zákaznický servis (8:00 - 17:00):  
+420 233 059 950

**Denmark:**

Medtronic Danmark A/S  
Tel: +45 32 48 18 00

**Deutschland:**

Medtronic GmbH  
Geschäftsbereich Diabetes  
Telefon: +49 2159 8149-370  
24-Std-Hotline: 0800 6464633

**Eire:**

Accu-Science Ltd.  
Tel: +353 45 433000

**España:**

Medtronic Ibérica S.A.  
Tel: +34 91 625 05 42  
24 horas: +34 900 120 330

**Estonia:**

AB Medical Group Eesti OU  
Tel: +372 6552310  
Helpline: +372 5140694

**Europe:**

Medtronic Europe S.A. Europe, Middle East  
and Africa HQ  
Tel: +41 (0) 21-802-7000

**France:**

Medtronic France S.A.S.  
Tel: +33 (0) 1 55 38 17 00

**Hellas:**

Medtronic Hellas S.A.  
Tel: +30 210677-9099

**Hong Kong:**

Medtronic International Ltd.  
Tel: +852 2919-1300  
To order supplies: +852 2919-1322  
24-hour helpline: +852 2919-6441

**India:**

India Medtronic Pvt. Ltd.  
Tel: (+91)-80-22112245 / 32972359  
Mobile: (+91)-9611633007  
Patient Care Helpline: 1800 209 6777

**Indonesia:**

Medtronic International Ltd.  
Tel: +65 6436 5090 or +65 6436 5000

**Israel:**

Medtronic  
Tel (orders): +9729972440, option 3 +  
option 1  
Tel (product support): +9729972440,  
option 2  
Helpline: (17:00 – 08:00 daily/weekends –  
Israel time): 1-800-611-888

**Italia:**

Medtronic Italia S.p.A.  
Tel: +39 02 24137 261  
Servizio assistenza tecnica:  
N° verde: 800 60 11 22

**Japan:**

Medtronic Japan Co. Ltd.  
Tel: +81-3-6776-0019  
24 Hr. Support Line: 0120-56-32-56

**Kazakhstan:**

Medtronic BV in Kazakhstan  
Tel: +7 727 311 05 80 (Almaty)  
+7 717 224 48 11 (Astana)  
Круглосуточная линия поддержки:  
8 800 080 5001

**Kosovo:**

Yess Pharma  
Tel: +377 44 999 900  
Helpline: +37745888388

**Latin America:**

Medtronic, Inc.  
Tel: 1(305) 500-9328

**Latvija:**

RAL SIA  
Tel: +371 67316372  
Helpline (9am to 6pm): +371 29611419

**Lithuania:**

Monameda UAB  
Tel: +370 68405322  
Helpline: +370 68494254

**Macedonia:**

Alkaloid Kons Doel  
Tel: +389 23204438

**Magyarország:**

Medtronic Hungária Kft.  
Tel: +36 1 889 0688

**Malaysia:**

Medtronic International Ltd.  
Tel: +603 7946 9000



**México:**

Medtronic Servicios S. de R. L. de C.V.  
Tel (México DF): +(11) 029 058  
Tel (Interior): +01 800 000 7867  
Medtronic Directo 24/7 (from México DF):  
+(55) 36 869 787  
Medtronic Directo 24/7:  
+01 800 681 1845

**Middle East and North Africa:**

Regional Office  
Tel: +961-1-370 670

**Montenegro:**

Glosarij d.o.o.  
Tel: +382 20642495

**Nederland, Luxembourg:**

Medtronic B.V.  
Tel: +31 (0) 45-566-8291  
Gratis: 0800-3422338

**New Zealand:**

Medica Pacifica  
Phone: 64 9 414 0318  
Free Phone: 0800 106 100

**Norge:**

Medtronic Norge A/S  
Tel: +47 67 10 32 00

**Philippines:**

Medtronic International Ltd.  
Tel: +65 6436 5090 or +65 6436 5000

**Россия:**

ООО «Медтроник»  
Tel: +7 495 580 73 77  
Круглосуточная линия поддержки:  
8 800 200 76 36

**Polska:**

Medtronic Poland Sp. z o.o.  
Tel: +48 22 465 6934

**Portugal:**

Medtronic Portugal Lda  
Tel: +351 21 7245100

**Puerto Rico:**

Medtronic Puerto Rico  
Tel: 787-753-5270

**Republic of Korea:**

Medtronic Korea, Co., Ltd.  
Tel: +82.2.3404.3600

**Romania:**

Medtronic Romania S.R.L.  
Tel: +40372188017  
Helpline: +40 726677171

**Schweiz:**

Medtronic (Schweiz) AG  
Tel: + 41 (0) 31 868 0160  
24-Stunden-Hotline: 0800 633333

**Serbia:**

Epsilon Research International d.o.o.  
Tel: +381 113115554

Medtronic Serbia D.o.o.  
Helpline: +381 112095900

**Singapore:**

Medtronic International Ltd.  
Tel: +65 6436 5090 or +65 6436 5000

**Slovenija:**

Zaloker & Zaloker d.o.o.  
Tel.: +386 1 542 51 11  
24-urna tehnična pomoč:  
+386 51316560

**Slovenská republika:**

Medtronic Slovakia, s.r.o.  
Tel: +421 26820 6942  
HelpLine: +421 26820 6986

**Sri Lanka:**

Swiss Biogenics Ltd.  
Mobile: (+91)-9003077499  
or (+94)-777256760

**Suomi:**

Medtronic Finland Oy  
Tel: +358 20 7281 200  
Help line: +358 400 100 313

**Sverige:**

Medtronic AB  
Tel: +46 8 568 585 20

**Taiwan:**

Medtronic (Taiwan) Ltd.  
Tel: 02-21836000  
Toll Free: +886-800-005285

**Thailand:**

Medtronic (Thailand) Ltd.  
Tel: +662 232 7400

**Türkiye:**

Medtronic Medikal Teknoloji  
Ticaret Ltd. Şirketi.  
Tel: +90 216 4694330

**Ukraine:**

Med Ek Service TOV  
Tel: +380 50 3311898  
or: +380 50 4344346  
Лінія цілодобової підтримки:  
0 800 508 300

**USA:**

Medtronic Diabetes Global  
Headquarters  
24 Hour HelpLine: +1-800-646-4633  
To order supplies: +1-800-843-6687

**United Kingdom:**

Medtronic Ltd.  
Tel: +44 1923-205167

**Österreich:**

Medtronic Österreich GmbH  
Tel: +43 (0) 1 240 44-0  
24 – Stunden – Hotline: 0820 820 190



# Medtronic

**Medtronic MiniMed**

18000 Devonshire Street  
Northridge, CA 91325 USA

1800 646 4633

+1 818 576 5555

[www.medtronicdiabetes.com](http://www.medtronicdiabetes.com)

EC	REP
----	-----

Medtronic B.V.  
Earl Bakkenstraat 10  
6422 PJ Heerlen  
The Netherlands



Sanmina Corporation Mexico  
Carretera Guadalajara-Chapala Km 15.5 No. 29  
Tlajomulco de Zuniga  
Jalisco, Mexico 45640

## C €0459

REF	MMT-7781
-----	----------

M10542489-021\_2