




印刷：连线哑油，双面四色

<h3>Troubleshooting</h3>																																			
<table> <tr> <th>Problem</th><th>Solution</th></tr> <tr> <td>No response after startup</td><td>Please use the incidental charging cable to charge. The charging time is 45mins</td></tr> <tr> <td>The app exits quickly</td><td>Delete the existing app and scan QR to reinstall it</td></tr> <tr> <td>Feel tricky with replacing ear tips</td><td>Align the replacement ear tip to camera and tighten clockwise</td></tr> <tr> <td>Indicator light flashes</td><td>Please use the incidental charging cable to charge. The charging time is 45mins</td></tr> <tr> <td>Abnormal display in app</td><td>Interrupt the connection, restart the product to connect again</td></tr> </table>	Problem	Solution	No response after startup	Please use the incidental charging cable to charge. The charging time is 45mins	The app exits quickly	Delete the existing app and scan QR to reinstall it	Feel tricky with replacing ear tips	Align the replacement ear tip to camera and tighten clockwise	Indicator light flashes	Please use the incidental charging cable to charge. The charging time is 45mins	Abnormal display in app	Interrupt the connection, restart the product to connect again																							
Problem	Solution																																		
No response after startup	Please use the incidental charging cable to charge. The charging time is 45mins																																		
The app exits quickly	Delete the existing app and scan QR to reinstall it																																		
Feel tricky with replacing ear tips	Align the replacement ear tip to camera and tighten clockwise																																		
Indicator light flashes	Please use the incidental charging cable to charge. The charging time is 45mins																																		
Abnormal display in app	Interrupt the connection, restart the product to connect again																																		
<p>If you still have problems with the above procedures, please contact our Customer Service Representatives directly for further assistance</p>																																			
<h3>Hazardous substance</h3> <table> <tr> <th rowspan="2">Part name</th><th colspan="6">Hazardous material</th></tr> <tr> <th>(Pb)</th><th>(Hg)</th><th>(Cd)</th><th>(Cr(VI))</th><th>(PBB)</th><th>(PBDE)</th></tr> <tr> <td>Printed Circuit Components</td><td>x</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td></tr> <tr> <td>Case</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td></tr> <tr> <td>Plastic rubber accessory</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td></tr> </table>		Part name	Hazardous material						(Pb)	(Hg)	(Cd)	(Cr(VI))	(PBB)	(PBDE)	Printed Circuit Components	x	○	○	○	○	○	Case	○	○	○	○	○	○	Plastic rubber accessory	○	○	○	○	○	○
Part name	Hazardous material																																		
	(Pb)	(Hg)	(Cd)	(Cr(VI))	(PBB)	(PBDE)																													
Printed Circuit Components	x	○	○	○	○	○																													
Case	○	○	○	○	○	○																													
Plastic rubber accessory	○	○	○	○	○	○																													
<p>This form is formulated in accordance with SJ/T 11364</p> <p>○ : represent the hazardous material ratio is below GB/T26572-2011</p> <p>x : represent the hazardous material ratio is over GB/T26572-2011</p>																																			
<h3>Product parameters</h3> <table> <tr> <td>Brand</td><td>Xlife</td></tr> <tr> <td>Product Name</td><td>Smart Visual Ear-Clean Rod</td></tr> <tr> <td>Model</td><td>Q3</td></tr> <tr> <td>Pixels</td><td>3 Mega-pixels</td></tr> <tr> <td>Lens diameter</td><td>3.5mm</td></tr> <tr> <td>Image transmission rate</td><td>20fps</td></tr> <tr> <td>Battery</td><td>130mAh</td></tr> <tr> <td>Battery life</td><td>45 minutes</td></tr> <tr> <td>Charging time</td><td>45 minutes</td></tr> <tr> <td>Input current</td><td>5V — 0.5A</td></tr> <tr> <td>Gyroscope</td><td>3-axis</td></tr> <tr> <td>Network standard</td><td>IEEE 802.11b/g/n</td></tr> <tr> <td>Working temperature</td><td>-10- 45°C</td></tr> <tr> <td>Rod weight</td><td>25g</td></tr> <tr> <td>Image sensor</td><td>CMOS</td></tr> </table> <div>    </div>		Brand	Xlife	Product Name	Smart Visual Ear-Clean Rod	Model	Q3	Pixels	3 Mega-pixels	Lens diameter	3.5mm	Image transmission rate	20fps	Battery	130mAh	Battery life	45 minutes	Charging time	45 minutes	Input current	5V — 0.5A	Gyroscope	3-axis	Network standard	IEEE 802.11b/g/n	Working temperature	-10- 45°C	Rod weight	25g	Image sensor	CMOS				
Brand	Xlife																																		
Product Name	Smart Visual Ear-Clean Rod																																		
Model	Q3																																		
Pixels	3 Mega-pixels																																		
Lens diameter	3.5mm																																		
Image transmission rate	20fps																																		
Battery	130mAh																																		
Battery life	45 minutes																																		
Charging time	45 minutes																																		
Input current	5V — 0.5A																																		
Gyroscope	3-axis																																		
Network standard	IEEE 802.11b/g/n																																		
Working temperature	-10- 45°C																																		
Rod weight	25g																																		
Image sensor	CMOS																																		
<h3>Attention</h3> <ol style="list-style-type: none"> When cleaning lens, please wipe it carefully with professional alcohol cotton swabs. Please observe your surroundings before use. Stay at least 2ft away from people, animals, or moving objects when using. Don't use in a moving vehicle or other unstable environments. Please don't move around when using. This product is unsuitable for children under 3 years old. This product is not suitable for children to use alone. Do not place in a humid environment for a long time or soak in liquids. There is a built-in rechargeable lithium-ion battery. Even if you do not use it, you should charge it once a month to avoid long-term battery damage. Use a 5V charger to charge the device. Do not use this product while being charged. Avoid direct exposure to the sun. Do not put heavy objects on the product. It is normal for the temperature to rise slightly when using (up to 95°F). 																																			
<h3>Warranty</h3> <table> <tr> <td colspan="2">User info</td></tr> <tr> <td>User name:</td><td>Telephone number:</td></tr> <tr> <td colspan="2">Address:</td></tr> <tr> <td>Seller:</td><td>Telephone number:</td></tr> </table>		User info		User name:	Telephone number:	Address:		Seller:	Telephone number:																										
User info																																			
User name:	Telephone number:																																		
Address:																																			
Seller:	Telephone number:																																		

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.