Implantable Insulin Pump System
Sistema de bomba de insulina implantable
Le système pompe à insuline implantable
Implantierbare Insulinpumpe
Pompa impiantabile per insulina
Implanteerbaar insulinepompstelsysteem van
Implanterbara insulinpumpssystem

options to be considered for diagnosing and restoring a failed or impaired MIP System.

Physicians must be completely familiar with the function of the MIP System prior to its implantation and use.

Patients must be fully trained on how to use the MIP System, and must demonstrate proficiency in programming their PPC and responding to alarms prior to using the MIP System on their own.

Individual components of the MIP System should NOT be used if damaged prior to or during the implantation procedure. The Pump should not be placed in direct contact with other metal implants.

When refilling the Pump, never push on the syringe plunger. When the Refill Needle is properly seated in the Pump fill port, vacuum in the Pump reservoir will automatically draw in insulin from the syringe. When refilling the Pump, never release the syringe-locking ring while there is a vacuum. The handle of the syringe may snap back forcibly and damage the syringe.

When measuring the Pump stroke volume, always check for leaks at the Pump - Pipette connection. If a leak occurs, the measurement is inaccurate and should be repeated.

ADVERSE EVENTS

Although the safety of the MIP System has been well demonstrated during clinical evaluations spanning more than 10 years, the following adverse events have been reported:

Abdominal Pain	Foreign Body Reaction
Inflammation at Refill Site	Skin Disorder
Infection	Urinary Disorder
Necrosis	Ketoacidosis
Retinal Disorder	Skin Erosion
Abnormal Liver Function	Kidney Disorder
Ileus	Pocket Lymph Edema
Hyperglycemia	Catheter Occlusion/Encapsulation
Hypoglycemia	Anti-Insulin Antibodies*

* A small number of patients may develop elevated titers of anti-insulin antibodies in their bloodstream. Clinical evaluations have shown that elevated anti-insulin antibodies alone do not interfere with the maintenance of glycemic control provided with this treatment.

STERILE DEVICE

The Pump and Infusion Catheter are sterilized by ethylene oxide and provided in tamper evident packaging. Do not use if the package has been opened or damaged. The MMT-4104 100 µl Pipette and MMT-4106 Refill Template are not supplied sterile. They should be steam sterilized prior to their use.

WIRELESS COMMUNICATION NOTICE

This device complies with part 15 of U.S. FCC regulations and European standards for the wireless communication of medical devices.

The MIP System is covered by one or more of the following U.S. patents. Comparable foreign patents also exist and other patent applications are pending:

[U.S.] 4,373,527	[U.S.] 4,395,259	[U.S.] 4,525,165
[U.S.] 4,568,250	[U.S.] 4,569,641	[U.S.] 4,573,994
[U.S.] 4,619,653	[U.S.] 4,636,150	[U.S.] 4,731,051
[U.S.] 4,776,842	[U.S.] 5,167,633	[U.S.] 5,176,644
[U.S.] 5,197,322	[U.S.] 5,217,442	[U.S.] 5,257,971
[U.S.] 5,460,618	[U.S.] 5,466,218	[U.S.] 5,514,103
[U.S.] 5,527,307	[U.S.] 5,559,828	[U.S.] 5,797,733
[U.S.] 5,915,929	[U.S.] 6,283,943	[U.S.] 6,427,088
[U.S.] 6,537,268	[U.S.] 6,562,001	[U.S.] 6,564,105
[U.S.] 6,571,128	[U.S.] 6,577,899	[U.S.] 6,585,644
[U.S.] 6,635,014	[U.S.] 6,648,821	[U.S.] 6,659,948
[U.S.] 6,687,546	[U.S.] 6,694,191	

INDICATIONS AND USAGE

The MIP System is indicated for the long-term, continuous, intraperitoneal administration of exogenous insulin, to maintain glycemic control in patients with Type 1 or insulin-requiring Type 2 diabetes mellitus. In order for the insulin therapy to be efficacious, patients must actively participate in their Pump treatment program.

CONTRAINDICATIONS

- Patients unwilling or unable to:
 - Test their blood glucose (BG) four or more times per day
 - Use their BG results to update Pump programming
 - Switch to subcutaneous insulin injections, if their Pump fails
 - Comply with the guidance of their primary care physician
- Return to their physician every 2-3 months for insulin refills
- Patients who reside or travel above 8,000 ft. (2,400 m.) elevation
- Patients with a history of hypersensitivity reactions to titanium alloy or silicone materials that are used in the MIP System

WARNINGS

The Pump can only be used with the MMT-3150 PPC and either the MMT-4024 or MMT-4027 Infusion Catheter. The Pump can only be refilled or rinsed using the MMT-4102 Refill Needle and MMT-4105 Refill Kit. Use of other needles may damage the Pump valve and septum. Only Aventis U-400 insulin can be used in the MIP System. Use of other insulins could cause damage to the Pump, resulting in inaccurate insulin delivery or Pump failure. Exposure to therapeutic ultrasound may damage the MIP System.

PRECAUTIONS

Patients must always have supplies available for taking subcutaneous insulin injections, in the event of a MIP System malfunction. If a malfunction is suspected, patients should immediately discontinue the use of their MIP System, and contact their primary care physician for instructions and MIP System diagnosis.

If insulin delivery is impaired, the Pump may need to be rinsed and/or the Infusion Catheter flushed to remove insulin precipitates. If impaired insulin delivery is suspected, patients should immediately contact their primary care physician for instructions and MIP System diagnosis.

Physicians should review the Physician's Manual for a full description of pump rinsing, catheter flushing and software reset

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