




Study Hub

Operator's manual

Warnings and precautions

- The device is not waterproof. Do not soak or submerge.
- The Study Hub may overheat at high temperatures. Keep the Study Hub out of direct heat and in a well ventilated environment.
- For individuals with a cardiac pacemaker, implantable defibrillator, medical pump, or other implantable medical electronic device: the device may use wireless signals like Bluetooth, Wi-Fi, and cellular connections that could potentially affect your implantable medical electronic device. The device is FCC certified for use with implantables. Please follow the directions for use for your implant. Keep the Study Hub at least 6 inches from your implantable device. If you still have concerns regarding your implanted device, please contact your healthcare provider.



Welcome to the Verily Study Hub

The Study Hub charges your sensor devices and syncs the data to a secure database over a cellular or Wi-Fi connection.


What's inside

- 1 Study Hub
- 2 Power cord

Set up the Study Hub



Plug the power cord into the round port on the back of the Study Hub.



Plug the power cord into an electrical outlet to provide power to the Study Hub. Outlets near an area with good cellular connection are best (for example, near a window).



The light next to the connection symbol (📶) will flash while establishing a connection. The Study Hub is ready to sync when the light changes to solid white, meaning a good connection has been made. The connection process may take up to 10 minutes depending on the connection quality.

Charge and sync the Study Watch with Study Hub


If you received a Study Watch as part of your study devices, you can charge and sync the Study Watch using the Study Hub.



Plug the Study Watch's charging dock into one of the USB ports, either on the side of the Hub or at the back.



You can put the dock either on top of the Study Hub or to the side.



Place the Study Watch on the charging dock. The time and date on the Study Watch should read the same direction as the “V” on the dock.

Connect the Study Hub using Wi-Fi or Bluetooth

As part of the study, your Study Hub may require the use of a Wi-Fi or Bluetooth connection to sync data. Follow the instructions provided as part of the study for additional instruction.

Study Hub lighted symbols

The Study Hub has lights next to symbols to help you understand what the device is doing:



CONNECTION

Flashing white: The Study Hub is looking for a connection. This may take up to 10 minutes.

Solid white: The Study Hub is connected and ready to sync.

Solid orange: The Study Hub is unable to find a connection. Try moving the Study Hub to another location. If that doesn't resolve the problem, contact support.



SYNC

Flashing white: The Study Hub is syncing and charging the sensor device.

Solid white. Sync is complete.



ERROR

Solid orange: Something's wrong. You can try unplugging the power cord and plugging it back in. If that doesn't resolve the problem, contact support.



WI-FI

Flashing white: The Study Hub is looking for a Wi-Fi connection. This can take a few minutes.

Solid white. The Study Hub is connected and ready to sync.

Solid orange: The Study Hub wasn't able to establish a connection over Wi-Fi. Contact support.




BLUETOOTH

Flashing white: The Study Hub is looking for a Bluetooth connection. This can take a few minutes.

Solid white. The Study Hub is connected via Bluetooth and ready to stream.

Solid orange: The Study Hub wasn't able to establish a Bluetooth connection. Contact support.

Turn lights on and off

After 5 minutes of inactivity, all the lighted symbols will go dark. You can turn the lights on and off yourself. On the top of the device, tap the light button  .

Supplies, contents, and storage

The devices are supplied in a cardboard box with appropriate labeling.

Depending on the configuration required for a study, the packaging may contain one or more devices including accessories and supplemental labeling. The devices should be stored in a climate controlled setting.

FEDERAL COMMUNICATION COMMISSION INTERFERENCE STATEMENT

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC CAUTION:

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

RADIATION EXPOSURE STATEMENT:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 20cm between the radiator & your body.

Note: The country code selection is for non-US model only and is not available to all US model. Per FCC regulation, all WiFi product marketed in US must fixed to US operation channels only.



Verily Life Sciences LLC
269 East Grand Ave
South San Francisco, CA 94080 USA

953-00710-01 Rev A