

说明书 516*430mm/80g双胶纸/4C/风琴折 折叠后尺寸为86*86mm



Keep your sound amplifier clean. Heat, moisture, and foreign substances can result in poor pe








HOW TO WEAR YOUR SOUND AMPLIFIER

- Hold the sound amplifier with your index finger and thumb. - Gently insert the sound amplifier into the ear canal, and then rotate it forward or backward until the product is stable in the ear and will not fall.



PAIRING-FIRST USE

1. Open the charging box, the sound amplifier will automatically enter the bluetooth pairing mode. 2. Open your phone settings, turn on the Bluetooth, search for device HAT607, and then click the circle icon to pair and connect.



Model	HAT607
Output.OSPL90 Max.	≤105dB±3dB
OSPL90 HFA Avg	≤101dB±3dB
FOG50 MAX Gain	≤30±3dB
EQ Input Noise	≤22dB+3dB
Frequency Range (Not narrower than)	220Hz-5500Hz
Working current	≤12mA
Total Harmonic Distortion	≤3%
Bluetooth version	5.0
Operating Voltage	DC 3.7V
play time	4H
The number of times the charging box can be charged	>2 times
Headphone charging time	about 2 hours
Charging time of charging box	about 2 hours
Input voltage	5V/500mA

Notice: The above parameter index testing is on basis of FONIX8000Hearing Aid Test System; Inspection standard: IEC 60118-7 2005, 2CC coupling cavity SB-1

- Use a cleaning brush or soft cloth to clean debris from around the user control and microphone; inspect the the receiver, earbud, or ear tip and wax guard regularly. - Never use water, solvents, cleaning fluids, or oil to clean your sound amplifier. - Charging contacts should always be kept clean and dry. Your hearing professional can provide further information on additional maintenance procedures for your sound amplifier, if needed.

MAINTENANCE AND CLEAN

1. Daily maintenance

It is important to keep your instrument clean and dry. Daily, cleanthe hearing instruments using a soft cloth or tissue. To avoid damage due to humidity or excessive perspiration, please use a drying kit.

2.Care and maintenance

Please put it in where children could not reach, dry and no exposed in the sun. Do not place sound amplifier near high-temperature heating areas such as stoves. Do not leave the sound amplifier in a humid place for a long time. Do not scratch the microphone with sharp objects to avoid damage to the microphone. Do not disassemble or attempt any repairs by yourself.

3.Cleaning

1) Cleaning the eartips Remove ear wax plugged in the earpiece by using a stick and brush. 2) Cleaning hook Using a soft cloth or tissue, wipe off any moisture, sweat, or debris from the surface of the device. Note: Do not use alcohol or other cleaning solvents to clean the hook as this could damage the protective covering.

SERVICE AND REPAIR

1. If for any reason, your sound amplifier and charger do not operate properly, do not attempt to fix it yourself. Not only are you likely to violate any applicable warranties or insurance, you may also cause further damage. 2. If your sound amplifier and charger fail or perform poorly, check the guide on the next page for possible solutions. If problems continue, contact your hearing professional for advice and assistance. Most of the problems can find in your hearing professional's office or clinic.

TROUBLESHOO			IMPORTANT SAFETY INSTRUCTIONS	IMPORTANT NOTICE FOR PROSPECTIVE SOUND AMPLIFIER USERS
Symptom	Possible Causes	Solutions	 A 1.CONTRAINDICATION The device is unsuitable in the following situations: Usage is prohibited for patients with chronic suppurative otitis media and congenital ear atresia. This sound amplifier can only compensate patients with hearing impairment and do not have hearing therapeutic value. 1) Congenital atresia of external ear or deformity of external auditory meatus, a large number of foreign bodies in external acoustic meatus; 2) Deafness is unstable, and the degree of hearing loss often fluctuates; 3) Complaint with headache, dizziness, earache, tinnitus, and other symptoms that are not suitable for adaptation; 4) Sudden deafness in the past three months; progressive rapid deafness; apparent unilateral deafness or other acute ear diseases; 5) Effusion in the middle ear and frequent purulence; 6) Central deafness; non-organic deafness; 7) Patients with acute otitis externa and tympanitis; 8) Patients with acute and chronic suppurative otitis media (in the period of purulent infection); 9) Patients who are allergic to the material of this product. 	medically treatable conditions that may affect hearing are identified and treated before the sound amplifier is purchased. Warning to sound amplifier dispensers A sound amplifier dispenser should advise a prospective sound amplifier user to consult promptly with a licensed physician (preferably an ear specialist) before dispensing a sound amplifier if the sound amplifier dispenser determines through inquiry, actual observation, or review of any other available information concerning the prospective user, that the prospective user has any of the following conditions: (i) Visible congenital or traumatic deformity of the ear. (ii) History of active drainage from the ear within the previous 90 days. (iii) History of sudden or rapidly progressive hearing loss within the previous 90 days. (iv) Acute or chronic dizziness.
	Blocked tubing/earbud	Clean or replace ear wax guard as needed.		
Not Loud Enough	Hearing change	Contact your hearing professional.		
	Too much dust	Clean both microphone and receiver with brush.		
Inconsistent Performance	Restart required	Put the sound amplifier into charging case,close the top cover to turn off. Then take out the sound amplifier.		
	Blocked tubing/ earbud	Clean or replace ear wax guard as needed.		
Unclear, Distorted	Restart required	Put the sound amplifier into charging case,close the top cover to turn off. Then take out the sound amplifier.		
Performance	Blocked tubing/ earbud	Clean or replace ear wax guard as needed.		
	Defective sound amplifier	Contact your hearing professional.		
	Charge required	Put the sound amplifier into charging case, close the top cover, and wait for charging.		

		,	
Dead sound amplifier	Blocked earbud/tubing	Clean or replace ear wax guard as needed.	 A 2.WARNINGS Keep small swallow-able accessories out of Dropping, immersing in liquid, and accessing sound amplifier. If a sound amplifier or accessories as possible, do not take it out by yourself. Warning not to use in explosive or oxygen- Warning that only the intended person can result in hearing damage. Statement about the special needs of partice Warning about common conditions that courties warning that results from the risk manager storage. Each warning/safety contains: nature of Hat a storage. Please Do not use it at high volume for a log uncomfortable, please stop using it. Please check the product every time before ear canal. It will affect the sound amplifier. Avoid shocking and dropping the product. If the sauna room). Do not throw the scrapped battery as standed the share with someone else will result in other share with someone else will result in other share with any unusual skin reactions.
ampimer	Ear canal fulfilled with oil wax	Contact your hearing professional.	
No LED indicator when sound amplifier are inserted into the charging port	Low battery	Connect the Type-c cable to your charging case and plug it into the wall. The on-board battery LEDs will flash for a few seconds to indicate the connection to the power source.	
Orange blinking LED	Low battery	Insert your type-c cord to charge the case	
sound amplifier whistling in the charging case	Dead battery	Put the sound amplifier into charging case, close the top cover, and insert your type-c cord to charge the case.	
If the sound amplifier	er Sound Amplifier earplug and ear canal sealing or fit is not good	1. Gently turn the sound amplifier, so that the sound amplifier eartips and ear canal sealing and fitting better, so that you can deal with the whistle phenomenon.	
whistles in the ear		2. The size of the eartips on the Sound Amplifier does not match your ear canal.Use earplug of a different size to improve the whistle phenomenon.	
Sound Amplifier	Low battery	Put the sound amplifier into the charging case to charging	
has no sound	Blocked up the Eartips' sound port or microphone port	Clean the sound port or microphone port	
	Low battery	Put the sound amplifier into the charging case to charging	the sound amplifier, and check with a medic - Do not place the sound amplifier near the m
The voice is very weak and unclear	The volume is too low	Click to adjust the volume level	 Do not place the sound amplifier hear the r Be sure to take off the sound amplifier whe Never pry the inside of the sound amplifier device. Do not put the sound amplifier near heating high moisture areas such as bathrooms or In a closed car in summer, the temperature
	Eartips loose	Re-wearing the sound amplifier	
	Sound amplifier	failure self-checking table	

Factory Reset Mode

1. Put the left sound amplifier and the right sound amplifier into the charging case bin, open the top cover of the charging case, long press the button at the bottom of the charging case for 5 seconds, and the Yellow indicator light on the charging case starts to blink quickly that enter TWS pairing mode or reset mode.

2. The Yellow indicator light on the charging case flashes quickly about 15 seconds and off; 3. You can connect your smart phone by "Pairing First Use" page;



Press and hold for 5 seconds

HELPFUL HINTS

- Do not take apart your sound amplifier or insert the cleaning tools inside them. (When not in use, place your sound amplifier in the charger and store them.)
- In a dry, safe place.
- Away from direct sunlight or heat to avoid extreme temperatures.
- Where you can easily find them.
- Safely out of reach from children and pets. Synergy Charger Care
- Keep your charger clean. Heat, moisture, and foreign materials can result in poor performance.
- Do not use water, solvents, or cleaning fluids to clean the charging ports.
- Keep lid closed as much as possible to avoid dust and debris build-up.
- Store your charger in a clean and dry location. - Fully charge sound amplifier batteries every night.
- INTENDED USE

PDF

- 1. Medical purpose - The sound amplifier intends to compensate for impaired hearing.
- 2. Intended patients and environment
- The sound amplifier uses for indoor and outdoor environment, and patients under 18 years old must equire the doctor's approval.

of children's reach. sing strong electromagnetic fields or over heat will damage essories remain in the ear, please go to see a doctor as soon -enriched atmospheres; it may damage the device. n use it; using a sound amplifier designed for someone else may ticular patient groups. ould damage the device. ment file permissible environmental conditions of transport and lazard; likely consequences; precautions

long time. Otherwise, it will injure your ear. If you feel re using it, and avoid wearing the earplug with foreign objects in Do not keep or use it in high humidity environments (such as dard rubbish. osive or dangerous environment. ers. Your sound amplifier setting is personal and individual, ers' hearing damage. . In some instances, wearing a sound amplifier can cause an

shes. If you notice any uncommon skin condition, stop wearing ical specialist.

magnetic objects. en getting a diagnosis or treatment.

with a needle or similar object; otherwise, it will damage the

g appliances or in direct sunlight, and you should also avoid near a humidifier.

- In a closed car in summer, the temperature can rise to very high levels. Do not leave your sound amplifier in a high-temperature car.

- Do not use the sound amplifier close to a cell phone, otherwise it may increase the noise and affect the sound quality.

- Please contact the MANUFACTURER:
- Any questions in setting up, using, or maintaining the sound amplifier;
- To report unexpected operations or events.

MATTERS NEEDING ATTENTION, WARNINGS, AND TIPS

1. Do not immerse the sound amplifier in any liquid. Remove the sound amplifier before swimming and shampooing. After swimming or washing, wait for the ear canal to dry before wearing the sound amplifier. 2. Please use the function and volume control switch correctly, open and close as lightly as possible to

- avoid excessive force.
- 3. Applicable objects: Suitable for people with air-conducting hearing loss of about 30-50dB. 4. Please use sound amplifier with sound amplifier charging cases.
- 5. When you do not use sound amplifier, please charge the product regularly to avoid battery damage; the recommended charging time interval of sound amplifier is ≤28 days;

MAINTENANCE, STORAGE AND DISPOSAL

1. Maintenance

- 1) We do not authorize any institution or individual to maintain and repair the product. If you suspect that
- the products have any issues, please contact the manufacturer or distributor to handle them. 2) The user must not attempt any repairs to the device or any of its accessories. Please contact the retailer for repair.
- 3) Opening of the equipment by unauthorized agencies is not allowed and will terminate any warranty claim. warning: No modification of this equipment is permitted!

2.Storage

- 1) If you Do not use it, and then put the sound amplifier into the charging case.
- 2) Do not put the device under the sunshine, high temperature, and moist environment or someplace,
- which maybe get in touch with fire or is vulnerable to vibration
- 3) Keep the device storage in following environmental ranges:-20~55 C, ≤95%RH, 700-1060hPa

3. Disposal

Please dispose of the device under the directive 2002/96/EC – WEEE (Waste Electrical and Electronic Equipment).and Don't dispose with domestic waste. If you have any queries, please (vii) Visible evidence of significant cerumen accumulation or a foreign body in the ear canal (viii) Pain or discomfort in the ear.

Special care should be exercised in selecting and fitting a sound amplifier whose maximum sound pressure level exceeds 132 decibels because there may be risk of impairing the remaining hearing of the sound amplifier user. (This provision is required only for those sound amplifier with a maximum sound pressure capability greater than 132 decibels (dB).)

FCC WARNING

RF exposure statement

This equipment complies with the FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

FCC Warning

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE 1: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures: - Reorient or relocate the receiving antenna.

- Increase the separation between the equipment and receiver.

-Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. -Consult the dealer or an experienced radio/TV technician for help.

NOTE 2: Any changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

WARM TIPS:

Due to the different actual situations, each customer's experience of using product is different. If you have any questions in the using process, please contact our professional customer service staff by our aftersales email (aftersalesmailbox@163.com) or by the platform where you purchased this product as soon as possible. Please believe that our customer service staff will offer you satisfactory assistance.

NORMALIZED SYMBOLS				
Ĩ	Consult instructions for use			
Ŕ	B type applied part			
F©	FCC ID: 2AGA6-HAT607			
	Disposal Please dispose of the device under the directive 2002/96/EC – WEEE (Waste Electrical and Electronic Equipment). and Don't dispose with domestic waste.			
	Manufacturer			
~	Date of manufacture			
SN	Serial number			
	Warning/ Danger: Improper use might cause serious injuries.			
IP41	IP code of the device: this device's grade of against ingress of solid foreign objects – the access probe of 2.5 mm \emptyset shall not penetrate; the grade of waterproof is vertically falling drops shall have no harmful effects.			

MANUFACTURE: JIANGXI YIANG MEDICAL TECHNOLOGY CO., LTD

Address: Floor 3, Building 3, Ganzhou National High-level Talents Science And Technology Innovation

3.Side effects - For patients with sensitive skin, the silicone material may irritate your ear canals. - First use may cause slight pain to the ear canal due to maladjustment.

refer to the local authorities responsible for waste dis-postal.

Park,No.1 Wudang M Ganzhou,Jiangxi,341000,CHINA