

Specification

Model: JB-318
Rated voltage: 5V-1A
Rated power: 5W
Material: ABS
Reliable, safe, comfortable and ergonomic design
Standard: GB4706.1-2005 GB4706.10-2008

Trouble and Solution

- Can't switch on the power source?
- low battery,
- please arrange charging.
 - incorrect connection for power adaptor,
- please connect it properly.
 - broken inner lithium battery,
- if can't charge, please change new battery.



Eye Massager

Instruction Manual



Improve blood circulation



Relieve eye fatigue

Use 10-15, every day

List of Contents

Thanks for purchasing this eye massager. Please read this instruction manual thoroughly before using, especially the precaution to ensure correct usage. And please keep this manual well for future use.

Brief Introduction

Product Feature	1
Component List	2
Preparation	2
Operation Method	3
Precaution	4
Specification	5

Product functions are subject to upgrade or change without notice, please refer to the actual product.

Brief Introduction

Product Principle

Adopting far infrared heating and microcomputer control technology through hot pack, vibration and relaxation to eyes, this eye massager helps eyes to increase oxygen and dredge main and collateral channels. It enhances crystalline lens regulation ability, restores the elasticity of ciliary muscle, improves eye blood circulation, alleviates visual pressure, relieves eye fatigue, reduces visual impairment, thus maintains eyes' health. This product setup includes relaxing music, brainwaves and nature music.

Application Area

It has an auxiliary massage effect on pseudomyopia, shallow myopia and amblyopia. It also has the health care function of alleviating eyestrain and preventing myopia. It is suitable for college, tertiary and primary school students, myopia and hyperopia people, eye overuse person, computer operator, accountant, designer, secretary, teacher, long time eye using worker, neuroticism, insomnia and old persons with uncomfortable eyes.

Do SPA for eyes at any time

Good helper for eye care



Right Person

Glasses Wearer

Computer User

Product Feature

Using the computer chip control technology to carry on the humanity simulation control, this product is a compound high-tech eye massager. The product adopts advanced high-tech infrared heating body which has the characteristics of uniform heating, low power consumption, high far-infrared density and light weight. Infrared ray can activate the irradiated cells, making the molecules produce resonance, enhance the intermolecular force, enhance the regeneration ability of tissues, improve immunity, increase oxygen and dredge collaterals and achieve the health effect.

Acupoint Function

This product successfully applies the new concept of air pressure to eye health care, making the product more humanized. Under the control of micro-computer, the air pressure can knead the eyes in an all-round way, relieve the pressure on the eyes, fully relieve the eye fatigue, restore the contraction ability of the ciliary muscle, regulate the physiological elasticity of the ciliary muscle and relieve the elongation or atrophy of the eye axis.

This product uses the built-in three-dimensional magnetic field and the moving when cutting the magnetic field line of the main acupoints of the face and the periorbital capillaries, so as to generate micro electric fluid inside the capillaries to achieve the function of activating meridians and collaterals and harmonizing Qi and blood.



Component List

1. Massager body
2. USB data cable (for charging)
3. Colour box
4. Operation manual
5. Qualified certificate

Product Operation

Preparation

1. Turn on the machine to check whether the power of the massager is sufficient
2. If the power is insufficient, please use DC5V stable power supply configured by our company to charge until the power is sufficient

Charge

1. Find the charging port on the left side of the massager body
2. Insert the output plug of the USB charging line into the charging socket of the product, and then connect the USB end to the charging adaptor of the home mobile phone.
3. During the charging process, the red charging indicator is always on, and the red indicator is off after the charging is completed (the blue indicator is on after charging is completed for voice model).
4. During the charging process, then machine can not work (self-protection function, automatic shut off for the machine during charging).
5. Charge for 3 hours.



Slow life, keep going

Key Layout



Operation Method

Press the switch for 2 seconds to turn on the machine. The default mode is automatic. When press key shortly, the key will enter the next mode. The massage mode will be cycled in turn. Press the button for 2 seconds to switch off the machine.

Mode cycle with one button control:



The default working time of the machine is 15 minutes after it switches on. If you want to extend the using time, you need to switch on the machine again(it is recommended that the using time should not exceed 15 minutes for each time)

Music control key+

Press this key for a long time to control the volume and sound to increase. Press this key for a short time to switch to the next music.

Music control key-

Press and hold this key for a long time to control the volume and sound to be smaller. Press this key for a short time to switch into the previous music.

Bluetooth connect

Switch on the bluetooth of the mobile phone, search for bluetooth devices, and connect bluetooth in pairs.

Call answer key

Connect bluetooth, press this key when call comes to answer.

Precaution

Please read thoroughly and understand all contents of the manual. Please operate according to the operation manual

1. It is forbidden to use for patients with eye diseases, such as eye surgery, glaucoma, cataract, retinal detachment and etc.

2. For patients with deep stroke, high blood pressure and hypotension, please consult a professional doctor. If you feel uncomfortable during use, please stop using immediately.

3. Those who are seriously sluggish in actions, have physical and mental discomfort or are mentally handicapped should use the machine with caution. Please use it under the supervision of the caregiver. If there is any discomfort during the use, please stop using immediately.

4. Please avoid using eye massagers after drinking alcohol.

5. Please make sure there is no contact lens in the eyes when using the eye massager.

6. Do not use this product in the environment with temperature over 40 °C and excessive humidity.

7. It is forbidden to wash the product with water. Please use cotton cloth soaked with alcohol for cleaning.

8. Do not disassemble or attempt to repair by yourself.

FCC Caution.

(1) § 15.19 Labelling requirements.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

§ 15.21 Changes or modification warning

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

§ 15.103 Information to the user:

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the

instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of

the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

RF warning for Portable device:

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

According to § 15.247(e)(1) and § 1.1307(b)(1), systems operating under the provisions of this section shall be operated in a manner that ensures that the public is not exposed to radio frequency energy level in excess of the Commission's guidelines.

According to KDB 447498 (2)(a)(i)