

封面logo尺寸：36*9.3mm或30*7.7mm



1. Scale is a high-precision measuring device. Never jump or stamp on the scale or disassemble. Please handle the scale carefully to avoid breakage.
2. Place the scale on a hard, flat surface during use. Using the scale on a soft surface like carpet or linoleum may affect the accuracy of the scale readings.
3. To ensure weighing accuracy, please gently step on the scale to wake it up and wait at least 10 seconds before weighing yourself on the scale. You must also wake up the scale if it has been moved. Take your measurements at the same time each day to ensure the most accurate results.

The Bluetooth Scale can connect with your mobile phone (Android 6.0+ & IOS 8.0+) or iPad via Bluetooth 4.0. Download the 'Aidot' app free from the App Store or Google Play.

Conductive areas

Click [+]

Click [Add Device]

Select [BLE Body Fat Scale]

A diagram illustrating the conversion of a bracketed expression to a decimal value. On the left, there is a bracketed expression $\left[\frac{1}{100} \right]$. An arrow points to the right, where the result is shown as 0.01 kg .

- Ensure personal data is entered.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

A Better Life Is Worth Living