

# A170 智能电子药盒使用说明书

A170 Smart Electronic Pill Box Instruction Manual

使用本产品前请仔细阅读说明书，并妥善保管。

Please read the instruction manual before using this product, and hold it for safekeeping.

## 目录

### Contents

应用下载方式 Downloading the App

包装内容 Package Contents

产品介绍 Product Description

显示介绍 UI Description

按键说明 Buttons Description

指示灯说明 Indicator Light Description

操作视频 Tutorials

产品使用 Instructions

产品规格 Specifications

产品清洁维护 Care and Maintenance

注意事项 Attentions

## 应用下载方式 Downloading the App

方法一：请扫码下载 App

Method 1: Scan the QR code to download the app

方法二：搜索下载

Method 2: Search and download

## 包装内容 Package Content

药盒* 1	Pill Box
说明书*1	Instruction Manual
数据线*1	Data Cable

注意：如果电池电量低，请充电后再使用。  
Note: If the battery power is low, please charge it before use.

产品介绍 Product Description

- 外盖 Outer Cover
- 内盖 Inner Cover
- 按键 Button
- 显示屏 Display
- 充电口 Charging Port
- 指示灯 Indicator Light

显示介绍 UI Description

- 电量图标 Power Level
- 设置的闹钟个数 The number of alarm clocks set
- 音量 Volume
- 时钟时间 Time
- 蓝牙 Bluetooth

按键说明 Buttons Description

状态 Status	操作 Operation	功能 Function
息屏时 Screen off	单击 Press once	唤醒屏幕，可连接蓝牙 Turn on the screen and Bluetooth
亮屏时 Screen on	单击 Press once	熄灭屏幕，蓝牙将断开 Turn off the screen and Bluetooth
闹钟响起时 Alarm goes off	单击 Press once	取消当前闹钟 Turn off the current alarm

指示灯说明 Indicator Light Description

蓝灯闪烁 Blue light flashes	闹钟响起 Alarm goes off
蓝灯常亮 Blue light always on	外盖打开 Outer cover opened
蓝灯灭 Blue light off	合上外盖 Outer cover closed

操作视频 Tutorials

请扫码观看操作视频  
Please scan the QR code to access tutorial videos.

## 产品使用 Instructions

### 1.补充药物 1. Replenish medicine

- ①打开外盖 Open the outer cover
- ②打开内盖 Open the inner cover
- ③填充药物 Replenish medicine
- ④合上盖子 Close the covers

### 2.药盒连接 APP Connect the pill box to the app

- ①唤醒屏幕 Turn on the screen
- ②iOS: 在 iOS 控制中心打开蓝牙  
iOS: Turn on Bluetooth in Settings  
    安卓: 在手机设置中打开蓝牙和定位  
    Android: Turn on Location and Bluetooth in the phone settings
- ③通过蓝牙连接药盒并在 APP 中设置提醒  
Connect the app to the pill box through Bluetooth, and set your alarms

### 3.取药 Take medicine

- ①闹钟响 Alarm goes off
- ②打开药盒 Open the pill box
- ③倒出药物 Pour out the medicine
- ④合上盖子 Close the cover

### 4.取消闹钟

#### Turn off the alarm

- (1)在药盒上取消 Turn off the alarm on the pill box
  - ①闹钟响起时 When the alarm goes off
  - ②单击按键, 取消当前闹钟 Press the button once to turn off the current alarm
- (2)在 APP 上取消 Turn off the alarm on the App
  - ①在蓝牙连接状态下 When the pill box is connected via Bluetooth
  - ②关闭闹钟或者删除闹钟 Turn off or cancel the alarm

### 5.铃声设置

#### Set the ringtone

- ①在蓝牙连接状态下 When the Pill Box is connected via Bluetooth
- ②打开/关闭铃声 Turn the ringtone on or off

## 产品规格 Specifications

名称 Name	智能电子药盒 Smart Electronic Pill Box
------------	-------------------------------------

型号 Model	A170 A170
尺寸 Dimension	218.2X49.3X37.9cm 218.2 x 49.3 x 37.9 cm
重量 Weight	167.60g 167.60 g
通讯接口 Connection	蓝牙 Bluetooth
供电方式 Power Port	Type-C USB-C
电源输入 Power Input	5V=1A 5V=1A
支持操作系统 Supported OS	Android Android
	iOS iOS

## 产品清洁维护 Care and Maintenance

①用柔软潮湿的布擦拭内格。

Use a soft damp cloth to clean the interior of the pill box.

注意：不要弄湿电子件部分，可能会发生短路。

Note: Do not get the electronic components wet. It may cause a short circuit.

②在补充药物前，请确保药盒的清洁干燥。

Before replenishing the medicine, please ensure that the pill box is dry and clean.

③电量低时，请使用配套的数据线充电。

When the battery power is low, please use the data cable supplied to charge it.

## 注意事项 Attentions

①不要让孩子使用或者玩这个设备。

Keep the device away from children.

②确保将药盒放在安全位置，坠落可能会损坏设备或里面的药片。

Make sure to keep the medicine box in a safe place. Dropping may damage the device or the pills inside.

③如果您因任何原因错过剂量，请联系您的医疗保健提供者或按照您的医生或药剂师提供的指示。

If you miss taking your medicine, please contact your healthcare provider or follow the instructions provided by your doctor or pharmacist.

④一些药物应严格密封，不要暴露在空气中。如果您不确定此药盒是否适合存储您的药物，请咨询您的医疗保健提供商或药剂。

Some medicine should be kept sealed. Do not expose them to the air. If you are unsure whether the pill box is suitable for storing your medicine, please contact your healthcare provider or pharmacist.

⑤当电量图标只剩一格时，请及时充电，以确保药盒的正常运行。

When there is only one bar of battery power left, please charge it in time to ensure the normal operation of the pill box .

⑥不要将设备放在明火或热源附近，火灾、触电、产品损坏或伤害可能发生。

Do not place the device near naked flames or heat sources to prevent fires, electrical shocks, damages to the product, and any other complications.

### **FCC Caution.**

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### **RF warning for Mobile device:**

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 0cm between the radiator & your body.