



125mm

100mm

FCC Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

USER MANUAL

Electronic scale CF516BLE

Thank you for choosing Electronic scale.

This model utilizes BIA (Bio-electrical Impedance Analysis) technology. Bio-electrical impedance technology makes the measurements more accurate by taking advantage of the conductive composition of muscle and water contents within the body. The device sends harmless signals through the body to measure the overall body-impedance which is increased due to fat composition and lowered by the conductive muscles and water in the body.

SPECIFICATION

Size: 280x280x25mm	Division: 0.05kg/0.1lb
Capacity: 0.1~180kg	Unit: kg /lb
Batteries: 3*1.5 V AAA	

ACCURACY

0.1-2Kg±0.00Kg;	2-5Kg±0.05Kg	5-50Kg ±0.1Kg;
50-100Kg ±0.2Kg;	100-150Kg ±0.3Kg;	150-182.5Kg ±0.5Kg

CAUTION

- This scale should not be used to diagnose or treat any medical condition. You should consult a your primary care doctor or physician before undertaking your diet or exercise program.
- DO NOT stand on the edge of the scale or jump on it.
- DO NOT overload the scale (Maximum 180kg).
- DO NOT immerse the scale in water.
- Please use a damp cloth or a glass cleaner when you clear the surface of scale.
- Always place scale on a hard, dry and flat surface before measurement.
- Connect the batteries in accordance with the correct polarities indicated.
- If you are not using the scale, it is recommended that you remove the batteries.

- Not recommended for any person using a pacemaker or other internal device.
- Not recommended for women who are pregnant, the measurements will be incorrect.
- Make sure your feet are bare and dry when you measure body fat and other data.

Frequently Asked Question & Answer:

Why doesn't the scale work? Why doesn't the display show data?
Why does the data on the screen disappear in a flash?
- Please check whether the insulation piece is taken out, batteries are installed or batteries are in low and need replacing.

Display indicates "Err" when measuring body fat:
- Your feet or scale are too wet.
- Overload the maximum capacity (180kg)
- Body parameter settings should be with Height (3'3"-7'3"/100-220cm), Aged (10-99 years old).