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EN

Smart digital bath scale:

- Shows body weight on the scale display
- 13 body composition metrics are available in the **Fitdays™** app.

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INTRODUCTION

Dear customer, Laica would like to thank you for choosing one of our products, designed according to strict performance and quality criteria to ensure total satisfaction.

TABLE OF CONTENTS

Introduction	1
Important note for users	1
Product specifications	1
For your safety	2
Safe use of batteries	2
Recommendations for getting an accurate reading	3
Preparations before use	3
Operations instructions	3
First time connecting the scale to your smartphone	4
After having connected the scale to smartphone	4
Body composition guide	5
Maintenance and troubleshooting	7
Warranty	8

IMPORTANT NOTES FOR USERS

READ CAREFULLY BEFORE USE - KEEP IN A SAFE PLACE FOR FUTURE REFERENCE

The instruction manual must be considered an integral part of the product and must be kept for its entire life. If the product is transferred to another owner, its documentation must also be transferred in its entirety. To ensure the safe and correct use of the product, the user must carefully read the instructions and warnings contained in the manual insofar as they provide important information regarding safety, user and maintenance instructions. If the instruction manual is lost or if you require additional information, please contact the company at the address shown below:


SPRINGFIELD INSTRUMENTS INC.
104 Barr Street, Saint Laurent, Quebec H4T 1Y4
1-844-687-4040
info@springfieldinstruments.com

This is a digital scale that measures weight and calculates body composition based on height, age and gender. Data from the scale is sent to your smartphone so you can monitor your progress at any time.


PRODUCT SPECIFICATIONS

Fig. 1

Scale - Platform side

Scale - underside

PRODUCT DESCRIPTION

1. Display
2. Electrodes
3. Battery compartment
4. Key to select unit of measurement (lb/kg)

TECHNICAL CHARACTERISTICS

- Capacity: 180 kg / 400 lb
- Division : 100 g / 0.2 lb
- 2 disposable 1.5V AAA batteries
- Operating frequency: 2426-2460 MHz, effective radiated power: <0.9 dBm
- Environmental conditions for use: +10°C ~ +40°C; RH ≤ 85%
- Environmental storage conditions: +10°C ~ +50°C; RH ≤ 85%

The scale calculates the body composition for individuals aged 12 to 99 measuring between 90 and 220 cm (1ft to 8ft 11in) in height.

FOR YOUR SAFETY

This section explains precautionary measures to take to avoid injury to the users of this device and others and to prevent damage to property.

Please read this information to ensure safe operations of this equipment.

Drawing**SYMBOLS**

Warning

Prohibition

SAFETY WARNINGS

- Before using, make sure the device is intact with no visible damage. If you have any doubts, do not use the device and contact the retailer.
- Keep the plastic bag used for packaging out of reach of children: suffocation hazard.
- This device must only be used for its intended purpose and in the manner indicated in the instructions. All other uses are considered improper and therefore hazardous. The manufacturer cannot be considered responsible for any damage due to improper or incorrect use.
- This device can be used by children on their own from the age of 8 upwards. Individuals with reduced physical, sensory or mental capacities or unskilled people must only use the product if previously instructed on its use, in safety and only if informed of the dangers connected to the product itself. Children must not play with the device.
- Cleaning and maintenance operations must not be performed by children without adult supervision.
- In the event of a failure and/or malfunction, do not tamper with the device. Always contact your retailer for repairs.

- Do NOT allow liquids to penetrate into the body of the scale. This appliance is not waterproof.
- Do NOT leave and Do NOT use the scale in water-saturated environments; humidity rates greater than 85% or contact with water or other liquids may cause corrosion or compromise proper scale operation.
- Do NOT step onto the scale with wet feet, or if the surface of the scale is wet: slipping hazard.

**ATTENTION! BEFORE USING THIS APPLIANCE**

All 13 data calculated by the scale includes: BMI, %body fat, muscle rate, fat-free body weight, %subcutaneous fat, visceral fat, % body water, skeletal muscle, muscle mass, bone mass, protein, BMR and body age must be considered as indicative data only. For further information, consult a doctor or dietician.

- **Self-measurement means checking, not diagnosis or treatment. Unusual readings must always be discussed with your physician. Under no circumstances should you alter the dosages of any medication prescribed by your physician.**
- **Through the electrodes located on the platform, the scale sends a weak, low intensity electrical current that passes through the body. The scale is not dangerous in any way.**
- **The scale is not recommended for the following individuals as the data detected by the scale may not be accurate:**
 - individuals wearing a pacemaker or other electronic medical device
 - people with fever
 - children under 12 years of age
 - pregnant women
 - individuals undergoing dialysis, with oedema in the legs or osteoporosis
 - people who take cardiovascular medication

SAFE USE OF BATTERIES

- Remove the batteries if the scale is not used for long periods of time and keep them in a well-ventilated and dry place, at room temperature.

- ⊗ Do NOT recharge the batteries if they are not rechargeable.
- ⊗ Do NOT recharge rechargeable batteries using methods other than those indicated in the manual or using unsuitable equipment.
- ⊗ NEVER expose the batteries to sources of heat or direct sunlight. Failure to observe these instructions may result in damage to and/or the explosion of the batteries:
- ⊗ Do NOT throw the batteries into an open flame.
 - The removal or replacement of batteries must always be performed by an adult.
 - Keep the batteries far from the reach of children: swallowing batteries carries the risk of death. If batteries are swallowed, immediately contact a doctor.
 - The acid contained in the batteries is corrosive. Avoid contact with the skin, eyes, or garments.

RECOMMENDATIONS FOR GETTING AN ACCURATE READING

Always weigh yourself on the same scale at the same time, preferably when you wake up and before getting dressed.

Wait about 15 minutes after getting up so the water can distribute itself through your body.

It is important to weigh yourself with dry, bare feet positioned correctly on the electrodes (see figure).

Avoid contact between the legs so the scale operates correctly.

Scale precision could vary in the following cases:

- a. high level of alcohol
- b. high level of caffeine or drugs
- c. after heavy meals
- d. after rigorous sports activities
- e. during illness
- f. during pregnancy
- g. when dehydrated
- h. with wet or unclean feet
- i. swollen legs

Figure.

PREPARATIONS BEFORE USE

Inserting the Batteries:

This scale works with 2 x 1.5V AAA batteries.

1. To insert the batteries, open the battery compartment by pulling up the tab on the cover.
2. Insert the batteries into the compartment according to the polarity markings.
3. After inserting the batteries, replace the battery cover.

Selecting the weighing unit of measurement:

Switch on the scale by stepping onto it and checking the unit of measurement indicated on the display: kg, lb and st:lb (UK measurement). To

OPERATIONS INSTRUCTIONS

This scale can be used in 2 ways:

- I. Body Weight Measurement only
- II. Body Weight and Body Composition Measurement with the **Fitdays™** app.

I. Body Weight Measurement only:

(does not require pairing to smartphone application)

Place the scale on a flat, sturdy surface.

Note: An uneven or disjointed floor (e.g., rug, carpet, linoleum) may compromise the accuracy of the reading.

2. With the scale off, step onto it with your feet parallel and your weight evenly distributed. The scale will begin measuring. Do not move until the scale has measured your weight. The weight will flash 3 times.

Note: Do not lean on anything while the scale measures your weight to avoid altering the readings.

Step off the scale. The display will continue to show the reading for a few seconds before automatically turning off.

II. Body Weight and Body Composition Measurement with the **Fitdays™** app.

(requires pairing to smartphone application)

When it is paired successfully, 13 body composition metrics will be available to be viewed on the **Fitdays™** app.

FIRST TIME CONNECTING THE SCALE TO YOUR SMARTPHONE

The scale sends body composition data to the smartphone and make it available on the **Fitdays™** app.

During the following operation, we recommend keeping the smartphone near the scale to keep the connection active.

List of compatible devices:

- IOS 8.0 and above
- Android™ version 4.3 (Jelly Bean) and above

1. Turn on bluetooth on your smartphone settings.
2. Download the free **Fitdays™** app from App Store or Play Store™.
You can also scan the QR code with your smartphone's camera to direct you to the app in the Apps store.
3. Launch the **Fitdays™** app.
4. Follow the instructions in the app to pair the LAICA scale and to create a new user.
If the scale is paired successfully with the app, you will see the weight measuring on the app.

If the weight is not measuring on the app, please make sure the device is paired by trying the following:
In the app, go to **Account** → **Device**. You should see an entry of "MY_SCALE". If not, tap on the + to add a device. Step on the scale and tap on **searching the device for pairing**. If the scale is found, it will appear with an option to connect. Tap on "connect" and it will confirm with "Device is paired". Tap back to the **Measure** screen to begin measuring.

AFTER HAVING CONNECTED THE SCALE TO SMARTPHONE

1. Ensure that you are on the **MEASURE** screen:
2. The scale will begin measuring. If the scale is paired successfully with the app, you will see the weight measuring on the app.

3. Stand still and do not move until the scale has measured your weight.
It is finished when the weight flashes 3 times.

At this time, the app will show:


- Your weight
- Target weight (if set in the app)

At-a-glance changes

Fluctuation in the following 3 health metrics will appear:

(You can set up to compare them to the most recent data or to a specified date where data is available.

- Weight
- BMI
- Body fat %

13 body composition metrics in addition to your weight are also available on the "Details" section on the app (click on "more" to expand and scroll down to see all the metrics). **The 13 body composition metrics are: BMI, % body fat, muscle rate, fat-free body weight, %subcutaneous fat, visceral fat, % body water, skeletal muscle, muscle mass, bone mass, protein, BMR and body age.**

- When you weigh yourself, make sure that the app is active and displayed on your smartphone so that you can check that the data has been promptly sent with your body composition. Before taking a new reading, always make sure that the connection on your smartphone is active.
- For Android™ devices it is necessary to activate the localisation services.
- Before weighing yourself for the first time, always select the user on your smartphone.
- Each time the scale is moved or lifted and then placed back on its support surface, the first weighing must not be considered due to the fact that the electronic board needs to readjust with respect to the load cells. Consider the second and successive weights, as in certain cases the first one could be altered for the reason mentioned above.

About Fitdays™ app:

Fitdays™ is the effort of ICOMON's application development team which has more than 30 years experience of the weighing measurement industry. The app has been upgraded iteratively for more than 5 years. Fitdays™ provides fast and accurate insights into body compositions and helps with management for a healthier life.

Here are some of Fitdays™ app features:**• Core Data**

For analyzing body compositions

• Girth Chart Record

Helping with better bodybuilding

• Trend Chart

Tracking changes of body compositions

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After getting your authorization, Fitdays™ could synchronize your body measurement data with these fitness applications.

BODY COMPOSITION GUIDE

The body composition is made up of:

- lean mass: all components without fat (water, muscle, connective tissue, and other non-fatty organic components), bone, bowels, organs, and blood;
- fatty mass: all fats that can be extracted from adipose tissue and other tissues/organs;
- bone mass.

Through the electrodes located on the platform, the scale sends a weak, low intensity electrical current that passes through the body. The different masses respond differently to the current: lean mass provides less resistance than the opposite. This body resistance is called bioelectrical impedance. The relative percentage of the various body masses varies

Below are some common body composition metrics that appear in app when the scale is connected to your smartphone (see section connecting the scale to your smartphone). The information to follow is to give you an indication only. Consult your doctor or dietician if you have any questions or concerns.

Considerations

When assessing your body composition trend, it is necessary to keep the following in mind:

- A decrease in weight accompanied by no change in body fat or even an increase in body fat means that your body only lost water (for example, after exercising, time in the sauna or a limited diet for rapid weight loss).
- An increase in weight accompanied by a decrease or no change in body fat indicates the development of muscle mass.
- When the weight and body fat decrease at the same time it means that your diet is working and volume is being lost.

Body Fat Percentage

Body fat percentage is the amount of body fat as a proportion of your body weight. The body fat percentage indicated in the table below represents the total fat that is calculated by adding essential fat and storage fat.

Essential fat is necessary for various bodily functions. The percentage is 3% in men and 12% in women. Storage fat gives energy protects against physical trauma and insulates the body to preserve heat.

These percentages vary with the age and gender of the individual as shown in the table below. In Fitdays™, there will be an automatic comparison of your

Bodyfat Table					
	Thim	Athietes	Fitness	Acceptable	Too high
Women	<14%	14~21%	21~25%	25~32%	>=32%
Men	<6%	6~13%	13~17%	17~25%	>=25%

*The data in the table is for references only.

Percentage of Body Water

The % of water in the body differs based on gender: in healthy adult women it ranges from 45% to 60% while in men it ranges from 50% to 65%.

The percentage of water is greater in childhood and decreases with age. Body weight variations can occur in a short period of time (hours), a medium period (days) or over a long period (weeks). Small variations in weight over a few hours is normal and due to water retention since extracellular water (which includes the interstitial fluid, When weight losses occur in just a few days, it is very likely that it is part of the total water. In **Fitdays™**, there will be an automatic comparison of your measured body water % against a range chart for your reference.

% Body Water		
	Men	Women
Normal	50-65%	45-60%

Body Mass Index

BMI Table	BMI Value
Underweight	<18.5
Normal	18.5-24.9
Overweight	25.0-29.9
Obese	>30

The body mass index (BMI) is the relation of body weight to height.

The classification given below is correlated to the risk that body weight exerts on health.

When BMI is between 18.5 and 24.9 it means that the weight does not statisti-

Any deviation, above or below this BMI interval should be considered a health risk. In this case, we recommend consulting a physician. Getting in the habit of following a balanced diet helps improve quality of life, promote health and prevent illness, increase vitality, and maintain ideal

It is important to know that BMI does not allow accounting for weight variations due to muscle mass, bone mass, water content or fat accumulation. Therefore, a simple BMI calculation alone should not be misinterpreted. The evaluation of BMI should be made by a doctor or other specialists such as a nutritionist or dietitian. In **Fitdays™**, there will be an automatic comparison of your measured BMI against a range chart for your reference.

Muscle Mass Percentage

Muscle mass is the weight of your muscles in you body. There are three different types of muscle: smooth, striated skeletal, and striated cardiac. Smooth muscle is made of smooth fibres. These are found in the walls of the digestive tract, respiratory tract, and genital tract and in the wall of the venous and arterial blood vessels. Human striated skeletal muscle is in the skeletal and facial muscles.

Cardiac striated muscle is found in the heart. Muscles play an important part in the overall health, including maintaining regular body functions. Muscle mass starts to decrease with age. Keeping track of muscle metrics overtime, along with other body composition metrics, can provide indications of your health fitness.

This scale captures and displays in the **Fitdays™** app 3 muscle metrics: muscle mass (lbs/kg), muscle rate (all muscle types) and skeletal muscle (%). On average, 40% of body mass is made up of skeletal muscles.

% MUSCLE MASS NORMAL RANGE		
Age	Men	Women
20-39	75-89%	63-75.5%
40-59	73-86%	62-73.5%
60-79	70-84%	60-72.5%

*The data in the table is for references only.

Bone Mass Information

Bone mass is the weight of bone mineral in your body. It is important that you develop and maintain healthy bones by having a balanced diet and plenty of exercises.

Keep track of your bone mass over time for any changes in the long term. Your bone mass tends to decrease after the age of 30. Below is a chart to show average of estimated bone mass based on weight and gender.

In **Fitdays™**, there will be an automatic comparison of your measured bone mass against a range chart for your reference also.

Bone Mass - Average bone mass			
Weight	Male	Weight	Female
<143 lbs	5.9 lb	<110 lbs	4.3 lb
143-209 lbs	7.3 lb	110-165 lbs	5.3 lb
>209 lbs	8.1 lb	>165 lbs	6.5 lb

*The data in the table is for references only.

MAINTENANCE

BATTERY REPLACEMENT

When the "LO" or battery symbol appears on the display, replace the batteries with the same type.

1. Open the battery compartment by pulling the tab on the cover.
2. Remove the worn batteries.
3. Insert the new batteries according to the polarity markings.
4. Close the battery compartment.

MAINTENANCE

- Store the device in a cool, dry area.
- Store the device in a horizontal position
- Only clean the scale with a soft, dry cloth.
- Never use chemicals or abrasive products

TROUBLESHOOTING

Problem	Possible cause	Solution
The scale doesn't turn on.	You have not stepped onto the scales properly.	Step onto the scale with your feet parallel. Note: the scale will turn on only after stepping onto it.
	The batteries have run out.	Replace the batteries.
"Lo" or the battery symbol is shown on the display.	The batteries are low.	Replace them.
The display shows "Err".	Weight greater than the scale capacity.	Step off the scale immediately so as not to damage it and check the maximum capacity.
The "Fitdays™" app doesn't display any data.	The connection has not been enabled in the smartphone settings.	Connect your smartphone.
	The smartphone went into standby during weighing.	Keep the app active on your smartphone during weighing to check that data is sent instantly.
The "Fitdays™" app only shows your body weight.	You did not step on the scale with bare feet.	Step onto the scale with bare feet.
	The skin on your feet is too dry.	Slightly dampen your feet.

All weights are approximate having been rounded up or down to the nearest whole number. We guarantee accuracy within $\pm 1\%$ + 1 division.

WARNING NOTE

Changes or modifications not expressly approved by the party responsible for compliance could avoid the user's authority to operate the equipment.

WARRANTY

This equipment is covered by a 5-year warranty from Springfield Instruments Inc. starting from the purchase date, which has to be proven by means of the stamp or signature of the dealer and by the fiscal payment receipt that shall be kept enclosed to this document.

This warranty period complies with the laws in force and applies only in case the consumer is an individual. Laica products are designed for domestic use only; therefore, their use for commercial purposes is forbidden. The warranty covers only production defects and does not apply, if the damage is caused by accidental events, misuse, negligence or improper use of the product. Use the supplied accessories only; the use of other accessories invalidate the warranty cover. Do not open the equipment for any reason whatsoever; the opening or the tampering of the equipment, definitely invalidate the warranty. Warranty does not apply to parts subject to tear and wear along the use, as well as to batteries, when they are included. Warranty decays after 5 years from the purchase date.

In this case, technical service interventions will be carried out against payment. Information on technical service interventions, under warranty terms or against payment, may requested by writing at info@springfieldinstruments.com or by contacting customer service at 1-844-687-4040. Corrective measures or replacement will be discussed with customer service representative if the product is deemed to be under warranty.

All interventions covered by warranty (including product replacement or the replacement of any of its parts) shall not be meant to extend the duration of the original warranty period of the replaced product.

The manufacturer cannot be made liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the

FCC STATEMENT

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference.
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

EN English

PS7011

YEARS GUARANTEE
ANS DE GARANTIE

LAICA S.p.A.
Viale del Lavoro, 10
36048 Barbarano Mossano (VI) - Italy
Tel. +39 0444.795314
info@laica.com
Made in China

www.laica.com

Distributed by / Distribué par: SPRINGFIELD INSTRUMENTS INC.

104 Barr Street, Saint Laurent, Quebec H4T 1Y4

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PÊSE-PERSONNE
NUMERIQUE INTELLIGENTE

LC 170/bTW - 06/22

The Bluetooth[®] word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Laica S.p.A. is under licence.
Other trademarks and trade name are those of their respective owners.
iOS, iPhone, iPad mini, iPod touch, App Store, Apple Health are trademarks of Apple Inc., registered in the U.S. and other countries - Android[™], Play Store[™], Google Play[™] and Google Health[™] are trademarks of Google Inc. - Fitbit[®] is a trademark of Fitbit Inc.
- Samsung Health is a trademark of Samsung Electronics Co.Ltd.

LAICA[®]

Distributed by/Distribué par:
SPRINGFIELD INSTRUMENTS INC
104 Barr Street, Saint Laurent, Quebec H4T 1Y4
Tel. 514-360-0223
Fax 514-360-0224
1-844-687-4040
info@springfieldinstruments.com

Made in China
Fabriqué en Chine

LAICA S.p.A.
Viale del Lavoro, 10
36048 Barbarano Mossano (VI) - Italy
Tel. +39 0444.795314
info@laica.com
Made in China

www.laica.com