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Product Introduction

This Smart Neck Massager adopts a streamlined circular design that conforms to the physiological curvature of the human cervical vertebra.

It adopts the five in one synchronous multi-channel health care functions of low-frequency electric pulse, magnetic effect, heat moxibustion, high-frequency vibration and circular traction to form a compound energy field.

During the whole working process, it always maintains the correct physiological posture, providing a platform for the health care of the cervical vertebra.

Our portable unit can help you to enjoy a relaxing, stimulating and comfortable neck massage anywhere!

Product Specifications

Product name: Smart Neck Massager

Charging input: DC5V / 1A

Rated power: 5W

Executive standard: GB4706.1-2005 GB4706.10-2008

Technical Performance Parameters

Pulse frequency: 1-1000Hz pulse

Pulse width: 20-400us

Pulse waveform: Square wave

Moxibustion temperature: + 38° C ~ 42° C (under 25° C environment)

Transportation and Storage Conditions

Ambient temperature range: - 20 ° C ~ + 55 ° C (60F-131F)

Relative humidity range: ≤ 80%

Atmospheric pressure range: 50kPa ~ 106kpa

Using the Remote Control

On the unit, turn on the power, and then press the power button on the remote control key once to light up the remote control screen.

Note: When using the remote control, the distance between the remote control and the unit should be kept within 5 meters(16.4 ft), otherwise the unit will not work.

1. Pulse mode adjustment (6 Modes: Beating / Massage / Acupuncture / Cupping / Scraping / Relaxation)
Press the power key to switch the pulse mode.
2. Pulse intensity adjustment (10 levels in total)
Press "+" "-" keys to adjust the intensity.
3. Time adjustment(10/15/20/30 minutes working time)
Press the "⌚" key to adjust the time.
4. Heat moxibustion regulation (2 gears in total)
Press the "🔥" key to switch on and off for the heat moxibustion function.

Tips:

1. The remote control uses the "CR2032" button battery, which can be used for about 60 days (30 minutes per day, under normal use).
2. When the brightness of the remote-control display is dimmed or the unit is not responding as usual, the remote-control battery is low. Please replace the battery with a new one.
3. If the Neck Massager is not to be used for a while, we recommend taking out the battery to prevent damage from a leaking battery.
4. Please do not use batteries other than the specified class.
5. Do not install the wrong polarity of the battery.
6. Do not use the battery in an environment above 45 °C[113F], otherwise it will affect the performance and life of the battery.
7. Please dispose of the battery after use according to the relevant environmental protection regulations of your city.

Charging

1. When the battery is low, you can use the 5V 1A charger to charge the unit. When charging, the power indicator light flashes white. When the battery is fully charged, the white light is always on. The charging time is about 2.5 hours.
2. After fully charged, it can be used continuously for about 5 hours, but the heating function will make the power consumption very fast, and it can be used continuously for 1 hour after turning on.
3. When the rechargeable battery is not used for a long time, the service life of the rechargeable battery will be shortened. It is recommended to charge at least once a month.

Precautions

1. The built-in intelligent recognition program of this machine will turn the unit off, if the machine is not used after 30 minutes.
2. Do not use Neck Massager with other medical electronic devices, such as cardiac pacemakers, artificial heart & lungs and electrocardiographs. Otherwise, it may cause medical problems or issues.

If the patient uses both high-frequency surgical equipment and the massager, it may cause burns in the massager's massage position and may damage the massager. If the massager is used in a shortwave or microwave therapy equipment accessory (1 meter), the output of the massager may be unstable.

3. Pregnant women, sensitive skin, heart disease, abnormal blood pressure, malignant tumors, cerebrovascular patients, patients with acute illness, or other patients who are being treated by a doctor must consult a doctor before using this product.
4. Do not use when bathing or sleeping.
5. Do not use it on unsupervised children, seniors or anyone without the ability to understand clearly instructions on how to use the machine.
6. Do not use Neck Massager if you suffer from skin irritation or are sensitive to heat.
7. Do not adjust the position of the Neck Massager when in use. Please turn off the machine and then proceed to make any adjustments. Otherwise, you may feel discomfort and or pressure from the unit.
8. Caution: when adjusting the pulse intensity, we advise to go slow and easy, carefully adjusting level as needed. Stronger does not necessarily mean better. The high-end level may be too strong, so go easy when starting out!
9. Do not under any circumstances disassemble this product to avoid any injuries and void the warranty.
10. Do not wash this product with water.

Cleaning and Maintenance

1. When cleaning the product, wipe it with a soft cloth or towel with little water, and do not use the volatility liquid.
2. Please keep it in a dry and ventilated place.
3. Please do not store in direct sunlight, high temperature and damp places.

Product Scrap Treatment

To help our planet and in order to avoid environmental pollution, when the Neck Massager and its accessories are damaged or can no longer be used, please make sure to dispose safely, according to local environmental requirements.

Name and content of harmful substances in products

Part Name	Hazardous substance					
	[Pb]	[Hg]	[Cd]	[Cr(VI)]	[PBB]	[PBDE]
Host	X	0	0	0	0	0
Parts	X	0	0	0	0	0
Packing	X	0	0	0	0	0

This form is prepared in accordance with SJ / T11364

0: It means that the content of the hazardous substance in all homogeneous materials of the part is within the limit requirements specified in GB / T26572.

X: It means that the content of the hazardous substance in at least one homogeneous material of the part exceeds the limit requirements specified in GB / T26572, and there is no mature alternative in the industry, meeting the requirements of GB / T26572.

The number in this sign indicates that the environmental protection service life of the product under normal use is 10 years. Some components may also have an environmental protection service life mark, and the number in the mark shall prevail. Due to different models, the product may not include all parts except the main engine. Please refer to the actual sales configuration of the product.

FAQs

Q: I'm using the Neck Massager, and I feel a funny sensation on my neck, is this normal ?

A: The reason for this tingling sensation is that your skin may be dry and the massager is not making complete contact with your neck. To avoid this, it is recommended to use a moist towel, allowing your neck to be in close contact with the massager. Start with a low-level pulse intensity and gradually adjust. If you feel uncomfortable, you can select other functions. (Electric pulse massage is to use Chinese traditional acupuncture to massage the Main and collateral channels of the neck)

Q: Can this unit cause a rash?

A: After using the electric pulse function, if you develop a rash, it may be that your skin is too dry, or the massager is not making contact directly and is loose. We recommend putting your neck with a wet towel to make the massager stick to the skin and make good contact. If a rash develops, it should go away on its own after 2 days.

Q: How long does the battery take to be fully charged? How long does a full battery charge last?

A: The neck massager has built in environmental protection lithium battery, and the full charging time is about 2 hours. Single pulse function based on full charge can last around 5 hours and the heating function can last about one hour.

Q: Will this Neck Massager cause pain?

A: The neck massager adopts ergonomic design, perfectly fits different neck sizes of all sizes. ABS anti-bacterial engineering plastics are used to protect the neck effectively; please wear this product correctly when using and keep the correct posture. To feel great results, we recommend using it twice a day for 15 minutes each time to relax neck.

Quality Assurance

Thank you for purchasing the intelligent cervical vertebra massager, please read the product manual carefully before use.

Under normal use, you can enjoy your Neck Massager that comes with lifetime warranty. If you need further assistance, please feel free to contact Rova customer service email address:

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This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Warning: changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment

Shenzhen Apeman Innovations Technology Co.,Ltd

1801, Heng Lu E Times Building, No.159, North Pingji Road, Hehua Community, Pinghu Street, Longgang District, Shenzhen, Guangdong, China

Lifetime Warranty

For any inquiry or customer service, please email us:
support.us@rovahealth.com
support.all@rovahealth.com

Made in China

version 1.0