0025 Scale Cover / Info Card

PMS 7464C 60% (Match Tray, Tray insert & Outer Sticker)	PMS 2965C (Match Tray insert & Outer Sticker)	PMS 178C
PMS 7464C (Match Tray insert & Outer Sticker)	PMS Cool Grey 428C	PMS 136C (Match Slip Sticker

Front Model No: welcome guide see reverse

Back **Quick Start:** WARNINGS Do not step on the edge of the scale while getting Remove the plastic pull tab from the back of the scale to activate the batteries. on or off-otherwise, it may tip. Do not jump on scale. Protect scale from hard knocks, temperature fluctuations, and heat sources that are too close (e.g. stoves, heating units). • Do not drop scale or drop any objects on it as this Set the scale on a hard, flat surface, then may damage the sensors. step on (not carpet). • Do not step or stand on the scale when your body and/or feet are wet-for example, after taking shower. • Do not step on the scale when the surface is wet. There is a danger of slipping. When you see the weight estimate icon or hear double beeps, step off. • Do not step on this scale if you weigh more than Clean the scale with a slightly damp cloth. DO NOT immerse the scale in water or use chemical/abrasive cleaning agents. • Do not heat or deform the batteries, or dispose of The green AccuCheck icon or a long beep means your weight has been confirmed. them in fire. Batteries should not be disposed of with household waste. Please check with your local authority for battery recycling advice. Please recycle where facilities exist. Your weight will be sending after your weigh-in, and "Sent" will appear when your weight has been transmitted. Visit www.me.health/help for more information. FCC Regulations:
This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation If this equipment does cause harmful interference by the sequipment of fand on, the user is determined by turning the equipment off and on, the user is edilowing measures:

• Reorient or relocate the receiving antenna.

• Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

• Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

• Connect the equipment and receiver.

• Changes or modifications no

310.0

Dielines