

User Manual

Baby Scale

0220



FCC ID:2ADUL0220

Greater Goods, LLC.
4427 Chouteau Ave. St. Louis, MO63110









- Thank you very much for selecting the Baby Scale 0220.
- Please do read the user manual carefully and thoroughly so as to ensure the safe usage of this product, and keep the manual well for further reference in case you have problems.

Table of Contents

Symbol Meaning	
Symbol Meaning	2
Safety Information	
Your Scale and Its Environment	3
Efficient Use of Your Scale	3
Overview	
General Instructions	4
Device Components	4
List	5
LCD Display	5
Initial Start-Up	
Insert the Batteries	6
Install APP and Pair-up	7
Start Measuring	
Select Measurement Unit	8
Daily Measurement	8
Measuring With ZERO Function	9
Data Management	
Data Transmission	10
Troubleshooting	
Error Prompt	11
When Measuring	11
When Data Transmitting	12
Specifications	13
Maintenance	14
Warranty	14
FCC Regulations	15

♥ Symbol Meaning

The warning signs and symbols are essential to ensure your correct and safe use of this product and protect you and others from injury. Please kindly find the meanings of the warning signs and symbols, which you may encounter in the label and user manual, as follows:

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "MANUFACTURER"
	The Bluetooth Combination Mark		Symbol for "COMPLIES WITH EU REQUIREMENTS"
	Symbol for "SERIAL NUMBER"		Symbol for "ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "DIRECT CURRENT"		

♥ Your Scale and Its Environment

To ensure your safety as well as the service life of your scale, please avoid using the scale under the following circumstances:

- Slippery floor such as tile floor
- Jumping onto the platform immediately after bath or with wet feet
- Near a cell phone or microwave oven

Avoid storage in the following locations:

- Where there is water
- Where the device may be exposed to extreme temperatures, humidity, moisture, direct sunlight, dust, or salt air
- Where there is risk of shock or drop
- Where you store chemicals or full of corrosive gases
- Where in reach of the infants or children

♥ Efficient Use of Your Scale


To ensure the accuracy of measurement, please follow below instructions when you start measurement.

- Place the scale on a flat, hard surface. Soft surface such as carpet will affect the performance of the scale.
- Start measurement at least two hours after Getting up or Dinning.
- Avoid measurement immediately after strenuous exercise, sauna or bath, drinking, and dining.
- Always start measurement in the same time slot and on the same scale located on the same flat, hard surface.
- Any information provided by this scale is in no way meant to treat, cure or prevent any disease or illness from happening. If in doubt, contact your physician.
- It is intended for use in the domestic setting only.

♥ General Instructions

The Baby Scale 0220 applies SENSE ON technology. It's not necessary to turn on or turn off the scale, just stand on the scale for measurement. It is delicate and trendy, convenience for using.


♥ Device Components



♥ List

- 1. Baby Scale 0220
- 2. Four AAA-size Batteries (1.5V per each)
- 3. User Manual

♥ LCD Display



kg	Kilogram		Bluetooth Transmitting
lb:oz	Ounce		Data pending to transmit to wireless wellness system
lb	Pound		Low battery

♥ Insert the Batteries

- Open the battery door in the back of the scale.
- Remove insulating strip from beneath the battery (if fitted) or insert batteries (4 x 1.5V AAA) observing the polarity signs (+ and -) inside the battery compartment.

* The digits "HELLO" will be shown on the LCD.



- Close the battery door and wait until the digits "0: 0.0lb:oz" are shown on the LCD.

* The scale will turn off about 120s if there is no operation after the batteries are installed.



CAUTION

- When the symbol "⚡" appears, the device will power off in three seconds. Then you shall replace with a new set of batteries. Please replace all four batteries at the same time. Do NOT mix the old batteries with the new one.
- Worn batteries are hazardous waste. Do NOT dispose of them together with the household garbage. Please refer to the local ordinances and recycling instructions regarding disposal of the worn batteries.


♥ Install App and Pair-Up

With the advanced Bluetooth 4.0 technology applied, the mobile or portable equipments, which are equipped with Bluetooth function in line with BLE Technical Specifications as well as BLP Protocol established by global organization Bluetooth SIG, are capable to receive your personal health data.

Just simply install the specially-designed app and pair up your scale with your mobile or portable equipments. Then you may enjoy the comprehensive health solution provided by Transtek.


Note: It is only available for the device which supports Bluetooth 4.0 technology.

- The APP "GG Baby Scale" can be freely downloaded at the APP Store, you may search and install it in your smartphone.
- Turn on Bluetooth and APP, then touch "SCAN" you will see a "BS1711-B" device name. The device will connect with your smartphone automatically.



♥ Select Measurement Unit

With batteries correctly installed, press "UNIT" button in the back of the scale to select measurement unit. The default measurement unit is "Ounce". You may press "UNIT" button to choose among ounce, kilogram and pound.



♥ Daily Measurement

- Turn on the scale, put the baby upon the scale until the LCD displays the weight and stop blinking.



- If the baby moves around and the weight is unstable, you can touch "HOLD" button to get a stable value.
- You can do the next weighing in a few seconds.

♥ Measuring with ZERO Function

ZERO function allows you to weigh different items consecutively with out removing the items from the scale.


1. When the scale is off, press the ON/OFF button to go to the weighing mode. You can start to weigh when "0: 0.0" is shown on the LCD.



2. Place the first item you wish to weight onto the scale, the LCD will show the weight of the first item.



3. Press the ZERO button to reset the display to zero.




4. Now add another item to weigh onto the scale. The LED will show the weight of the newly added item.





5. If there is no operation, the scale will switch off automatically in about 2 minute. Or you can press and hold the ON/OFF button for about 2s to switch the scale off.

♥ Data Transmission


- Once the Baby Scale 0220 is successfully connected to your smartphone, please open the APP and turn on your Bluetooth. Keep smartphone and scale in transmission distance. Once the measurement is finished you can touch "GET DATA" to get your measurement data.



- The symbol  only appears when the data is transferred, and disappears after transmission. Then you may check your health data stored in your smartphone.
- If the data transmission fails, the symbol  will remain. The pending measurement data will be temporarily kept in the scale. But the data is not saved after shutdown.






CAUTION

- Interference may occur in the vicinity of equipment marked with the following symbol . And the Analyzer may interfere vicinity electrical equipment.
- To enable the data transmission function, this product should be paired to a Bluetooth 4.0 end.

How to mitigate possible interference?

- The range between the Analyzer and the Bluetooth end should be reasonably close, from 1 meter to 10 meters. Please ensure no obstacles between the Analyzer and the Bluetooth end so as to obtain quality connection.
- To avoid interference, other electronic devices (particularly those with Bluetooth transmission / Transmitter) should be kept at least 1 meter away from the Analyzer.

♥ Error Prompt

Error	Description	Solution
	Overload. The device will power off.	Stop using this scale for measurement.
	Low Battery. The device will power off in three seconds.	Replace all four batteries at the same time. Please purchase the authorized batteries for replacement.
	There has data doesn't upload. (The notebook icon will be displayed.)	Upload the data.

♥ When Measuring ...

Problem	Root Cause	Solution
No display on LCD when the device powers on.	Batteries not yet installed.	Install the batteries. (Please refer to Insert the Batteries)
	Worn batteries.	Replace all four batteries at the same time. Please purchase the authorized batteries for replacement.
The device powers off automatically.	Low battery.	Replace all four batteries at the same time. Please purchase the authorized batteries for replacement.

♥ When Data Transmitting ...

Problem	Root Cause	Solution
Data transmission failed.	Bluetooth is OFF.	Turn ON the Bluetooth via Settings.
	APP is OFF.	Press the icon to turn ON your APP.
	Out of range of Bluetooth transmission.	Place your smartphone closer to the scale.

♥ Specifications

Product Name	Baby Scale
Dimension	Scale: Approx. 609.7x294.9x84.8mm
Net Weight	Approximately 2.3kg (Excluding the dry cells)
Display	Digital LCD
Measurement Unit	Ounce / Kilogram / Pound
Measurement Range	66.138lb / 30kg
Division	5g
Accuracy	Temperature requirement: 25±5 °C 10kg: ±0.04kg; 20kg: ±0.07kg; 30kg: ±0.1kg;
Working Environment	Temperature: 5 °C to 40 °C Relative Humidity: ≤90% RH Atmospheric pressure: 86kPa to 106kPa
Storage Environment	Temperature: -20 °C to 60 °C Relative Humidity: ≤90% RH Atmospheric pressure: 50kPa to 106kPa
Power Source	6V (Four AAA-size Batteries)
Auto-OFF	About 2 minutes while showing 0.0; About 2 minutes after the weight data is locked.
Accessories	1. Four AAA-size Batteries 2. User Manual
Mode of Operation	Continuous Operation

About the Accuracy of This Product

- This product passes strict inspection before delivery and therefore its accuracy is guaranteed by the manufacturer. Please refer to the above table for the descriptions on accuracy.
- This product is specially designed for body scale as well as weight measurement. It should NOT be used by anyone during the process of transaction for verification of goods' weight.

♥ Maintenance

When carrying out usual maintenance, please ensure practice of the following Do's and Don'ts:

- DO use a dry soft cloth to wipe the dust.
- DO use a wet soft cloth, dipped into water and wrung out, to wipe the dirt. Then use a dry soft cloth to dry up the device.
- DON'T wash the device with water or immerse it in water.
- DON'T use propellant, abrasive or other chemicals to wipe the dirt in avoidance of discolor or malfunction.
- DON'T disassemble this device. If you have any problems, please contact us. *(Please refer to [Warranty](#) for contact information)*

♥ Warranty

- Transtek warrants its products free of defects in materials and workmanship in normal use for a period of FIVE years from the date of retail purchase.
- This warranty does NOT cover damages caused by misuse or abuse, including but not limited to:
 - Failure caused by unauthorized repairs or modifications;
 - Damage caused by shock or drop during transportation;
 - Failure caused by improper operation inconsistent with the instructions stated in this user manual;
 - Malfunction or damage from failure to provide the recommended maintenance;
 - Damage caused by improper use of power supply.
- Should this device require maintenance (or replacement at our option) under warranty, please deliver the original package to **Greater Goods, LLC** prepaid. Please return the store receipt (with the retail purchase date) and a note with reasons to return on it as well.

Greater Goods, LLC.
4427 Chouteau Ave. St. Louis, MO63110

♥ FCC Regulations

FCC User Guide Information

Radio Frequency Interface Requirements - FCC

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna;
- Increase the separation between the equipment and receiver;
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected;
- Consult the dealer or an experienced radio / TV technician for help.

Radio Transmitters (Part 15)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

FCC RF Exposure Guidelines

Safety Information

Reducing RF Exposure - Use Properly

Only operate the device in accordance with the instructions supplied.

This device complies with FCC radiation exposure limits set forth for an uncontrolled environment.