

FOREO

UFO

A WHOLE NEW WAY TO MASK IN ONLY 90 SECONDS!

[insert image]

FULL USER MANUAL

GETTING STARTED

Congratulations on taking the first step toward out of this world skincare by acquiring UFO. Before you begin to enjoy all the benefits of sophisticated mask technology in the comfort of your home, please take a few moments to carefully read the instructions in this manual.

Please READ ALL INSTRUCTIONS BEFORE USE and utilize this product only for its intended use as described in this manual.

WARNING: NO MODIFICATION OF THIS EQUIPMENT IS ALLOWED.

UFO OVERVIEW

The world's most sophisticated smart mask, UFO combines highly effective skincare technologies with exclusive mask formulas for a professional-level skincare experience in seconds. Featuring enhanced Hyper-Infusion Technology, UFO integrates Thermo-Therapy and Cryo-Therapy to push mask essence deep below skin's surface and lock in moisture while reducing the appearance of pores. Full-spectrum LED light therapy completely rejuvenates skin, unmasking a refined, radiant complexion.

GETTING TO KNOW YOUR UFO

[insert diagram of device front & back - then picture of app]

FEATURES

Medical Grade Silicone - Ultra-hygienic and super soft

T-Sonic Pulsations - Gently massage mask essence into skin

RGB LED Light- Red, green and blue LED lights offer a range of skincare benefits

Thermo/Cryo Therapy - Heating and cooling to open pores and firm the skin

UFO Activated Mask - Soft microfiber masks ensure revitalized and refreshed skin

Smart Swedish Design - Ergonomic and lightweight, so you can mask on the go

FOREO App- Syncs pre-programmed customized beauty routine to device automatically

HOW TO USE:

PREP FOR USE

1. REMOVE THE PLASTIC RING FROM UFO.
2. REMOVE UFO ACTIVATED MASK FROM SACHET
3. SECURE MASK ON DEVICE BY CLIPPING THE PLASTIC RING IN PLACE
4. PRESS ONCE TO START THE SMART BEAUTY TREATMENT

[how to use section images****]

ACTIVATING YOUR TREATMENT

MAKE MY DAY

- Press the button to turn on the device.
- The first quadrant will be illuminated by a pink light.
- After 10 seconds, the device will be pre-heated and your treatment will begin.

CALL IT A NIGHT

- Press the button once more.
- The first two quadrants will be illuminated by pink lights.
- After 10 seconds, the device will be pre-heated and your treatment will begin.

SKIN REHAB COLLECTION

- Press and hold the button for 3 seconds. A blinking white light indicates UFO is in pairing mode.
- A solid white light indicates UFO has been successfully connected to your smartphone or tablet.
- Use the FOREO app to sync your selected treatment routine.

#

*Indicator lights encircling UFO let you know which treatment you're about to enjoy! (Image example)

#

APPLYING YOUR UFO TREATMENT

1. Clean and dry your face and neck. Gently glide UFO across your skin, distributing the essence evenly.
2. Using the UFO device, massage the formula into your skin with a circular motion. Move the device in a clockwise circle around the face, starting on the chin, then moving up to the right cheek and the forehead. Continue down the left cheek to the chin and neck. Repeat this movement until treatment ends.
3. When the treatment ends, gently pat any remaining essence into your skin, or remove excess with a cotton pad.
4. Follow the mask with the moisturizer of your choice.
5. Remove the plastic ring from device and throw away the mask. Then, rinse both UFO and plastic ring under warm running water. Pat dry, re-attach ring to the device and place back on the plastic stand.

CLEANING UFO

Always clean UFO thoroughly after use. Wash the device and plastic ring with soap and then rinse with warm water. Pat dry with a lint-free cloth or towel. After use, we recommend placing device on plastic stand to dry.

NOTE: Never use cleaning products containing alcohol, petrol or acetone, as they may irritate the skin and damage the silicone.

IMPORTANT

FOR OPTIMUM SAFETY

- If you have a skin condition or any medical concerns, please consult a dermatologist before use.
- Using UFO should be comfortable – if you experience discomfort, discontinue use immediately and consult a physician.
- Exercise particular care when cleansing the under-eye regions and do not bring the device into contact with the eyelids or eyes themselves.
- For reasons of hygiene, we do not recommend sharing your UFO with anyone else.
- Avoid leaving your UFO in direct sunlight and never expose it to extreme heat or boiling water.
- Close supervision is necessary when this device is used by, on, or near children, as well as those with reduced physical and mental abilities.
- Discontinue use if this product appears damaged in any way. This product contains no serviceable parts.
- Use this device only for its intended use as described in this manual. If you do not find the answer to your specific question, or if you have any other questions regarding the device's operation, please visit foreo.com/support.

TROUBLESHOOTING

Precautions to be taken in the event of changes in the performance of UFO.

If UFO is not activated when pressing the universal button:

- **Battery is empty. Charge your UFO via USB charger until the light emits a steady glow. A full 2.5 hour charge will last up to 40 treatments.**

If UFO cannot be switched off and/or interface buttons do not respond:

- Microprocessor is temporarily malfunctioning. Press and hold power button to restart the device.

If UFO won't sync to the FOREO app:

- Switch your Bluetooth off and then on again to try reconnecting
- Close the FOREO app and then reopen it to start the process over

WARRANTY TERMS & CONDITIONS

REGISTER WARRANTY

To activate your 2-Year Limited Warranty and 10-Year Quality Guarantee, register the number provided on the FOREO Magnetic Scratch Card at www.foreo.com/support.

2-YEAR LIMITED WARRANTY

FOREO warrants this device for a period of TWO (2) YEARS after the original date of purchase against defects due to faulty workmanship or materials arising from Normal Use of the device. The warranty covers working parts that affect the function of the device. It does NOT cover cosmetic deterioration caused by fair wear and tear, or damage caused by accident, misuse or neglect. Any attempt to open or take apart the device (or its accessories) will void the warranty.

If you discover a defect and notify FOREO during the warranty period, FOREO will, at its discretion, replace the device free of charge. Claims under warranty must be supported by reasonable evidence that the date of the claim is within the warranty period. To validate your warranty, please keep your original purchase receipt together with these warranty conditions for the duration of the warranty period.

To claim your warranty, you must log in to your account at www.foreo.com and then select the option to make a warranty claim. Shipping costs are nonrefundable. This undertaking is in addition to your statutory rights as a consumer and does not affect those rights in any way.

10-YEAR QUALITY GUARANTEE

As a continuation of the terms of FOREO's TWO (2) Year Limited Warranty, FOREO's TEN (10) Year Quality Guarantee makes owners eligible to purchase a new device from www.foreo.com for 50% of the listed price.

DISPOSAL INFORMATION

Disposal of old electronic equipment (applicable in the EU and other European countries with separate waste collection systems).

[insert all symbols here]

The crossed-out dustbin symbol indicates that this device should not be treated as household waste, but rather be brought to the appropriate collection point for recycling of electrical and electronic equipment. By ensuring this device is disposed of correctly, you will help prevent the potential negative consequences for the environment and human health which could be caused by inappropriate waste handling of the product. The recycling of materials will also help conserve natural resources.

For more information about the recycling of your device, please contact your local household waste disposal service or your place of purchase.

[INSERT UFO BATTERY REMOVAL ILLUSTRATIONS]

SPECIFICATIONS

MATERIALS: Body-safe silicone & PC + ABS

COLOR: Pearl Pink/Mint/Fuchsia

SIZE: 72mm(L)* 72mm(W)* 31 mm(H)

WEIGHT: 146.2g (TBC on UFO mini) gr

BATTERY: Li-Ion 1200mAh 3.7V

USAGE: up to 40 uses per charge

STANDBY: 180 days

FREQUENCY: 180 Hz

MAX NOISE LEVEL: <50 dB

INTERFACE: 1-button

Disclaimer: Users of this device do so at their own risk. Neither FOREO nor its retailers assume any responsibility or liability for any injuries or damages, physical or otherwise, resulting, directly or indirectly, from the use of this device. Further, FOREO reserves the right to revise this publication and to make changes from time to time in the contents thereof without obligation to notify any person of such revision or changes.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

RF exposure statement: This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Model may be changed for improvements without notice.

©2018 FOREO AB. ALL RIGHTS RESERVED.

IMPORTER AND DISTRIBUTOR IN EU: FOREO AB, BIRGER JARLSGATAN 22, 114 34 STOCKHOLM, SWEDEN IMPORTER AND DISTRIBUTOR IN US: FOREO INC., 3565 S. LAS VEGAS BLVD. #159, LAS VEGAS, NV 89109, USA MANUFACTURED FOR FOREO AB DESIGNED AND DEVELOPED BY FOREO SWEDEN.

Note: Don't touch the Top Panel gold of the UFO directly when it is powered on.

WWW.FOREO.COM