


### Genesis AI Custom ITE or ITC Rechargeable Left Hearing Aid

The custom hearing aid is custom manufactured to the patient's ear canal. Whether the aid is called an ITE or ITC depends on the size of the patient's ear:

- ITE (In-the-ear) devices reside primarily in the bowl and past the 1<sup>st</sup> bend of the canal.
- ITC (In-the-canal) devices reside partially in the bowl and into the 2<sup>nd</sup> bend of the ear canal.

These photos show the interior circuitry that is the same for an ITE or ITC




left side



right side



Front view



Back view