

-No Automatic Measurements: Confirm that measurement intervals are correctly set in the app.

-Low Battery: Recharge the ring using the supplied charging dock.

9. Maintenance and Care

- Clean the ring with a soft, dry cloth.
- Avoid exposure to chemicals or abrasive materials.
- Store in a cool, dry place when not in use.

10. Regulatory Information



CE Certification: [CE ID]

11. Warranty and Support

Warranty Period: 1 year from the date of purchase  
Customer Support: [viv@vivring.ai](mailto:viv@vivring.ai)  
Support Phone Number: +82 2 6487 6878

Thank you for choosing VIV Ring. Stay connected to your

# VIV Ring

Model: VIV-Ring-6

Input: 5V~0.4A (Max)

VIV HEALTH INC

FCC ID:2BN8V-VIVRING

Made In China

size: 50mm\*30mm

Due to the small size of the product, the FCC ID number is located on the instruction manual