


2:1 view

View only - use layout file for print.



Printed subjects:

- Use layout file (1:1) 83-03743 as input;
- Use CE mark (4mm height) as reference to ensure proper size;
- Bright for product name and dimmed for TA marks & others;
- Current production day "WXXY", capital letter "W" follow by week/year code.
Year code using base36 coding. (i.e.: "W35H" for week 35 year 2017)

CONFIDENTIALITY NOTICE:

This document is the property of GN A/S and must be treated as confidential. The document must not be used, or disclosed to any third party, without written permission of GN A/S. Copyright © 2007 GN A/S. All rights reserved.

Process control dimensions are marked: ☐

IQC/OQC Characteristics for Evaluation:

1, 2

Date of Change:

Description of Change:

Responsible:



Matr. Spec: (See B.O.M.)

Matr. type:

3D file name:

Cosmetic Standard: GSC.ST.00002
Part Marking: GSC.PR.00102

Issued by:

ECAI

Issued:

04.09.17

Drawing Scale:

2:1

Size:

A4

Approved by:

Date Approved:

Projection:



Description:

Laser Placement, HS LW, Rio

Document No:

50-07917

Doc. Revision:

C

Units:

mm

Page No:

1 of 1

A

B

C

D

E

F