



AirSense™10
AirCurve™10

Addendum

ENGLISH

This addendum replaces the wireless module information found in the User and Clinical Guides.

Wireless module

Technology used:	Frequencies (MHz)*	Max RF power output (dBm)
2G	850/900/1800/1900	33.5
3G	850/900/1900/2100	33.5
4G	700/800/850/900/1700/1800/1900/2100/2300/2600	33.5

*Bands may not be available in all regions.

FCC ID: 2ACHL-AIR104GU

IC ID: 9103A-AIR104GU

This device complies with FCC Rules and Industry Canada rules.

This device should be used at a minimum distance of 0.8" (2 cm) from the body during operation.

Additional information regarding the FCC Rules and IC compliance for this device can be found in

www.resmed.com/downloads/devices

FRANÇAIS

Le présent addenda remplace les informations relatives au module sans fil figurant dans le Guide de l'utilisateur et le Guide clinique.

Module sans fil

Technologie utilisée:	Fréquences (MHz)*	Puissance de sortie RF max. (dBm)
2G	850/900/1 800/1 900	33,5
3G	850/900/1 900/2 100	33,5
4G	700/800/850/900/1 700/1 800/1 900/2 100/2 300/2 600	33,5

*Les bandes sont susceptibles de ne pas être disponible dans toutes les régions.

ID FCC : 2ACHL-AIR104GU

ID IC : 9103A-AIR104GU

Cet appareil est conforme aux règlements de la FCC et d'Industrie Canada.

Cet appareil doit être utilisé à une distance minimale de 2 cm (0,8 po) du corps pendant son fonctionnement.

Le site www.resmed.com/downloads/devices contient des renseignements supplémentaires au sujet de l'observance des règlements de la FCC et d'IC applicables à ce dispositif.



ResMed Pty Ltd

1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia

See ResMed.com for other ResMed locations worldwide. AirSense and AirCurve are trademarks and/or registered trademarks of the ResMed family of companies. For patent and other intellectual property information, see ResMed.com/ip.

© 2022 ResMed. 3781050/1 2022-05

ResMed.com

Rx Only

