

CARETALK®
Wireless
Pain Relief Device

CARETALK®
Wireless
Pain Relief Device

6 preset TENS programs
5 preset EMS programs
3 preset massage programs
for targeted pain relief
Does not contain natural
rubber latex
TENS therapy has been proven
effective by the American
Physical Therapy Association.

Includes:

Wireless device with
rechargeable lithium ion battery
USB power cord
Remote with 2 AAA batteries
included
2 adhesive electrode pads


- Neck, Lback, Arm, Shoulder,
Leg, Joint treatment areas
- 60 levels of intensity
- Bluetooth connectivity
- Effective in temporarily relieving
pain and stimulating muscles

This device complies with Part
15 of the FCC Rules. Operation
is subject to the following two
conditions: (1) this device may
not cause harmful interference,
and (2) this device must accept
any interference received,
including interference that may
cause undesired operation.

CARETALK®
Wireless
Pain Relief Device



Neck



Lback



Arm



Shoulder



Leg



Joint

Used with the Remote or the APP

CARETALK®
Wireless
Pain Relief Device

TENS for Pain Relief
EMS for Muscle Recovery
Massage for Relaxation

What is TENS?

TENS stands for Transcutaneous Electrical Nerve Stimulation, a drug-free method of controlling pain. TENS uses tiny electrical impulses sent through the skin to the nerves to modify pain perception. It is effective in temporarily relieving pain, allowing for a return to normal activity.

What is EMS?

EMS is Electronic Muscle Stimulation, which uses electrical impulses to contract muscles in order to improve and facilitate muscle performance. This aids in strength training, recovery from minor injuries or post exercise, as well as helping to temporarily relieve pain and has been shown to enhance blood circulation.


What is Electrotherapy Massage?

The massage setting provides low level muscle vibration to loosen tight muscles. Its main function is to mimic manual massage to promote relaxation and reduce muscle tension. The massage function is great for use before or after a workout session to help keep muscles from tightening.

Compatibility: iOS and Android™ compatible



Shenzhen Dongdixin Technology Co., LTD.
No.3 Building XiliBaimang Xusheng Industrial
Estate, 518108 Nanshan, Shenzhen China
Tel: 0086-755-27652316



Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg Germany
Tel: 0049-40-2513175 Fax: 0049-40-255726

**Surround your pain
using two electrodes**



55°C
-10°C

90%
10%

1060hPa
700hPa

