
 Medtronic	MEDTRONIC CONFIDENTIAL THIS DOCUMENT IS THE PROPERTY OF MEDTRONIC, INC. AND MUST BE ACCOUNTED FOR. INFORMATION HEREIN IS CONFIDENTIAL. REVISION LEVELS ARE ELECTRONICALLY CONTROLLED. DO NOT REPRODUCE IT, REVEAL IT TO UNAUTHORIZED PERSONS OR SEND IT OUTSIDE MEDTRONIC WITHOUT PROPER AUTHORIZATION.	p. 1 of 6
		Doc. No.
		M976704A
TITLE DEVICE REGULATORY LABEL MODEL PA97000 PATIENT ASSISTANT		REV B

REVISION	DESCRIPTION
A	Initial release
B	Revise section 2 to reflect metric values, add tabs 002 through 006 to view 2, update to new template.

1.0	TABLE OF CONTENTS	
1.0	Table of Contents	2
2.0	Scope	3
3.0	Component Identification	3
4.0	Handling Requirement	5
5.0	Notification of Change	6
6.0	Completion	6
7.0	Notes	6

	DEVICE REGULATORY LABEL MODEL PA97000 PATIENT ASSISTANT	M976704A	REV. B	p. 3 of 6
---	--	-----------------	---------------	------------------


2.0 SCOPE

This specification shall identify the procurement source for Device Regulatory Label for the Model PA96000 Patient Assistant.

3.0 COMPONENT IDENTIFICATION

1. Material: (.125 mm) Thick velvet matte finish polycarbonate with PSA adhesive. Overall thickness excluding release liner (.185 mm) nominal.
2. Adhesive: (.06 mm) thick Acrylic PSA, 3M 467MP or Medtronic approved equivalent.
3. Colors: Background – Dark Grey, Pantone 431. Graphics – Black, Pantone Process Black.
4. See graphics file for artwork. Refer to associated artwork file (.eps file).
5. Packaging Kiss-cut on sheets with waste removed or die cut individual labels with release liner pull tab. Bulk package as to prevent damage and packaging waste.
6. All dimensions in View 1 are in millimeters unless otherwise noted. Tolerances are +/- .254 mm.
7. Radius tolerances are +/- .3 mm unless otherwise noted.
8. Identify each label package with the part number at a minimum.
9. All borders and window/s are represented by line-art. See the supplied artwork. Cut back the adhesive layer approx. (.508 mm) REF. from window border. No adhesive allowed in window area.
10. This component to be RoHS compliant with restrictions on hazardous substances.


View 1 – Dimensions – (not to scale)




	DEVICE REGULATORY LABEL MODEL PA97000 PATIENT ASSISTANT	M976704A	REV. B	p. 4 of 6
---	--	----------	--------	-----------

View 2 – Artwork Tabs (Reference use only not to scale)


TAB -001 (DEMO)




TAB -002 AUS / NZ




TAB -003 China




 Medtronic	DEVICE REGULATORY LABEL MODEL PA97000 PATIENT ASSISTANT	M976704A	REV. B	p. 5 of 6
--	--	----------	--------	-----------


TAB -004 Japan



TAB -005 US / Canada




TAB -006 OUS



4.0 HANDLING REQUIREMENT

The manufacturer shall package and ship the component parts in a manner designed for protection from physical damage and, if electrically susceptible, from electrostatic discharge.

This document is electronically controlled. Printed copies are considered uncontrolled.

 Medtronic	DEVICE REGULATORY LABEL MODEL PA97000 PATIENT ASSISTANT	M976704A	REV. B	p. 6 of 6
--	--	-----------------	---------------	------------------

5.0 NOTIFICATION OF CHANGE

Upon approval by Medtronic, Inc. of the initial design, any part, process changes or deviations considered by the manufacturer must be submitted to Medtronic, Inc. In writing for review, the information submitted should include a complete description of the change and effect the change will have on all characteristics of the device, material and /or process.

6.0 COMPLETION

This paragraph concludes this specification.

7.0 NOTES

Not applicable to this document.