



Antenna Info

List of Antennas:

Coil loop antenna located in programmer head

- Manufacturer: Inspire Medical Systems, Inc.
- Frequency: 175.8 kHz
- Number of Turns: 36 (see attached drawing)
- Description: Loop antenna

Drawing or Photograph – See “Telemetry Head Measurements.pdf”.

FCC 15.203 Compliance – This antenna is permanently attached, thus complying with §15.203.

Two ceramic loop antennae located in the programmer cable controller

- Gain: Peak 2.5 dBi
- Description: Loop Antenna

Drawing or Photograph – See “Antenna 1+2 Datasheet.pdf”.

FCC 15.203 Compliance – These antennae are permanently attached, thus complying with §15.203.