



*Title:* **SIMILARITIES AND DIFFERENCES -  
US/CAN VARIANTS**

*Project/Product* **AIR 10**

## 1. Purpose

This document discusses the similarities and differences between US/CAN variants of the Air10 project applicable to the following device models:

Model Name / Product Code	Marketing Name
37158	AirSense 10 CPAP
37159	AirSense 10 Elite
37160	AirSense 10 AutoSet
37161	AirSense 10 AutoSet FH
37162	AirCurve 10 ASV
37163	AirCurve 10 S
37164	AirCurve 10 VAuto
37165	AirCurve 10 ST

## 2. Equivalence

All the device models mentioned in [Section 1](#) have the following **same** components/parts:

- Hardware
  - Cellular
  - Flow Generator
- Software
  - Cellular
  - Flow Generator

## 3. Differences

Different models are derived for marketing purposes. The difference between these models are the medical therapy modes and region/customer centric features supported by the device. These settings are configured at manufacturing via a configuration file and is independent of the cellular and flow generator software.

In summary, the differences between the device models mentioned in [Section 1](#) are:

- Therapy modes
- Data settings
- UI language
- Name plaque
- Device colour