

BIOTRONIK GmbH & Co. KG • Postfach 47 02 55 • 12311 Berlin

Federal Communications Commission  
Equipment Authorization Branch  
7435 Oakland Mills Road  
Columbia, MD 21046  
USA

Berlin, September 20, 2006

Re: Label Sample FCC ID: QRILUMAXT

# Lumax 300 VR-T

**BIOTRONIK**

IS-1  
DF-1  
DF-1



REF 355270  
Cat. No.

SN 12345678  
Serial No.

2007-09  
Use Before

41-113 °F  
Star. Temp.



[01]04035479094454(17)070920(21)12345678

Implantable Cardioverter Defibrillator



# Lumax 300 VR-T

PID 89  
Prod. ID

Caution: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician (or properly licenced practitioner)



FCC ID: QRILUMAXT



For single use only, do not resterilize. Do not use if packaging is damaged. Refer to technical manual for instructions. ICD is shipped with all therapy disabled.

BIOTRONIK GmbH & Co. KG  
Weermarkstraße 1  
12359 Berlin - Germany  
Tel: +49 (0) 30 68905-0  
Fax: +49 (0) 30 6852804  
sales@biotronic.com

**BIOTRONIK**

# Lumax 300 DR-T

**BIOTRONIK**



REF 355266  
Cat. No.

SN 12345678  
Serial No.

2007-09  
Use Before

41-113 °F  
Stor. Temp.




[01]04035479094416[17]070920[21]12345678

Implantable Cardioverter Defibrillator



## Lumax 300 DR-T

PID 89  
Prod. ID

Caution: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician (or properly licenced practitioner)



FCC ID: QRILUMAXT



For single use only, do not resterilize. Do not use if packaging is damaged. Refer to technical manual for instructions. ICD is shipped with all therapy disabled.

BIOTRONIK GmbH & Co. KG  
Wormannstraße 1  
12359 Berlin - Germany  
Tel +49 (0) 30 68906-0  
Fax +49 (0) 30 6852904  
sales@biotronic.com

**BIOTRONIK**

# Lumax 300 HF-T

**BIOTRONIK**

IS-1  
IS-1  
IS-1  
DF-1  
DF-1



REF 355262  
Cat. No.

SN 12345678  
Serial No.

2007-09  
Use Before

41-113 °F  
Star. Temp.




10110403547909437911707092021112345678

Cardiac Resynchronization Therapy Defibrillator



## Lumax 300 HF-T

PID 89  
Prod. ID

Caution: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician (or properly licenced practitioner)



FCC ID: QRILUMAXT



For single use only, do not resterilize. Do not use if packaging is damaged. Refer to technical manual for instructions. ICD is shipped with all therapy disabled.

BIOTRONIK GmbH & Co. KG  
Wormannstraße 1  
12359 Berlin - Germany  
Tel. +49 (0) 30 68905-0  
Fax +49 (0) 30 6852804  
sales@biotronic.com

**BIOTRONIK**

# Lumax 340 VR-T



IS-1  
DF-1  
DF-1



REF 355271

Can. No.

SN 12345678

Serial No.

2007-09

Use Before

41-113 °F

Stor. Temp.



(01)04035479094461(17)070920(21)12345678


Implantable Cardioverter Defibrillator

# Lumax 340 VR-T


PID 89

Prod. ID

Caution: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician (or properly licenced practitioner)



FCC ID: QRILUMAXT



For single use only, do not resterilize. Do not use if packaging is damaged.  
Refer to technical manual for instructions. ICD is shipped with all therapy disabled.

BIOTRONIK GmbH & Co. KG  
Wormannstraße 1  
12359 Berlin - Germany  
Tel: +49 (0) 30 68905-0  
Fax: +49 (0) 30 6852804  
sales@biotronic.com



# Lumax 340 DR-T

**BIOTRONIK**

IS-1  
IS-1  
DF-1  
DF-1



REF 355267

Cat. No.

SN 12345678

Serial No.

2007-09

Use Before

41-113 °F

Stor. Temp.



[01]04035479094423(17)070920(21)12345678


Implantable Cardioverter Defibrillator

## Lumax 340 DR-T



PID 89

Prod. ID

Caution: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician (or properly licenced practitioner).



FCC ID: QRILUMAXT



BIOTRONIK GmbH & Co. KG  
Wormannkehl 1  
12359 Berlin - Germany  
Tel +49 (0) 30 68985-0  
Fax +49 (0) 30 6852804  
sales@biotronic.com

**BIOTRONIK**

For single use only, do not resterilize. Do not use if packaging is damaged.  
Refer to technical manual for instructions. ICD is shipped with all therapy disabled.

# Lumax 340 HF-T



IS-1  
IS-1  
IS-1  
DF-1  
DF-1



REF 355263  
Cat. No.

SN 12345678  
Serial No.

2007-09  
Use Before

41-113 °F  
Stor. Temp.




[01]04035479094386[17]070920[21]12345678

Cardiac Resynchronization Therapy Defibrillator



# Lumax 340 HF-T

PID 89  
Prod. ID

Caution: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician (or properly licenced practitioner)



FCC ID: QRILUMAXT



BIOTRONIK GmbH & Co. KG  
Wormannkehl 1  
12359 Berlin - Germany  
Tel: +49 (0) 30 68905-0  
Fax: +49 (0) 30 6852604  
sales@biotronic.com



For single use only, do not resterilize. Do not use if packaging is damaged.  
Refer to technical manual for instructions. ICD is shipped with all therapy disabled.