

#whdentalwerk



iodent.wh.com

Now at
your dealer
or **wh.com**



Quick Installation Guide – ioDent[®] system



1

Connect the ioDent[®]
gateway mini to
power supply



2

Plug in the
ioDent[®] Wi-Fi dongle
into the device



3

Connect your device
to power supply



4

Follow the link sent to
your e-mail address or
follow-up ioDent.com



Enter e-mail
address and
password

5

Run the initial
set up on the
ioDent.com platform



6

You have successfully
connected your device

W&H Dentalwerk Bürmoos GmbH
Ignaz-Glaser-Straße 53, Postfach 1
5111 Bürmoos, **Austria**

t +43 6274 6236-0
f +43 6274 6236-55
office@wh.com **wh.com**



51026 AEN Rev. 000 / 26.03.2019
Subject to alterations