

④ Enter the Wi-Fi settings menu of the

phone to find and connect with the wifi named Smart Oscscope.

④ Open the app once the wifi is connected. Observe the condition within the ear through the otoscope and APP in a real-time manner.

Image Library
Device controls

IMPORTANT NOTE:

For a Class B digital device, pursuant to part 25 of the FCC Rules, These I are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with instructions, may cause harmful interference to radio communication. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to you,

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from the one to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

— 60 —