

OneTouch
Verio Flex®

⬇ + ⌂ = Edit date/time • Editar la fecha/hora

⬆ + ⬇ = Edit range • Editar rango

⬆ + ⌂ = ⚡ On/Off • Encendido/
apagado



AW 06932105A

||||| SN: XXXXXXXX

Blood Glucose Meter
Medidor de glucosa en la sangre

1 888 567-3005 (English) • 1 888 567-3010 (Español)

In vitro diagnostic. For self-testing. • Uso diagnóstico
in vitro. Para autodiagnóstico.

LifeScan Europe
6300 Zug, Switzerland

Made in China •
Fabricado en China

FCC ID: 2AC15-202 Rev. Date: /
CR0002, SV  Fecha rev: 06/2012

© 2008-2012 LifeScan IP Holdings, LLC
AW 06932105A