

TEST REPORT

Test Report No.: UL-RPT-RP12761350-116C

Customer Remote Diagnostic Technologies Ltd

Model No. / HVIN 00-1026-R

PMN Tempus Pro

FCC ID Contains FCC ID: NCMOMO6012 & PV7-WIBEAR11N-DF2

ISED Certification No. Contains IC: 2734A-M06012 & 7738A-WB11NDF2

Technology WLAN

Test Standard(s) FCC Parts 15.209(a) & 15.247(d)

Innovation, Science and Economic Development Canada

RSS-247 Issue 2 February 2017 Section 5.5 & RSS-Gen Issue 5 April 2018 Section 6.13

This test report shall not be reproduced except in full, without the written approval of UL VS LTD. 1.

- The results in this report apply only to the sample(s) tested. 2.
- The sample tested is in compliance with the above standard(s). 3.
- The test results in this report are traceable to the national or international standards. 4.

Version 1.0 5.

> Date of Issue: 02 May 2019

Checked by:

Ben Mercer

Senior Test Engineer, Radio Laboratory

Company Signatory:

Sarah Williams

Senior Test Engineer, Radio Laboratory

UL VS LTD



This laboratory is accredited by UKAS. The tests reported herein have been performed in accordance with its' terms of accreditation.

4) Blang

Telephone: +44 (0)1256 312000 Facsimile: +44 (0)1256 312001

Customer Information

Company Name:	Remote Diagnostic Technology Ltd	
Address:	Pavilion C2 Ashwood Park Ashwood way Basingstoke RG23 8BG United Kingdom	

Report Revision History

Version Number	I ISSUA I ISTA RAVISION I ISTAIIS		Revised By
1.0	02/05/2019	Initial Version	Ben Mercer

Page 2 of 19 UL VS LTD

Table of Contents

Customer Information	2
Report Revision History	2
1. Attestation of Test Results	4
1.1. Description of EUT	4
1.2. General Information	4
1.3. Summary of Test Results	4
1.4. Deviations from the Test Specification	4
2. Summary of Testing	5
2.1. Facilities and Accreditation	5
2.2. Methods and Procedures	5
2.3. Calibration and Uncertainty	6
2.4. Test and Measurement Equipment	7
3. Equipment Under Test (EUT)	8
3.1. Identification of Equipment Under Test (EUT)	8
3.2. Modifications Incorporated in the EUT	8
3.3. Additional Information Related to Testing	9
3.4. Description of Available Antennas	9
3.5. Description of Test Setup	10
4. Radiated Test Results	16
4.1. Transmitter Radiated Emissions <1 GHz	16
4.2. Transmitter Radiated Emissions >1 GHz	18

UL VS LTD Page 3 of 19

1. Attestation of Test Results

1.1. Description of EUT

The equipment under test was a medical vital signs monitor that contains FCC / ISED Canada certified GSM/UMTS, *Bluetooth* and 2.4 GHz WLAN radio modules (FCC ID: NCMOMO6012 & PV7-WIBEAR11N-DF2, ISED Certification No. IC: 2734A-M06012 & 7738-WB11NDF2).

1.2. General Information

Specification Reference:	47CFR15.247	
Specification Title:	Code of Federal Regulations Volume 47 (Telecommunications): Part 15 Subpart C (Intentional Radiators) - Section 15.247	
Specification Reference:	47CFR15.209	
Specification Title:	Code of Federal Regulations Volume 47 (Telecommunications): Part 15 Subpart C (Intentional Radiators) - Sections 15.209	
Specification Reference:	RSS-Gen Issue 5 April 2018	
Specification Title:	General Requirements for Compliance of Radio Apparatus	
Specification Reference:	RSS-247 Issue 2 February 2017	
Specification Title:	Digital Transmission Systems (DTSs), Frequency Hopping Systems (FHSs) and Licence-Exempt Local Area Network (LE-LAN) Devices	
FCC Test Firm Registration No.:	621311	
ISED#:	20903	
CAB Identifier:	UK0001	
Location of Testing:	UL VS LTD, Unit 3 Horizon, Wade Road, Kingsland Business Park, Basingstoke, Hampshire, RG24 8AH, United Kingdom	
Test Dates:	16 March 2019 to 22 March 2019	

1.3. Summary of Test Results

FCC Reference (47CFR)	ISED Canada Reference	Measurement	Result
15.247(d) & 15.209(a)	RSS-Gen 6.13 & RSS-247 5.5	Transmitter Radiated Emissions	②
Key to Results			•
= Complied	= Did not comply		

1.4. Deviations from the Test Specification

For the measurements contained within this test report, there were no deviations from, additions to, or exclusions from the test specification identified above.

Page 4 of 19 UL VS LTD

2. Summary of Testing

2.1. Facilities and Accreditation

The test site and measurement facilities used to collect data are located at Unit 3 Horizon, Wade Road, Kingsland Business Park, Basingstoke, Hampshire, RG24 8AH, United Kingdom. The following table identifies which facilities were utilised for radiated emission measurements documented in this report. Specific facilities are also identified in the test results sections.

Site 1	
Site 2	
Site 17	Х

UL VS LTD is accredited by UKAS. The tests reported herein have been performed in accordance with its terms of accreditation.

2.2. Methods and Procedures

Reference:	ANSI C63.10-2013
Title:	American National Standard of Procedures for Compliance Testing of Unlicensed Wireless Devices

UL VS LTD Page 5 of 19

2.3. Calibration and Uncertainty

Measuring Instrument Calibration

In accordance with UKAS requirements all the measurement equipment is on a calibration schedule. All equipment was within the calibration period on the date of testing.

Measurement Uncertainty

No measurement or test can ever be perfect and the imperfections give rise to error of measurement in the results. Consequently the result of a measurement is only an approximation to the value measured (the specific quantity subject to measurement) and is only complete when accompanied by a statement of the uncertainty of the approximation.

The expression of uncertainty of a measurement result allows realistic comparison of results with reference values and limits given in specifications and standards.

The uncertainty of the result may need to be taken into account when interpreting the measurement results.

The reported expanded uncertainties below are based on a standard uncertainty multiplied by an appropriate coverage factor such that a confidence level of approximately 95% is maintained. For the purposes of this document "approximately" is interpreted as meaning "effectively" or "for most practical purposes".

Measurement Type	Range	Confidence Level (%)	Calculated Uncertainty
Radiated Spurious Emissions	30 MHz to 1 GHz	95%	±4.65 dB
Radiated Spurious Emissions	1 GHz to 25 GHz	95%	±2.94 dB

The methods used to calculate the above uncertainties are in line with those recommended within the various measurement specifications. Where measurement specifications do not include guidelines for the evaluation of measurement uncertainty the published guidance of the appropriate accreditation body is followed.

Page 6 of 19 UL VS LTD

2.4. Test and Measurement Equipment

Test Equipment Used for Transmitter Radiated Emissions

Asset No.	Instrument	Manufacturer	Туре No.	Serial No.	Date Calibration Due	Cal. Interval (Months)
M2003	Thermoygrometer	Testo	608-H1	45046641	06 Jan 2020	12
K0017	3m RSE Chamber	Rainford	N/A	N/A	16 Feb 2020	12
M1995	Test Receiver	Rohde & Schwarz	ESU40	100428	10 Aug 2019	12
A3167	Pre-Amplifier	Com-Power	PAM-103	18020010	11 Feb 2020	12
A2948	Pre-Amplifier	Com-Power	PAM-118A	551087	12 Feb 2020	12
A3142	Pre-Amplifier	Schwarzbeck	BBV 9718 B	00020	12 Feb 2020	12
A2893	Pre-Amplifier	Schwarzbeck	BBV 9721	9721-021	15 Feb 2020	12
A3161	Antenna	Teseq	CBL611D	50859	17 Dec 2019	12
A2889	Antenna	Schwarzbeck	BBHA 9120 B	BBHA 9120 B 563	12 Feb 2020	12
A2890	Antenna	Schwarzbeck	HWRD 750	014	12 Feb 2020	12
A2892	Antenna	Schwarzbeck	BBHA 9170	9170-727	16 Feb 2020	12
A3113	Attenuator	AtlanTecRF	AN18-06	219706#3	17 Dec 2019	12
A2916	Attenuator	AltalTecRF	AN18W5-10	832827#1	20 Feb 2020	12
A2131	Low Pass Filter	AtlanTecRF	AFL-02000	JFB1004- 002	20 Feb 2020	12
A2914	High Pass Filter	AltalTecRF	AFH-03000	2155	20 Feb 2020	12
A2947	High Pass Filter	AltalTecRF	AFH-07000	1601900001	20 Feb 2020	12

UL VS LTD Page 7 of 19

3. Equipment Under Test (EUT)

3.1. Identification of Equipment Under Test (EUT)

Brand Name:	Remote Diagnostic Technology Ltd	
Model No. / HVIN:	00-1026-R	
PMN:	Tempus Pro	
Test Sample Serial Number:	602429	
Hardware Version:	Trizeps VII	
Software Version:	V7.01	
FCC ID:	Contains FCC ID: NCMOMO6012 & PV7-WIBEAR11N-DF2	
ISED Certification Number:	Contains IC: 2734A-M06012 & 7738A-WB11NDF2	

3.2. Modifications Incorporated in the EUT

No modifications were applied to the EUT during testing.

Page 8 of 19 UL VS LTD

3.3. Additional Information Related to Testing

Technology Tested:	WLAN (IEEE 802.11b,g,n) / Digital Transmission System				
Type of Unit:	Transceiver	Transceiver			
Modulation Type:	DQPSK				
Data Rates:	802.11b	11 Mbps			
Power Supply Requirement(s):	Nominal 12.0 VDC via 120 VAC 60 Hz				
Channel Spacing:	20 MHz				
Transmit Frequency Range:	2412 MHz to 2472 MHz				
Transmit Channels Tested:	Channel ID Channel Number Channel Frequency (MHz)				
	Bottom 1 2412				
	Middle	6	2437		
	Тор	11	2462		

3.4. Description of Available Antennas

The radio utilises an integrated antenna with the following maximum gain:

Manufacturer	Model	Туре	Frequency Range (MHz)	Antenna Gain (dBi)
Antenova	A10470-BLADE	Monopole	2400 to 2500	Not specified

UL VS LTD Page 9 of 19

3.5. Description of Test Setup

EUT Accessories

Description: Power Supply Brand Name: ASTEC Model Name or Number: DPS53-M / RDT Tempus PSU 01-2049 Serial Number: Not marked or stated Description: Earpiece. Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Pulse Oximeter (Masimo rainbow finger sensor). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: End Tidal CO2 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Serial Number: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Model Name or Number: Not marked or stated	The following accessories were used to exercise the EUT during testing:		
Model Name or Number: DPS53-M / RDT Tempus PSU 01-2049 Serial Number: Not marked or stated Description: Earpiece. Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: Pulse Oximeter (Masimo rainbow finger sensor). Quantity 1 Brand Name: Not marked or stated Description: Pulse Oximeter (Masimo rainbow finger sensor). Quantity 1 Brand Name: Not marked or stated Not marked or stated Description: End Tidal C02 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Description: Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Not marked or stated Model Name or Number: Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated	Description:	Power Supply	
Serial Number: Description: Earpiece. Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: Pulse Oximeter (Masimo rainbow finger sensor). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Description: End Tidal C02 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Model Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Model Name: Not marked or stated Model Name: Not marked or stated Model Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Model Name: Not marked or stated	Brand Name:	ASTEC	
Description: Earpiece. Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Pulse Oximeter (Masimo rainbow finger sensor). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Berial Number: Not marked or stated Description: End Tidal C02 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name: Not marked or stated Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated	Model Name or Number:	DPS53-M / RDT Tempus PSU 01-2049	
Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Pulse Oximeter (Masimo rainbow finger sensor). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: End Tidal CO2 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Model Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Description: Not marked or stated Description: Not marked or stated	Serial Number:	Not marked or stated	
Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Pulse Oximeter (Masimo rainbow finger sensor). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: End Tidal CO2 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Model Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Description: Not marked or stated Description: Not marked or stated	B	Francisco O continua	
Model Name or Number: Not marked or stated Pulse Oximeter (Masimo rainbow finger sensor). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Not marked or stated Pescription: End Tidal C02 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Pescription: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Pescription: Brand Name: Not marked or stated Not marked or stated Not marked or stated Model Name or Number: Not marked or stated Not marked or stated Not marked or stated Pescription: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	•		
Description: Pulse Oximeter (Masimo rainbow finger sensor). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: End Tidal CO2 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: Description: Description: Not marked or stated Not marked or stated Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Model Name: Not marked or stated Description: Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated			
Description: Pulse Oximeter (Masimo rainbow finger sensor). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: End Tidal C02 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Description: Not marked or stated Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Not marked or stated Description: Not marked or stated Not marked or stated Not marked or stated Description: Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated			
Brand Name: Model Name or Number: Not marked or stated Not marked or stated Description: End Tidal C02 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: Not marked or stated Description: Not marked or stated Not marked or stated Description: Not marked or stated	Serial Number:	Not marked or stated	
Model Name or Number: Not marked or stated Not marked or stated Description: End Tidal C02 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Not marked or stated Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Description:	Pulse Oximeter (Masimo rainbow finger sensor). Quantity 1	
Description: End Tidal C02 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Description: Not marked or stated Description: Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Brand Name:	Not marked or stated	
Description: End Tidal C02 Monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Pescription: Brand Name: Not marked or stated ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Description: Not marked or stated Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Model Name or Number:	Not marked or stated	
Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: Not marked or stated Serial Number: Not marked or stated Description: Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Serial Number:	Not marked or stated	
Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: Not marked or stated Serial Number: Not marked or stated Description: Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated			
Model Name or Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Model Name or Number: Not marked or stated Description: Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Description:	End Tidal C02 Monitor (plastic hose). Quantity 1	
Serial Number: Not marked or stated Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Brand Name:	Not marked or stated	
Description: ECG (RDT 12-lead cable 01-2073). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Description: Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Model Name or Number:	Not marked or stated	
Brand Name: Model Name or Number: Not marked or stated Not marked or stated Not marked or stated Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Serial Number:	Not marked or stated	
Brand Name: Model Name or Number: Not marked or stated Not marked or stated Not marked or stated Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated			
Model Name or Number: Not marked or stated Not marked or stated Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Description:	ECG (RDT 12-lead cable 01-2073). Quantity 1	
Serial Number: Not marked or stated Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Model Name or Number: Not marked or stated Serial Number: Not marked or stated Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Brand Name:	Not marked or stated	
Description: Non-invasive blood pressure monitor (plastic hose). Quantity 1 Brand Name: Not marked or stated Not marked or stated Serial Number: Not marked or stated Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Model Name or Number:	Not marked or stated	
Brand Name: Model Name or Number: Not marked or stated Serial Number: Not marked or stated Not marked or stated Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	Serial Number:	Not marked or stated	
Brand Name: Model Name or Number: Not marked or stated Serial Number: Not marked or stated Not marked or stated Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	B	No. 12 and a law and a second a	
Model Name or Number: Not marked or stated Not marked or stated Description: Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated	•		
Serial Number: Not marked or stated Description: Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated			
Description: Contact temperature thermocouple. Quantity 2 Brand Name: Not marked or stated			
Brand Name: Not marked or stated	Serial Number:	Not marked or stated	
	Description:	Contact temperature thermocouple. Quantity 2	
Madel News or Newsham	· · · · · · · · · · · · · · · · · · ·		
Model Name or Number: Not marked or stated	Model Name or Number:	Not marked or stated	
Serial Number: Not marked or stated	Serial Number:	Not marked or stated	

Page 10 of 19 UL VS LTD

EUT Accessories (continued)

Description:	USB Cable (RDT data transfer cable 01-2243). Quantity 1
Brand Name:	Not marked or stated
Model Name or Number:	Not marked or stated
Serial Number:	Not marked or stated

Description:	Invasive blood pressure; internal (RDT 2-channel cable 01-2108 & 01-2113 with transducers)
Brand Name:	Not marked or stated
Model Name or Number:	Not marked or stated
Serial Number:	Not marked or stated

Support Equipment

The following support equipment was used to exercise the EUT during testing:

Description:	Laptop PC
Brand Name:	Lenovo
Model Name or Number:	L430
Serial Number:	R9-Z2L03 13/06

Description:	USB Cable. Length 2 metres. Quantity 1.
Brand Name:	Not marked or stated
Model Name or Number:	Not marked or stated
Serial Number:	Not marked or stated

Description:	Ethernet Cable. Length 2 metres. Quantity 1.
Brand Name:	Not marked or stated
Model Name or Number:	Not marked or stated
Serial Number:	Not marked or stated

Description:	USB Hub
Brand Name:	НАМА
Model Name or Number:	00078498
Serial Number:	Not marked or stated

Description:	Ethernet Hub
Brand Name:	Netgear
Model Name or Number:	DG834G
Serial Number:	1JX167B008C4A

UL VS LTD Page 11 of 19

Support Equipment (continued)

Description:	Wideband Radio Comms Tester
Brand Name:	Rohde & Schwarz
Model Name or Number:	CMW 500
Serial Number:	145923

Operating Modes

The EUT was tested in the following operating mode(s):

• Continuously transmitting with a modulated carrier at maximum power on the relevant channels as required using the supported data rates/modulation types.

Configuration and Peripherals

The EUT was tested in the following configuration(s):


- WLAN was configured using a command prompt application installed on the laptopc PC supplied by the customer. The application was used to enable continuous transmission and to select the test channels as required.
- Transmitter radiated spurious emissions were performed with the EUT transmitting with a worst case data rate of 11 Mbps as declared by the client.
- The EUT was powered from a 120 VAC 60 Hz single phase supply via a 12V power adaptor.
- All ports were terminated with typical end-user hardware.

Page 12 of 19 UL VS LTD

Test Setup Diagrams


Radiated Tests:

Test Setup for Transmitter Radiated Emissions




UL VS LTD Page 13 of 19

Test Setup for Transmitter Radiated Emission (continued)



Page 14 of 19 UL VS LTD

Test Setup for Transmitter Radiated Emission (continued)



UL VS LTD Page 15 of 19

4. Radiated Test Results

4.1. Transmitter Radiated Emissions <1 GHz

Test Summary:

Test Engineer:	Mark Perry	Test Date:	22 March 2019
Test Sample Serial Number:	602429		

FCC Reference:	Parts 15.247(d) & 15.209(a)
ISED Canada Reference:	RSS-Gen 6.13 & RSS-247 5.5
Test Method Used:	ANSI C63.10 Sections 6.3 and 6.5
Frequency Range	30 MHz to 1000 MHz

Environmental Conditions:

Temperature (°C):	21 to 24
Relative Humidity (%):	37 to 51

Note(s):


- 1. The final measured value, for the given emission, in the table below incorporates the calibrated antenna factor and cable loss.
- The preliminary scans showed similar emission levels below 1 GHz, for each channel of operation. Therefore final radiated emissions measurements were performed with the EUT set to the middle channel only.
- 3. Measurements below 1 GHz were performed in a semi-anechoic chamber (Asset Number K0017) at a distance of 3 metres. The EUT was placed at a height of 80 cm above the reference ground plane in the centre of the chamber turntable. Maximum emission levels were determined by height searching the measurement antenna over the range 1 metre to 4 metres.
- 4. Pre-scans were performed with the EUT transmitting on middle channel with a data rate of 802.11b 11 Mbit/s and markers placed on the highest measured levels. The test receiver resolution bandwidth was set to 120 kHz and video bandwidth 500 kHz. A peak detector was used, sweep time was set to auto and trace mode was Max Hold.
- 5. Final measurements were performed on the marker frequencies and the results entered into the table below. The test receiver resolution bandwidth was set to 120 kHz, using a CISPR quasi-peak detector and span wide enough to see the whole emission.

Page 16 of 19 UL VS LTD

Transmitter Radiated Emissions (continued)

Results: Middle Channel / 802.11b / 11Mbps

Frequency (MHz)	Antenna Polarity	Level (dBμV/m)	Limit (dBμV/m)	Margin (dB)	Result
112.840	Vertical	30.3	43.5	13.2	Complied
250.119	Horizontal	30.1	46.0	15.9	Complied
331.780	Vertical	30.6	46.0	15.4	Complied
609.043	Horizontal	31.1	46.0	14.9	Complied
613.903	Horizontal	36.4	46.0	9.6	Complied
990.024	Vertical	37.2	54.0	16.8	Complied



Note: This plot is a pre-scan and for indication purposes only. For final measurements, see accompanying table.

UL VS LTD Page 17 of 19

ISSUE DATE: 02 MAY 2019

VERSION 1.0

4.2. Transmitter Radiated Emissions >1 GHz

Test Summary:

Test Engineer:	Mark Perry	Test Dates:	16 March 2019 & 17 March 2019
Test Sample Serial Number:	602429		

FCC Reference:	Parts 15.247(d) & 15.209(a)	
ISED Canada Reference:	RSS-Gen 6.13 & RSS-247 5.5	
Test Method Used:	ANSI C63.10 Sections 6.3, 6.6, 11.11 & 11.12	
Frequency Range	1 GHz to 25 GHz	

Environmental Conditions:

Temperature (°C):	21 to 24
Relative Humidity (%):	37 to 43

Note(s):

- 1. The final measured value, for the given emission, in the table below incorporates the calibrated antenna factor and cable loss.
- 2. The emission shown on the 1 GHz to 3 GHz plot is the EUT fundamental.
- 3. No spurious emissions were detected above the noise floor of the measuring receiver therefore the highest peak and average noise floor readings of the measuring receiver were recorded as shown in the tables below.
- 4. Pre-scans were performed in a fully anechoic chamber (Asset Number K0017) at a distance of 3 metres. The EUT was placed at a height of 1.5 metres above the test chamber floor in the centre of the chamber turntable. All measurement antennas were placed at a fixed height of 1.5 metres above the test chamber floor, in line with the EUT.
- 5. Pre-scans were performed with the EUT transmitting on middle channel with a data rate of 802.11b 11 Mbit/s and a marker placed on the highest measured level of the appropriate plot. The test receiver resolution bandwidth was set to 1 MHz and video bandwidth 3 MHz. The sweep time was set to auto.

Results: Peak


Frequency (MHz)	Antenna Polarity	Peak Level (dBμV/m)	Peak Limit (dBμV/m)	Margin (dB)	Result
2881.500	Vertical	54.9	74.0	19.1	Complied


Results: Average


Frequency (MHz)	Antenna Polarity	Average Level (dBμV/m)	Average Limit (dBμV/m)	Margin (dB)	Result
1033.000	Vertical	47.8	54.0	6.2	Complied


Page 18 of 19 UL VS LTD

Transmitter Radiated Emissions (continued)









--- END OF REPORT ---

UL VS LTD Page 19 of 19