

Electronic Blood Pressure Monitor Model: HWA11
Power: DC 3V(AAAx2) 1.5W
Measuring Range: 0~295 mmHg
FCC ID:2ADYL-HWA11
This device complies with Part 15 of the FCC Rules.
Operation is subject to the following two conditions:
(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Shenzhen Jumper Medical Equipment Co., Ltd.
D Building, No. 71, Xintian Road, Fuyong
Street, Baoan, Shenzhen, Guangdong, China

EC REP MedPath GmbH
Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany