## **Declaration on model difference**

We the undersigned hereby confirm that these models are identical, and all models have the same RF module, PCB layout, Antenna, schematics and component. Only different on model name.

The layout of the left and right earphones is the same, but the materials used are slightly different. The left side of the right earphone has a 1K resistor more than that of the left earphone, so as to generate a pull-down signal when the software is connected. The side of the earphone that generates the signal is the right earphone by default, and this component will not affect the test. The left earphone has one more output filter capacitor than the right earphone. Since the left earphone is the host, the audio signal is first sent to the left earphone, and then sent to the right earphone after filtering by the filter capacitor. The right earphone just receives the signal without transmitting it, so there is no filter capacitor, which has no impact on the test.

Production name	Trademark	FCC	Model no.
Vivitar Pro True Wireless	N/A	2AMD8-PROW	TWSPRO-WHT,
Earbuds			SKU#2198530

Confirmed by:

James Cao

James Cao Shenzhen Ground Enterprises Co., Ltd 2023-05-12