

| | |
|----------|-------------|
| SPEC NO. | eLB1392-074 |
| REV | 5.0 |
| PAGE | 1 OF 4 |

Nisha Rajbhandari
27 February 2021

Common labels/moulded in artwork




BRAND
ResMed branding
chrome or colour foiled
on front of product



“BF APPLIED PART”
Body Fibrillation Icon
Moulded into plastic case
using an insert
on rear of device,
adjacent to the device outlet



“R390-411 Label - MR Unsafe”
White pre-printed label stuck next to
power connector on rear of assembly



“R390-429 Label - FCC & IC NIR M1”
White pre-printed label
stuck on bottom of assembly

ResMed Pty Ltd
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia
See www.resmed.com for other
ResMed locations worldwide

ResMed SAS
Parc Technologique de Lyon
292 Allée Jacques Monod
69791 Saint Priest Cedex France

STORAGE & TRANSPORT RANGES

-25C

-13F

70C

158F

5

95

0

3010m

9870ft

ALTITUDE LIMITATION

5C

41F

35C

95F

10

95

OPERATING RANGES

RTCA/DO-160 Section 21

Category M Compliant

& FAA Compliant

MD

i

IP22

For patent information, see www.resmed.com/ip
See symbols glossary at www.resmed.com/symbols

A390-7374/2

A390-7374 Bottom Case Assembly
Laser-Etched on bottom of chassis

| | |
|----------|-------------|
| SPEC NO. | eLB1392-074 |
| REV | 5.0 |
| PAGE | 2 OF 4 |


Nisha Rajbhandari
27 February 2021

Product Specific Labels for FCC




39001 AirSense 11 AutoSet USA CO
(Base Code: 39485)

R390-7396




39002 AirSense 11 CPAP USA CO
(Base Code: 39486)

R390-7397



39001 AirSense 11 Elite USA CO
(Base Code: 39487)

R390-7398




PRODUCT NAME PLAQUES
Foiled product variant badges
“permanently” attached to front of product

SERIAL NUMBER LABELS
Thermally printed labels
attached to rear of product




39001 AirSense 11 AutoSet USA TRI
(Base Code: 39485)

R390-7396




39002 AirSense 11 CPAP USA TRI
(Base Code: 39486)

R390-7397



39001 AirSense 11 Elite USA TRI
(Base Code: 39487)

R390-7398



NOTE:

- The device ID is the base code, which can be packaged in different boxes of accessories
- CO means it is packed with a humidifier.
- TRI means it is packed with a humidifier and heated tube.

| | |
|----------|-------------|
| SPEC NO. | eLB1392-074 |
| REV | 5.0 |
| PAGE | 3 OF 4 |


Nisha Rajbhandari
27 February 2021

Product Specific Labels for ISED




39001 AirSense 11 AutoSet CAN CO
(Base Code: 39488)

R390-7396




39002 AirSense 11 CPAP CAN CO
(Base Code: 39489)

R390-7397



39001 AirSense 11 Elite CAN CO
(Base Code: 39490)

R390-7398




PRODUCT NAME PLAQUES
Foiled product variant badges
“permanently” attached to front of product

SERIAL NUMBER LABELS
Thermally printed labels
attached to rear of product




39001 AirSense 11 AutoSet CAN TRI
(Base Code: 39488)

R390-7396




39002 AirSense 11 CPAP CAN TRI
(Base Code: 39489)

R390-7397



39001 AirSense 11 Elite CAN TRI
(Base Code: 39490)


R390-7398




NOTE:
- The device ID is the base code, which can be packaged in different boxes of accessories
- CO means it is packed with a humidifier.
- TRI means it is packed with a humidifier and heated tube.

| | |
|----------|-------------|
| SPEC NO. | eLB1392-074 |
| REV | 5.0 |
| PAGE | 4 OF 4 |

Nisha Rajbhandari
27 February 2021



REAR VIEW



BOTTOM VIEW