



Electronic Blood Pressure Monitor
Model: HWA21
Power: DC 3V(AA x2) 1.5W IP22
Measuring Range: 0~295 mmHg
FCC ID: 2ADYL-HWA21

This device complies with Part 15 of the FCC Rules.
Operation is subject to the following two conditions:
(1) this device may not cause harmful interference, and
(2) this device must accept any interference received,
including interference that may cause undesired operation.

Shenzhen Jumper Medical Equipment Co., Ltd.
D Building, No. 71, Xintian Road, Fuyong
Street, Baoan, Shenzhen, Guangdong, China

MedPath GmbH
Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany

MD V1.0
CE 0482
2022-05
Warning icons: crossed-out trash can, exclamation mark in triangle, person with heart, person with pulse line, and person with heart and pulse line.